



# DHX Heat Exchanger

A new solution to heat transfer

# Sterile, efficient and modular

The Thermo Scientific™ DHX™ Heat Exchanger is a modular system that employs a plate-and-frame concept, using a single-use BioProcess Container (BPC) as the sterile fluid path. The BPCs fit tightly between five dimpled, stainless-steel plates efficiently transferring heat in a counter-current flow path. The system provides efficient, sterile heat transfer that easily integrates into any new or existing process.

## Key features and benefits

- Completely isolated flow paths for process and heat transfer fluid
- Counter-current, serpentine flow patterns
- Dimpled jacketing on the plates provides turbulent flow
- BPCs fill in place with no operator interaction
- Modular design and small overall footprint allows for changing process needs
- Reduced infrastructure requirements
- Reduced processing time
- Improved product consistency


## Applications

### Upstream

- Media hold
- Mixing and fermentation
- Cell separation/protein harvest
- Harvest cooling
- Harvest hold

### Downstream

- Harvest hold
- Buffers
- Protein purification
- Bulk drug substance



# Technical specifications

## DHX BPC assembly and plate system

### Sterile fluid path for process fluid

Each single-use DHX BPC fits tightly between the stainless-steel plate assembly. In a completely isolated flow path, the process fluid flows through the BPCs counter-currently to the heating/cooling fluid within the plates. Once the BPCs are loaded in place, no further operator interaction is required.

DHX BPCs	
General specifications	
<b>Material of construction</b>	Low-density polyethylene ASI™ 26/77 film
<b>Interconnecting tubing</b>	C-Flex
<b>Connections</b>	ReadyMate™ DAC 500 as standard; custom tubing and connections upon request
<b>Flow rate capacity</b>	Up to 15 L/min
<b>Pressure/temperature rating</b> (Installed in DHX plates)	20 psig at 122°F (50°C)


BPC ordering information	Cat. No.
<b>One BPC-26/77 film</b> (DAC connections on outlet ports)	DX00006-I
<b>Two BPC assembly-26/77 film</b> (DAC connections on outlet ports)	DX00007-I
<b>Three BPC assembly-26/77 film</b> (DAC connections on outlet ports)	DX00008-I
<b>Four BPC assembly-26/77 film</b> (DAC connections on outlet ports)	DX00009-I
<b>One BPC-26/77 film</b> (DAC connections and drain tubing)	DX00010-I
<b>Two BPC assembly-26/77 film</b> (DAC connections and drain tubing)	DX00011-I
<b>Three BPC assembly-26/77 film</b> (DAC connections and drain tubing)	DX00012-I
<b>Four BPC assembly-26/77 film</b> (DAC connections and drain tubing)	DX00013-I

DHX bioprocess equipment: plate system	
General specifications	
<b>Material of construction</b>	316L stainless-steel
<b>Effective heat transfer area</b>	Up to 27 sq. ft.
<b>Overall dimensions (WxDxH)</b>	51 x 74 x 69 cm (20 x 29 x 27 in)
<b>Number of plates/BPCs</b>	5 plates/up to 4 BPCs
<b>Dry weight</b>	150 kg (331 lbs)
<b>Full weight</b> (includes 4 BPCs)	190 kg (419 lbs)
<b>Pressure/temperature rating</b>	FV/140 psig at 150°F
<b>Pressure vessel code</b>	ASME U-1
<b>Connections</b>	1/2" compression

Equipment ordering information	Cat. No.
<b>DHX stainless steel unit</b>	DHX1001

Please contact your Thermo Fisher BioProduction sales representative for more information regarding customizations.






## How it works

### Efficient heat transfer

#### Single-use solution for heat transfer

The BPC and the plates each have a matching serpentine pattern. The heat transfer fluid flows through the plate system while the BPCs provide the sterile fluid path for the process fluid, each flowing counter-currently to the other. The temperature differential of the heat transfer fluid provides efficient heat transfer to the process fluid.

Tabs secure  
BPCs in place




# Heat transfer efficiency

## Typical single-pass application

### DHX cooling efficiency - from 37°C (single pass)


- Process fluid: Water at 37°C
- Heat transfer fluid: 30% propylene glycol at 2°C
- Heat transfer flow rate: 15 L/min
- Number of BPCs: 4




Represents the ΔT of process fluid (with a variable flow rate) as a result of a constant 15 L/min flow rate of the cooling fluid.

### DHX heating efficiency - from ambient temperature (single pass)


- Process fluid: water at 20°C
- Heat transfer fluid: 30% propylene glycol at 42°C
- Heat transfer flow rate: 15 L/min
- Number of BPCs: 4



Represents the ΔT of process fluid (with a variable flow rate) as a result of a constant 15 L/min flow rate of the heating fluid.



Cooling efficiency measured by temperature vs. time



Heating efficiency measured by temperature vs. time

## Integrated solutions for bioproduction

### Single-Use Mixers (S.U.M.s)

A variety of options up to 5,000 L for both upstream and downstream applications



### Liquid- and dry-format media

We offer both custom manufacturing and a full range of chemically defined performance media and supplement products



### BioProcess Containers (BPCs)

A variety of configurations up to 2,000 L for liquid harvest, storage, and transportation



### Sera

Our sera are the industry standards for consistent quality and reliability



### Single-Use Bioreactors (S.U.B.s)

50–2,000 L bioreactors capable of integrating with an existing control system



### Buffers and process liquids

Custom and standard buffers and process liquids, including Gibco™ Water for Injection (WFI)—quality water



### Integrity testing systems

A true point-of-use integrity testing system to confirm the integrity of BPCs before use



Find out more at [thermofisher.com/sut](http://thermofisher.com/sut)

Logan, Utah 84321  
United States of America  
+(1) 435-792-8500

Cramlington NE23 1WA  
United Kingdom  
+44 (0) 1670-734093

Millersburg, Pennsylvania 17061  
United States of America  
+(1) 717-692-2104

Matamoros, Tamaulipas 01152  
Mexico  
+011 (52) 868-8109700

**ThermoFisher**  
SCIENTIFIC