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

Biotinylated Human CD20, His, Avitag (HEK293) (Cat. No. CD0-H82E5) captured on Biotin CAP-Series S Sensor Chip can bind Rituximab with an affinity constant of 1.73 nM as determined in a SPR assay (in presence of DDM and CHS) (Biacore T200).

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Immobilized Human CD133 Full Length, His Tag (Nanodisc) (Cat. No. CD3-H52H1) at 1 µg/mL (100 µL/well) can bind Monoclonal Anti-Human CD133 (133-3) Antibody, Human IgG1 with a linear range of 0.2-4 ng/mL.

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Acro BIOSYSTEMS

- Closest to Natural Conformation and Modification
- High Bioactivity & Detection Sensitivity
- High Batch-to-Batch Consistency

MABSOL® Biotinylated Proteins

Bring pre-labeled biotinylated proteins directly to your bench

Chemically Biotin-labeled Proteins
@UltraLys Series

AviTag™ Biotin-labeled Proteins
@PrecisionAvi Series
Suitable for phage display

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

The purity of Biotinylated Human TNF-alpha, epitope tag free, primary amine labeling (active trimer) (MALS verified) (Cat. No. TNA-H8211) is more than 95% in HP-SEC, and the molecular weight of this protein is around 48-65 kDa verified by SEC-MALS.

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
In the above ELISA experiments, three different batches of biotin-labelled human TNF-alpha (Cat. No. [TNA-H82E3](#)) were tested for binding to Adalimumab, and the results showed that the batch-to-batch variability of the tested samples was small.

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FACS analysis shows that the binding of Biotinylated Human PD-L1, Fc, Avitag, His Tag (Cat. No. PD1-H82F3) to 293 overexpressing PD-1 was inhibited by increasing the concentration of neutralising anti-PD-1 antibody. The concentration of PD-L1 used is 2 µg/mL. The IC50 is 0.64 µg/mL (Routinely tested).

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Acro *Star Staining* **Shining Your CAR Detection**

- ★ Using new-generation site-specific labeling technology
- ★ Maintain natural conformation and modification
- ★ High batch-to-batch consistency and uniformity
- ★ Uniquely designed for CAR detection
- ★ More fluorescent labels available

FITC	Alexa Fluor 555
PE	Alexa Fluor 647
APC	Alexa Fluor 488

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

Binding activity of FITC-Labelled Human BCMA and CD19 protein from two different vendors were evaluated by the FACS analysis. The result showed that ACRO's Star Staining FITC-Labelled Human BCMA (Cat. No. [BCA-HF2H3](#)) and CD19 (Cat. No. [CD9-HF2H3](#)) protein have a much higher binding activity than that of the other competitor.

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
No non-specific binding to non-transduced PBMCs between FITC-Labelled Human BCMA Protein of ACRO and competitor. 5e5 of non-transduced PBMCs were stained with FITC-Labelled Human BCMA Protein and anti-CD3 antibody, washed and then analysed with FACS. PE signal was used to evaluate the expression of CD3+ T cells in non-transduced PBMCs, and FITC signal was used to evaluate the non-specific binding activity to non-transduced PBMCs.

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Binding activity of the Human BCMA before and after FITC labelling was evaluated in the above FACS analysis. The result shows that FITC-Labelled BCMA (Cat. No. [BCA-HF2H3](#)) and unconjugated Human BCMA have almost the same level of binding activity.

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Lambda | C.37

Critical Spike Single Mutation

High Frequency Nucleocapsid Mutations

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

The purity of SARS-CoV-2 S protein, His Tag, Super stable trimer (Cat. No. [SPN-C52H9](#)) is more than 90% verified by SDS-PAGE under reducing (R) condition. The molecular weight was around 550-660 kDa confirmed by SEC-MALS. The particles are similar in size and appearance to SARS-CoV-2 trimers reported in published literature verified by negative stain electron microscopy.

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
Different SARS-CoV-2 S proteins were used to detect the sera of recovered patients, among which SARS-CoV-2 S protein, His Tag, Super stable trimer (Cat.No. SPN-C52H9) had a high signal-to-noise ratio.

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

Immobilized SARS-CoV-2 Spike RBD (Omicron, Cat. No. [SPD-C522e](#)) can bind Anti-SARS-CoV-2 Spike RBD Antibody, Mouse IgG1 (Cat. No. [SPD-M305](#)) with a linear range of 0.4-12.5 ng/mL. The antibody does not bind Spike RBD of WT (Cat. No. [SPD-C52H1](#)), Alpha (Cat. No. [SPD-C52Hn](#)), Beta (Cat. No. [SPD-C52Hp](#)), Delta (Cat. No. [SPD-C52Hh](#)) and Gamma (Cat. No. [SPD-C52Hr](#)) (Routinely tested).

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
Immobilized SARS-CoV-2 Spike RBD (Delta, Cat. No. [SPD-C52Hh](#)) can bind Anti-SARS-CoV-2 Spike RBD Antibody, Mouse IgG1 (Cat. No. [SPD-M110](#)) with a linear range of 0.1-3 ng/mL. The antibody does not bind Spike RBD of WT (Cat. No. [SPD-C52H1](#)), Alpha (Cat. No. [SPD-C52Hn](#)), Beta (Cat. No. [SPD-C52Hp](#)), Gamma (Cat. No. [SPD-C52Hr](#)) and Omicron (Cat. No. [SPD-C522e](#)) (QC tested).

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Immobilized SARS-CoV-2 spike protein (Beta, Cat. No. [SPN-C52Hk](#) and Gamma, Cat. No. [SPN-C52Hg](#)) can bind Anti-SARS-CoV-2 Spike Antibody-H6, Mouse IgG1 (AM416) with a linear range of 6-8 ng/mL. The antibody does not bind spike protein of WT (Cat. No. [SPN-C52H7](#)), Alpha (Cat. No. [SPN-C52H6](#)) and Delta (Cat. No. [SPN-C52He](#)) (QC tested).

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Cross-validation data for antibody subtype detection

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


The purity of Human CD19 (20-291), Fc Tag (Cat. No. [CD9-H5251](#)) is more than 95% and the molecular weight of this protein is around 140-160 kDa as verified by SEC-MALS.

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

1e6 of the anti-CD19 CAR-293 cells were stained with 100 μ L of 1:50 dilution (2 μ L stock solution in 100 μ L FACS buffer) of PE-Labeled Human CD19 (20-291), His Tag (Cat. No. [CD9-HP2H3](#)). PE Streptavidin was used as negative control (QC tested).

Binding activity of FITC-Labeled Human CD19, His Tag from two different vendors were evaluated in the ELISA analysis against FMC63 Mab. The result showed that ACRO's FITC-Labeled Human CD19, His Tag has a much higher binding activity than that of the other vendor.

Binding activity of FITC-Labeled Human CD19, His Tag from two different vendors were evaluated in the flow cytometry analysis against anti-CD19-CAR-293 cells. The result showed that ACRO's FITC-Labeled Human CD19, His Tag has a much higher binding activity than that of the other vendor.

Binding activity of Human CD19, His Tag from two different vendors were evaluated by SPR assay against FMC63 Mab. The result showed that ACRO's Human CD19, His Tag can bind FMC63 Mab with an affinity constant of 2.95 nM which is much higher than that of the other vendor.

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


1e6 of the Anti-CD19 CAR-293 cells were stained with 100 μ L of 1:50 dilution (2 μ L stock solution in 100 μ L FACS buffer) of PE-Labelled Monoclonal Anti-FMC63 scFv Antibody, Mouse IgG1 (Cat. No. [FM3-HPY53](#)) and negative control antibody respectively. PE signal was used to evaluate the binding activity (QC tested).

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
5e5 of PBMCs were stained with PE-Labeled Monoclonal Anti-FMC63 scFv Antibody, Mouse IgG1 (Cat. No. [FM3-HPY53](#)) and anti-CD3 antibody, washed followed by FACS analysis. FITC signal was used to evaluate the expression of CD3+ T cells in PBMCs, and PE signal was used to evaluate the non-specific binding activity to PBMCs.

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Immobilized FMC63 scFv at 2 $\mu\text{g}/\text{mL}$ (100 $\mu\text{L}/\text{well}$) can bind Rabbit Anti-Mouse FMC63 scFv Polyclonal Antibody (Cat. No. [FM3-S93](#)) with a linear range of 0.098-0.78 ng/mL (QC tested).

Immobilized anti-CD19 antibody (Clone: FMC63) at 1 $\mu\text{g}/\text{mL}$, add increasing concentrations of Rabbit Anti-Mouse FMC63 scFv Polyclonal Antibody (Cat. No. [FM3-S93](#)) and then add Biotinylated anti-CD19 antibody (Clone: FMC63) at 2 $\mu\text{g}/\text{mL}$. Detection was performed using HRP-conjugated streptavidin with sensitivity of 78 ng/mL (QC tested).

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製品番号	説明
FM3-Y45A1	Monoclonal Anti-FMC63 scFv Antibody, Mouse IgG1 (Y45) (Carrier-free) (recommended for ADA assay) DMF Filed
FM3-Y45	Monoclonal Anti-FMC63 scFv Antibody, Mouse IgG1 (Y45) DMF Filed
FM3-Y45P1	Monoclonal Anti-FMC63 scFv Antibody, Mouse IgG1 (Y45) (HEK293) DMF Filed
FM3-FY45	FITC-Labeled Monoclonal Anti-FMC63 scFv Antibody, Mouse IgG1 (Y45) DMF Filed^{New}
FM3-BY54	Biotinylated Monoclonal Anti-FMC63 scFv Antibody, Mouse IgG1, Avitag™ (Y45) DMF Filed
FM3-HPY53	PE-Labeled Monoclonal Anti-FMC63 scFv Antibody, Mouse IgG1 (Y45) (Site-specific conjugation) DMF Filed^{New}
FM3-S93	Rabbit Anti-Mouse FMC63 scFv Polyclonal Antibody (recommended for ADA assay) (MALS verified)

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

Human CD3E&CD3D Heterodimer Protein (Cat. No. CDD-H52Wa) on SDS-PAGE under reducing (R) and non-reducing (NR) condition and the purity of the protein is greater than 95%. The purity of the protein is more than 85% and the molecular weight of this protein is around 80-90 kDa verified by SEC-MALS.

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
ELISA verified that Human CD3 epsilon & CD3 delta Heterodimer Protein, His Tag & Tag Free (MALS verified) (Cat. No. [CDD-H52W1](#)) can specifically bind to anti-BCMA x CD3 scFv-IgG1 Fc Tag with a linear range of 0.08- 3 ng/mL, does not bind to the isotype control Human IgG1 Fc Protein, Tag Free (MALS verified) (Cat. No. [FCC-H5214](#)).

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BIOSYSTEMS Acro

Antibody-Drug Conjugates (ADCs): the Magic Bullets for Treatment

Power your Magic Bullets

- Target** • 50+ high-quality ADC target proteins
- Linker** • Proteases for cleavable linkers
- PK** • Tools for ADCs PK/PD analysis
- Service** • SPR/BLI analytical & ADA service


Diagram labels: Linker, Payload, Target

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
Anti-LIV-1 mAb captured on CM5 chip via anti-human IgG Fc antibody can bind Human LIV-1, His Tag (Cat. No. [LV1-H5223](#)) with an affinity constant of 5.24 nM as determined in a SPR assay (Biacore T200) (Routinely tested).

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
Competitive experiment of neutralizing SIRP-a antibody

FACS analysis shows that the binding of Human CD47 to ACHN expressing SIRP-a was inhibited by increasing concentration of neutralising SIRP-a antibody. The concentration of Human CD47 used is 20 µg/mL. IC50=9.334 µg/mL.

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


Validated by SDS-PAGE, Human CD32a (R167), His Tag (Cat. No. [CDA-H5221](#)) were more than 95% pure; validated by SEC-MALS, Human CD32a (R167), His Tag (Cat. No. [CDA-H5221](#)) H5221) is more than 90% pure and has a molecular weight of 25-35 kDa.

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SPR verified that Human Fc gamma R1IB / CD32b Protein (Cat. No. [CDB-H5228](#)) immobilised on a CM5 chip could bind to Rituximab (Rituximab) with an affinity constant of 10 μ M (Biacore T200).

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


The purity of Human IL-2RB&IL-2RG Heterodimer Protein, Fc Tag&Fc Tag (MALS verified)(Cat. No. [ILG-H5254](#)) is more than 90% and the molecular weight of this protein is around 140-160 kDa verified by SEC-MALS.

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
Immobilized Human TGFβ1, His Tag (Cat. No. [TGL-H524x](#)) at 1 μg/mL (100 μL/well) can bind Human TGFBR2, Fc Tag (Cat. No. [TG2-H5252](#)) with a linear range of 0.3-2 ng/mL (QC tested).

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
Erbix (Cetuximab) captured on CM5 chip via anti-human IgG Fc antibodies surface, can bind Human EGFR, His Tag (Cat. No. [EGR-H5222](#)) with an affinity constant of 1.3 nM as determined in a SPR assay (Biacore T200).

Loaded Anti-IL23A & IL12B MAb (P40 domain) on Protein A Biosensor, can bind Human IL-12B&IL-12A Heterodimer Protein, His Tag&Flag Tag (Cat. No. [IL2-H4210](#)) with an affinity constant of 1.07 nM as determined in BLI assay (ForteBio Octet Red96e).

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
Bioactivity of three different lots of Human IL-2 Protein, Tag Free (Cat. No. [IL2-H4113](#)) is verified by binding with Human IL-2 R alpha, Fc Tag (Cat. No. [ILA-H5251](#)), and the result shows very high batch-to-batch consistency.

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

TSLPR (Luc) HEK293 Reporter Cell

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
Co-expression analysis of human TSLPR and IL7Rα on TSLPR (Luc) HEK293 Reporter Cell by FACS. Cell surface staining was performed on TSLPR (Luc) HEK293 Reporter Cell or negative control cell using PE-labeled anti-TSLPR antibody and FITC-labeled anti-IL7Rα antibody.


Response to human TSLP protein (RLU). The TSLPR (Luc) HEK293 Reporter Cell was stimulated with serial dilutions of human TSLP protein (ACROBiosystems, Cat.No.TSP-H52Hb). The EC50 was approximately 0.1392 µg/mL.


Response to human TSLP protein (Fold). The TSLPR (Luc) HEK293 Reporter Cell was stimulated with serial dilutions of human TSLP protein (ACROBiosystems, Cat.No.TSP-H52Hb). The max induction fold was approximately 45.15.

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Passage stability analysis of receptors expression by FACS. Flow cytometry surface staining of human TSLPR and IL7Rα on TSLPR (Luc) HEK293 Reporter Cell demonstrates consistent mean fluorescent intensity across passage 10-26. (A) Human IL7Rα expression analysis. (B) Human TSLPR expression analysis.


Passage stability analysis by Signalling Bioassay. The continuously growing TSLPR (Luc) HEK293 Reporter Cell was stimulated with serial dilutions of human TSLP protein. Human TSLP protein stimulated response demonstrates passage stabilization (fold induction and EC50) across passage 10-26.

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製品番号	ホスト	詳細	使用例
CHEK-ATF044	HEK293	VEGFR2 (Luc) HEK293 Reporter Cell	Screen for anti-human VEGF or anti-human VEGFR neutralizing antibody
CHEK-ATF045	HEK293	TSLPR (Luc) HEK293 Reporter Cell	Screen for anti-human TSLP or anti-human TSLPR neutralizing antibody
CJUR-STF046	Jurkat	NFAT (Luc) Jurkat Reporter Cell	Screen for T cell activators
CHEK-ATF047	HEK293	STAT3 (Luc) HEK293 Reporter Cell	Screen for STAT3 activators or inhibitors
CHEK-ATF048	HEK293	NF-κB (Luc) HEK293 Reporter Cell	Screen for NF-κB inhibitors and activators
CHEK-ATF049	HEK293	EGFR (Luc) HEK293 Reporter Cell	Screen for anti-human EGFR or anti-human EGF neutralizing antibody or small molecule inhibitor
CHEK-ATF050	HEK293	NFAT (Luc) HEK293 Reporter Cell	Screen NFAT activators or inhibitors

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

Pre-Coupling Magnetic Beads

- Ready to use
- Minimal non-specific binding
- High-efficiency target capture

SARS-CoV-2 S1	CD3	BCMA
CD19	HER2	PSMA
SARS-CoV-2 RBD	CD33	ROR1
CD38	EpCAM	EGFR

Inhibitor Screening Kits

- High Throughput
- High Specificity
- Easy Process

ELISA Assay Kits for Quantitative Analysis of Therapeutic Antibodies in Serum Samples

PD-1	HER-2	CTLA-4
------	-------	--------

- Low Background
- High Batch-to-Batch Consistency
- Application
- High Stability

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診断用の試薬			
Amyloid Precursor Protein (APP)	Tau-441	α-synuclein	Glial Fibrillary Acidic Protein (GFAP)
Neurofilament Light (NfL)	Neurofilament Heavy (NfH)	S100B	Neuron-specific enolase (NSE)

BIOSYSTEMS Acro

Proteins for Neuroscience

Aneuro

Advance Neuroscience Research

- Neurodegenerative diseases**
 - Alzheimer's disease
 - Parkinson's disease
 - Amyotrophic lateral sclerosis
 - Huntington's disease
- Structural brain diseases**
 - Brain tumors
 - Spinal muscular atrophy
- Neurovascular diseases**
 - Stroke
 - Subarachnoid hemorrhage
- Functional Neurological Disorders**
 - Epilepsy
 - Neuropathic pain
- Neuropsychiatric disorders**
 - Depression
 - Autism spectrum disorder
- Neuroinflammatory diseases**
 - Encephalitis
 - Meningitis

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