Hyperion™

Hyperion XTi Imaging System

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Hyperion XTi Imaging System specifications

Instrument Specifications

Description	Specification
Channels	135
Mass range	75–209 amu
Abundance sensitivity	<0.3% 141Pr (M+1 for 140Ce) <0.3% 174Hf (M–1 for 175Lu)
Frequency	800 pixels/sec
Detection limit	<250 copies per cell
Dynamic range	4.5 orders of magnitude
Calibration	Automated
Operating system	Windows 10 Enterprise LTSC
Data storage	3.6 TB
Cross-cell contamination (redeposition)	≤2% (selected direction)
Crosstalk pixel to pixel	≤25%
Maximum addressable sample size	16 x 58 mm
Brightfield field of view	≥250 x 250 μm
File type	TXT, multipage TIFF, OME-TIFF, OME-Zarr, MCD
Scan rate	Cell Mode: ≥2 mm ² /hr Tissue Mode: >50 mm ² /hr Preview Mode: >700 mm ² /hr
Optimized resolution	Cell Mode: 1 μm Tissue Mode: 5 μm Preview Mode: 1 μm
Lowest step size	Cell Mode: 0.33 μm Tissue Mode: 5 μm Preview Mode (fixed): 25 μm

Data Size

File Type	Size*
MCD	4–5 MB/channel

* 1,000 x 1,000 μm region of interest (ROI) in Cell Mode
* 5,000 x 5,000 μm ROI in Tissue Mode

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FLDM-01130 Rev 03 062025

Hyperion XTi Unlocking New Views Brochure JPN

For Research Use Only. Not for use in diagnostic procedures.

Patent and License Information: www.standardbio.com/legal/notices. Trademarks: www.standardbio.com/legal/trademarks.
Any other trademarks are the sole property of their respective owners. ©2025 Standard BioTools Inc. All rights reserved.

XTi Tissue Imager Dimensions*

Description	Specification
Dimensions	Width: 58 cm (23 in) Height: 132 cm (52 in) Depth: 58 cm (23 in)

Weight	192 kg (423 lb)
--------	-----------------

* The XTi Tissue Imager must be placed directly in front of the CyTOF™ XT instrument.

Slide Loader (installed)

Description	Specification
Dimensions	Width: 47 cm (19 in) Height: 45 cm (18 in) Depth: 19 cm (7 in)

Weight	20 kg (43 lb)
--------	---------------

Capacity	40 slides (20 slides per cassette)
----------	---------------------------------------

Compatible slide dimensions	L = 76 mm +0/–1 mm (75–76 mm) W = 26 mm +0/–1 mm (25–26 mm) T = 1.1 mm +0.1/–0.2 mm (0.9–1.2 mm)
-----------------------------	--

Note that beveled edge slides and chamber slides are not supported.

Supported barcodes	QR, Data Matrix
--------------------	-----------------

Minimum barcode dimensions	5 x 5 mm
----------------------------	----------

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