

# Maxicaps® MR

## Unique Large Scale Single-Use Filter Device

### Product Information

Maxicaps® MR is a fully contained single-use assembly with up to 27 m<sup>2</sup> filtration area, designed for large scale filtration in biopharmaceutical applications. The compact and ready-to-use Maxicaps® MR comes pre-sterilized and pre-assembled with 90 % less tubing and connectors compared to standard multi-capsule assemblies. Maxicaps® MR is the only logical choice for the lowest total cost of ownership in large-scale single-use processes.



## Introduction

Single-use filter capsules have been systematically replacing stainless steel housings and filter cartridges as a highly economical and risk-adverse choice for the biopharmaceutical industry. From capsules to complex custom assemblies, implementation of single-use filter systems reduces the time it takes for equipment setup and virtually eliminates the need for cleaning.

Conventional multi-round filter housings have now evolved into single-use Maxicaps® MR systems to meet today's advanced requirements. Until Maxicaps® MR, there has been no single-use equivalent to large-scale, multi-round filter configurations provided by stainless steel systems. Maxicaps® MR is the first ready-to-use, fully self-contained, single-use filtration unit featuring a wide choice of configurations. With 90% less tubing and only two connections, Maxicaps® MR reduces the installation time and the risk of operating errors significant.

## Single-Use Applications

- Media & feeds filtration
- Post cell harvest bioburden reduction for mAb's
- Clarification of vaccines
- Capture-column guard filtration
- Large-scale buffer preparation
- Virus filtration upstream and downstream
- Adsorptive virus pre-filtration

## Features

- Filtration area of up to 27 m<sup>2</sup>
- Complete device integrity testable as a single unit
- Large variety of pre-, sterile- and virus filters
- Flexible connections: Opta®, 1.5" Tri-Clamp, AseptiQuik®\* for virus filters or weldable tubing
- One single air filter for easy system venting

## Benefits

- Ready-to-use → pre-sterilized & pre-assembled
- Certified safety → Sterile & Sanitary delivery option
- Risk mitigation → 90% less tubings & connectors
- Space saving → Compact and organized design
- Time saving → 90% less test time – saves up to 4 hours

## Delivery Conditions

### Sterile

- For all gamma stable filter materials
- Assembled in a classified clean room, complete device gamma irradiated in a validated sterilization procedure

### Sanitary

- For non-gamma stable filter materials
- All fluid contact materials are sterilized in validated sterilization procedures and assembled in a classified clean room following specific hygienic measures and rules of conduct

### Non-Sterile

- For non-gamma stable filter materials
- Assembled in a classified clean room

## Validation

Maxicaps® MR have been qualified applying the most comprehensive and innovative test regimes. Biological, chemical and physical tests combined with extensive extractable testing. A sterilization validation in order to obtain a 10<sup>-6</sup> Sterility Assurance Level was performed to demonstrate the effectiveness of the gamma sterilization method for configurations with gamma stable filter material. The Maxicaps® filter capsules of the Sanitary delivery option are sterilized by autoclaving using a validated process following DIN|EN ISO 17665-1 regulations.

\*AseptiQuik® is a registered trademark of the Colder Products Company.

# Technical Data

Filter Families	Filtration Area	Materials	Max. Diffusion	Delivery Condition
<b>Sartopore® Platinum</b>				
MR3	9 m <sup>2</sup>   96.9 ft <sup>2</sup>	Polyethersulfone, surface modified	225 ml/min at 2.5 bar   36 psi	Gamma Irradiated
MR6	18 m <sup>2</sup>   193.8 ft <sup>2</sup>	Polyethersulfone, surface modified	450 ml/min at 2.5 bar   36 psi	Gamma Irradiated
MR9	27 m <sup>2</sup>   290.7 ft <sup>2</sup>	Polyethersulfone, surface modified	675 ml/min at 2.5 bar   36 psi	Gamma Irradiated
<b>Sartopore® 2 0.2 µm</b>				
MR3	5.4 m <sup>2</sup>   58.2 ft <sup>2</sup>	Polyethersulfone	162 ml/min at 2.5 bar   36 psi	Gamma Irradiated
MR6	10.8 m <sup>2</sup>   116.4 ft <sup>2</sup>	Polyethersulfone	324 ml/min at 2.5 bar   36 psi	Gamma Irradiated
MR9	16.2 m <sup>2</sup>   174.6 ft <sup>2</sup>	Polyethersulfone	486 ml/min at 2.5 bar   36 psi	Gamma Irradiated
<b>Sartopore® 2 XLG</b>				
MR3	7.2 m <sup>2</sup>   77.4 ft <sup>2</sup>	Polyethersulfone	207 ml/min at 2.5 bar   36 psi	Gamma Irradiated
MR6	14.4 m <sup>2</sup>   154.8 ft <sup>2</sup>	Polyethersulfone	414 ml/min at 2.5 bar   36 psi	Gamma Irradiated
MR9	21.6 m <sup>2</sup>   232.2 ft <sup>2</sup>	Polyethersulfone	621 ml/min at 2.5 bar   36 psi	Gamma Irradiated
<b>Sartopore® 2 XLI</b>				
MR3	7.2 m <sup>2</sup>   77.4 ft <sup>2</sup>	Polyethersulfone	189 ml/min at 2.5 bar   36 psi	Gamma Irradiated
MR6	14.4 m <sup>2</sup>   154.8 ft <sup>2</sup>	Polyethersulfone	378 ml/min at 2.5 bar   36 psi	Gamma Irradiated
MR9	21.6 m <sup>2</sup>   232.2 ft <sup>2</sup>	Polyethersulfone	567 ml/min at 2.5 bar   36 psi	Gamma Irradiated
<b>Sartopore® 2 XLM</b>				
MR3	7.2 m <sup>2</sup>   77.4 ft <sup>2</sup>	Polyethersulfone	180 ml/min at 2.5 bar   36 psi	Gamma Irradiated
MR6	14.4 m <sup>2</sup>   154.8 ft <sup>2</sup>	Polyethersulfone	360 ml/min at 2.5 bar   36 psi	Gamma Irradiated
MR9	21.6 m <sup>2</sup>   232.2 ft <sup>2</sup>	Polyethersulfone	540 ml/min at 2.5 bar   36 psi	Gamma Irradiated
<b>Sartoguard PES 0.1 µm nom.</b>				
MR3	7.2 m <sup>2</sup>   77.4 ft <sup>2</sup>	Polyethersulfone	225 ml/min at 1.5 bar   22 psi	Gamma Irradiated
MR6	14.4 m <sup>2</sup>   154.8 ft <sup>2</sup>	Polyethersulfone	450 ml/min at 1.5 bar   22 psi	Gamma Irradiated
MR9	21.6 m <sup>2</sup>   232.2 ft <sup>2</sup>	Polyethersulfone	675 ml/min at 1.5 bar   22 psi	Gamma Irradiated
<b>Sartoguard PES 0.2 µm nom.</b>				
MR3	7.2 m <sup>2</sup>   77.4 ft <sup>2</sup>	Polyethersulfone	162 ml/min at 1.2 bar   17.5 psi	Gamma Irradiated
MR6	14.4 m <sup>2</sup>   154.8 ft <sup>2</sup>	Polyethersulfone	324 ml/min at 1.2 bar   17.5 psi	Gamma Irradiated
MR9	21.6 m <sup>2</sup>   232.2 ft <sup>2</sup>	Polyethersulfone	486 ml/min at 1.2 bar   17.5 psi	Gamma Irradiated
<b>Sartopure® GF Plus 0.65 &amp; 1.2 µm nom.</b>				
MR3	3.6 m <sup>2</sup>   38.7 ft <sup>2</sup>	Glass Fiber		Sanitary or Non-Sterile
MR6	7.2 m <sup>2</sup>   77.4 ft <sup>2</sup>	Glass Fiber		Sanitary or Non-Sterile
MR9	10.8 m <sup>2</sup>   116.1 ft <sup>2</sup>	Glass Fiber		Sanitary or Non-Sterile

Filter Families	Filtration Area	Materials	Max. Diffusion	Delivery Condition
<b>Sartopure® PP3 0.45 µm nom.</b>				
MR3	3.6 m²   38.7 ft²	Polypropylene		Sanitary or Non-Sterile
MR6	7.2 m²   77.4 ft²	Polypropylene		Sanitary or Non-Sterile
MR9	10.8 m²   116.1 ft²	Polypropylene		Sanitary or Non-Sterile
<b>Sartopure® PP3 0.65, 1.2 &amp; 3.0 µm nom.</b>				
MR3	4.05 m²   43.5 ft²	Polypropylene		Sanitary or Non-Sterile
MR6	8.1 m²   87 ft²	Polypropylene		Sanitary or Non-Sterile
MR9	12.15 m²   130.5 ft²	Polypropylene		Sanitary or Non-Sterile
<b>Sartopure® PP3 5.0, 8.0, 20.0 &amp; 50.0 µm nom.</b>				
MR3	5.85 m²   63 ft²	Polypropylene		Sanitary or Non-Sterile
MR6	11.7 m²   126 ft²	Polypropylene		Sanitary or Non-Sterile
MR9	17.55 m²   189 ft²	Polypropylene		Sanitary or Non-Sterile
<b>Virosart® HF 20 nm nominal hollow fibre</b>				
MR2	4.8 m²   51.7 ft²	Polyethersulfone surface modified	≤ 41 ml/min at 2.5 bar   36 psi	Gamma Irradiation
MR3	7.2 m²   77.5 ft²	Polyethersulfone surface modified	≤ 60 ml/min at 2.5 bar   36 psi	Gamma Irradiation
MR4	9.6 m²   103.3 ft²	Polyethersulfone surface modified	≤ 79 ml/min at 2.5 bar   36 psi	Gamma Irradiation
MR5	12 m²   129.2 ft²	Polyethersulfone surface modified	≤ 99 ml/min at 2.5 bar   36 psi	Gamma Irradiation
MR6	14.4 m²   155 ft²	Polyethersulfone surface modified	≤ 117 ml/min at 2.5 bar   36 psi	Gamma Irradiation
<b>Virosart® Media 20 nm nominal hollow fibre</b>				
MR3	3 m²   32.3 ft²	Polyethersulfone surface modified	≤ 48 ml/min at 2.5 bar   36 psi	Gamma Irradiation
MR6	6 m²   64.6 ft²	Polyethersulfone surface modified	≤ 97 ml/min at 2.5 bar   36 psi	Gamma Irradiation
<b>Virosart® Max 0.1 µm</b>				
MR3	6.3 m²   68 ft²	Polyamide	≤ 16 ml/min at 2.0 bar   29 psi	Sanitary or Non-Sterile
MR6	12.6 m²   136 ft²	Polyamide	≤ 31 ml/min at 2.0 bar   29 psi	Sanitary or Non-Sterile
MR9	18.9 m²   203 ft²	Polyamide	≤ 46 ml/min at 2.0 bar   29 psi	Sanitary or Non-Sterile

## Max.Differential Pressure

2.5 bar | 36 psi at 20°C

## Accessoires (Reusable – Need to Be Ordered Separately)

SU Valve Actuator\*      Order Code: BPR0202  
Pressure Safety Device      Order Code: 26787---PS

\*3 reusable actuators are needed for each Maxicaps® MR

## Delivery Condition

Sterile, for gamma stable filter material  
Sanitary, for non-gamma stable filter material  
Non-Sterile, for non-gamma stable filter material

## Materials

### Filter Material

Refer to the technical reference of the respective filter.

### Maxicaps® Housing and Distribution Manifold Pipes

Polypropylene (PP)

### Inlet | Outlet Tubing

Silicone (reinforced)  
Thermoplastic Elastomer (TPE)

### Rack

Polypropylene (PP), Polyethylene (PE)

### Mounting Parts

Screws, Washer, Threaded Rod: Stainless Steel  
Gaskets: Silicone  
Tri Clamp: Polyamide (PA)

### Venting

Sartopore® Air with hydrophobic Polyethersulfone (PES)  
Pure-Fit TCL Clamp: Polyvinylidenfluorid (PVDF)  
Inspection Glass: Polyethylenterephthalat (PET)

## Technical References

🌐 For further information regarding pre-, sterile- & virus filters please click [here](#).

For further information on Maxicaps® MR & Virosart® Validation Guides, please see references below:

### Maxicaps® MR

Validation Guide Maxicaps® MR      2646224

### Virosart® HF

Datasheet      SPK2180-e  
Validation Guide      SPK5801-e

### Virosart® Media

Datasheet      DIR 2650737  
Validation Guide      SPK5812-e

### Virosart® Max

Datasheet      DIR 2650739  
Validation Guide      DIR 2650008

## Regulatory Compliance

- Each individual Maxicaps® element is tested for integrity (membrane filters only).
- Fully validated as sterilizing grade filters according to current ASTM F838 guideline for Sartopore® filter family.
- Designed, developed and manufactured in accordance with ISO 9001 certified Quality Management System.
- Non pyrogenic according to USP Bacterial Endotoxins.
- All assembled filters and tubing meet the requirements of the current USP Class VI Biological reactivity tests.
- Non-fiber releasing: This product is manufactured with membranes which meet the criteria for a “non-fiber releasing” filter as defined in 21 CFR 210.3 (b) (6) and 211.72.
- This product is conform to Pressure Equipment Directive 2014/68/EU.

# Ordering Information

## Sartopore® Platinum

549	73	07H	3	G	-	
-----	----	-----	---	---	---	--

**Pore Size**  
07H: 0.45 µm + 0.2 µm


**Number of Filter Elements per Device**  
MR3: 3 Filter Elements  
MR6: 6 Filter Elements  
MR9: 9 Filter Elements

## Sartopore® 2

544	73		3	G	-	
-----	----	--	---	---	---	--

**Pore Size**  
07H: 0.45 µm + 0.2 µm  
07G: 0.8 µm + 0.2 µm (XLG)  
07I: 0.35 µm + 0.2 µm (XLI)  
58M: 0.2 µm + 0.1 µm (XLM)

**Number of Filter Elements per Device**  
MR3: 3 Filter Elements  
MR6: 6 Filter Elements  
MR9: 9 Filter Elements



## Sartoguard PES

547	73		3	G	-	
-----	----	--	---	---	---	--

**Pore Size**  
07F: 0.2 µm nominally  
58G: 0.1 µm nominally

**Number of Filter Elements per Device**  
MR3: 3 Filter Elements  
MR6: 6 Filter Elements  
MR9: 9 Filter Elements

## Sartopure® PP3

505	73		3	--		
-----	----	--	---	----	--	--

**Retention Ratings**  
50P: 50 µm  
20P: 20 µm  
01P: 8 µm  
42P: 5 µm  
02P: 3 µm  
03P: 1.2 µm  
05P: 0.65 µm -- Non-Sterile  
06P: 0.45 µm C- Sanitary

**Number of Filter Elements per Device**  
MR3: 3 Filter Elements  
MR6: 6 Filter Elements  
MR9: 9 Filter Elements



## Sartopure® GF Plus

555	73		3	--		
-----	----	--	---	----	--	--

**Retention Ratings**  
03P: 1.2 µm  
05P: 0.65 µm  
-- Non-Sterile  
C- Sanitary

**Number of Filter Elements per Device**  
MR3: 3 Filter Elements  
MR6: 6 Filter Elements  
MR9: 9 Filter Elements

## Virosart® HF

3VI--	28-	M	C	G-	MR2
-------	-----	---	---	----	-----

Number of Filter Elements  
per Device

MR2: 2 Filter Elements

MR3: 3 Filter Elements

MR4: 4 Filter Elements

MR5: 5 Filter Elements

MR6: 6 Filter Elements



MR6

## Virosart® Media

3V2--	28-	I	V	G-	MR3
-------	-----	---	---	----	-----

Number of Filter Elements  
per Device

MR3: 3 Filter Elements

MR6: 6 Filter Elements



MR9

## Virosart® Max

54A	73	58	N3	--	MR3
-----	----	----	----	----	-----

-- Non-Sterile  
C- Sanitary

Number of Filter Elements  
per Device

MR3: 3 Filter Elements

MR6: 6 Filter Elements

MR9: 9 Filter Elements

### Germany

Sartorius Stedim Biotech GmbH  
August-Spindler-Strasse 11  
37079 Goettingen  
Phone +49 551 308 0

### USA

Sartorius Stedim North America Inc.  
565 Johnson Avenue  
Bohemia, NY 11716  
Toll-Free +1 800 368 7178

For further contacts, visit  
[www.sartorius.com](http://www.sartorius.com)

Specifications subject to change without notice.

© 2021 Sartorius Stedim Biotech GmbH, August-Spindler-Strasse 11, 37079 Goettingen, Germany

DIR No.: 2649560-000-01 | Order No.: 85037-559-85 | Status: 03 | 03 | 2021