4Cell® Nutri-T Media

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
Figure 1: (A) 0.5M healthy donor PBMCs were seeded in 24 w/p (2 mL media/well). Cells were activated with TransAct 1:100 and 600 IU/mL IL-2. Cells were split and media renewed every 2-3 days. Fold expansion (FE) and cell viability were measured at Day 11.

(B) 0.2M PBMCs from healthy donors were seeded in 24 w/p (2 mL media/well). Cells were activated with TransAct 1:100 and 600 IU/mL IL-2. 24 hr. After seeding cells were transduced with a lentiviral vector expressing an EGFR-CAR-T. Cells were split and media renewed every 2-3 days. FE and cell viability were measured at Day 11.

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
Figure 2: (A) TILs were isolated from a melanoma patient. 15,000 cells were seeded in a 24 w/p (2 mL/well) with PBMCs (1:100). Cells were activated with IL-2 (3,000 IU/mL) and OKT-3 (50 ng/mL). 2 mL and 4 mL of fresh medium + IL2 were added at Days 5 and 7 respectively (total volume of 8 mL). Fold expansion was measured at 14 days. Inherent variations among primary T lymphocyte donor populations may result in varying outcomes.

(B) PBMCs were separated from peripheral blood of a lymphoma patient. Tested mediums were supplemented with 50 ng/mL OKT3 and 300 IU/mL IL2. At Day 2 post seeding, 2-3M cells for the G-Rex24 were transduced with a CD19-CAR lentiviral vector in 6w/p pre-coated with RTN. Post transduction the cells were collected and reseeded. At day 4, 4 mL fresh medium +IL2 were added and at Day 6, 50% medium was replaced with fresh medium + IL2. At Day 9 transduction efficiency was evaluated and at Day 10 Fold expansion was measured.

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
🌐 For further contacts, visit
www.sartorius.com

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Sartorius Stedim Biotech GmbH
Sartorius Stedim North America Inc.
Biological Industries Israel Beit Haemek Ltd.