

# Steamer Line Hollow Fiber Membrane Products

Hollow Fiber Modules  
Designed for Perfusion and  
Aseptic Applications.



## Executive Summary

Steamer Line hollow fiber membrane products incorporate our latest glycerin free low extractables, heat resistant modified polyethersulfone (m-PES) membrane technology. All Steamer Line modules are gamma irradiated and ready to use without any tedious pre-rinse. The extractables level for the Steamer Line modules is approximately 80x less than a glycerin conditioned membrane, see Extractables Flush Analysis. After a quick buffer conditioning the module is ready to be used or autoclaved. For autoclaving directions, please see the Steamer Line Validation Guide. Steamer Line modules are fully scalable from batch volumes from 10 mL up to 1500 L with corresponding membrane surface areas from 0.056 ft<sup>2</sup> (0.0052 m<sup>2</sup>) up to 166.0 ft<sup>2</sup> (15.42 m<sup>2</sup>).

## Unique Selling Points

The Steamer Line modules offer true ready to use crossflow filtration


- Self-contained units
- No humectants - 80x less extractables
- No pre-rinse required (glycerin free)
- Ready to use
- Easy and reliable scale up and scale down
- Membrane offering from 3kD to 750kD MWCO and from 0.1 µm to 0.65 µm pore size
- Gamma irradiated
- Cost efficiency
- Free of time-consuming cleaning validation
- Integrity tested and individual lot number for easy traceability

## Relevant Applications

- Rapid clarification of volumes from 25 L to 100 L of cell culture, fermentation solutions and virus | vaccine suspensions
- Production scale aseptic operations or any other sterile application requiring autoclaving
- Cell perfusion with a conventional pump operated cross-flow perfusion system. Bioreactor and Steamer hollow fiber assembly can be autoclaved simultaneously
- Concentration and diafiltration of gene therapy products as well as monoclonal antibodies, enzymes, blood components and other proteins

## Relevant Process Steps

### Extractables Flush Analysis - Glycerin Free vs Glycerin Stored



**Membrane**  
m-PES 30kD, 1mm Fiber ID

**Extractant**  
Water (20 °C)

■ GF Gamma Irradiated (Glycerin Free)  
Autoclaveable Steamer Line

■ Glycerin Stored (No Gamma)

## Technical Specifications

Attribute	Steamer Line Modules
Bio compatibility	The module as assembled is USP Class VI compliant
Membrane	Modified Polyethersulfone (m-PES)
Housing	White Polysulfone
Encapsulant	Medical Grade Epoxy

Name	Scale	Fiber Length	Dimensions in inches (cm)	Membrane surface area ft <sup>2</sup> (m <sup>2</sup> )*	Recommended batch volume per module	Connections	
						Feed   Retentate	Feed   Retentate
Discover	Lab	12 inch	0.38 × 12.88 (0.95 × 32.72)	0.056 (0.0052)	10 - 250 mL	Luer Lock	Luer Lock
		24 inch	0.38 × 24.38 (0.95 × 61.93)	0.115 (0.0107)	50 - 400 mL		
		41 inch	0.38 × 41.8 (0.95 × 106.2)	0.208 (0.0193)	80 - 850 mL		
Explorer	Lab	12 inch	0.5 × 12.3 (1.3 × 31.2)	0.17 (0.0155)	150 - 175 mL	½-inch TC	⅜-inch Hose Barb
		24 inch	0.5 × 23.8 (1.3 × 60.5)	0.35 (0.0321)	250 - 1,500 mL		
		41 inch	0.5 × 41.8 (1.3 × 106.2)	0.62 (0.0579)	300 - 3,000 mL		
Researcher	Lab   Pilot	12 inch	0.75 × 12.3 (1.91 × 31.2)	0.48 (0.0444)	400 - 2,000 mL	¾-inch TC	⅜-inch Hose Barb
		24 inch	0.75 × 23.66 (1.91 × 60.1)	1.01 (0.0940)	700 - 4,000 mL		
		41 inch	0.75 × 41.8 (1.91 × 106.1)	1.85 (0.1716)	1,000 - 8,000 mL		
Researcher XL	Lab   Pilot	12 inch	0.75 × 12.3 (1.91 × 31.2)	0.66 (0.0617)	0.56 - 2.8 L	¾-inch TC	⅜-inch Hose Barb
		24 inch	0.75 × 23.66 (1.91 × 60.1)	1.41 (0.1305)	0.98 - 5.6 L		
		41 inch	0.75 × 41.8 (1.91 × 106.1)	2.57 (0.2383)	1.4 - 11.2 L		
Investigator	Lab   Pilot	12 inch	1.32 × 12.01 (3.34 × 30.5)	1.4 (0.13)	1 - 6 L	1.5-inch TC	½-inch TC
		24 inch	1.32 × 23.5 (3.34 × 59.7)	3.0 (0.28)	2 - 12 L		
		41 inch	1.32 × 41.5 (3.34 × 105.4)	5.4 (0.51)	3 - 25 L		
mini- Bioproducer	Pilot   Production	12 inch	2.70 × 15.0 (6.86 × 38.1)	6.6 (0.61)	5 - 50 L	1.5-inch TC	0.75-inch TC
		24 inch	2.70 × 26.5 (6.86 × 67.3)	14.5 (1.35)	10 - 100 L		
		41 inch	2.70 × 44.5 (6.86 × 113)	26.9 (2.50)	25 - 250 L		
Bioproducer	Production	12 inch	3.50 × 15.0 (8.89 × 38.1)	13.1 (1.21)	25 - 100 L	1.5-inch TC	1-inch TC
		24 inch	3.50 × 28.5 (8.89 × 67.3)	29.0 (2.69)	50 - 250 L		
		41 inch	3.50 × 44.49 (8.89 × 113)	53.8 (5.00)	100 - 500 L		
Maximizer	Production	24 inch	4.62 × 28.8 (11.73 × 73.15)	58.4 (5.42)	100 - 500 L	2-inch TC	1-inch TC
		41 inch	4.62 × 46.5 (11.73 × 118.1)	104.9 (9.74)	250 - 1,000 L		
Grand XL (only with 1.0 and 2.0 mm Fiber IDs)	Production	24 inch	6.64 × 32.65 (16.88 × 82.93)	87.2 (8.10)	250 - 800 L	2.5-inch TC or 6-inch TC	1-inch TC
		43 inch	6.015 × 51.65 (16.88 × 131.19)	166.0 (15.42)	400 - 1,500 L		

\*Membrane surface valid for all fiber IDs. For number of fibers per module please consider the Validation Guide.

# Available MWCO and Pore Sizes for Different Fiber IDs


The tables stated below are valid for all sizes and lengths, if not stated otherwise.

Fiber ID	MWCO								
	3 kD	5 kD	10 kD	30 kD	50 kD	100 kD	300 kD	500 kD	750 kD
0.5 mm*	■	■	■	■	■	■	■	■	■
1.0 mm	■	■	■	■	■	■	■	■	■
2.0 mm			■		■		■	■	■

Fiber ID	Pore size			
	0.1 µm	0.2 µm	0.45 µm	0.65 µm
0.5 mm*	■	■	■	■
1.0 mm	■	■	■	■
2.0 mm	■	■	■	■

\* Not available for Grand XL

## Ordering Information




Please consult the specification table on page 3 to learn more about the possible combinations.

# Overview


## Discover

Exact length for different Fiber length: 12, 24 or 41 inch




## Explorer

Exact length for different Fiber length: 12, 24 or 41 inch




## Researcher & Researcher XL

Exact length for different Fiber length: 12, 24 or 41 inch




## Investigator

Exact length for different Fiber length: 12, 24 or 41 inch




## mini-Bioproducer

Exact length for different Fiber length: 12, 24 or 41 inch




## Bioproducer

Exact length for different Fiber length: 12, 24 or 41 inch




## Maximizer

Exact length for different Fiber length: 24 or 41 inch



## Grand XL

Exact length for different Fiber length: 24 or 43 inch





**Germany**

Sartorius Stedim Biotech GmbH  
August-Spindler-Strasse 11  
37079 Goettingen  
Phone +49 551 308 0

**USA**

Sartorius Stedim North America Inc.  
565 Johnson Avenue  
Bohemia, NY 11716  
Toll-Free +1 800 368 7178

 For further contacts, visit  
[www.sartorius.com](http://www.sartorius.com)