

# Reuse Line Hollow Fiber Modules

Low Binding m-PES  
Membrane Enables Reliable  
Crossflow, High Product Flux  
and Easy Cleaning.

## Executive Summary

Reuse Line hollow fiber modules offers a modified polyethersulfone (m-PES) membrane which is gentle on your cells, biomolecules and viruses. Generate high yields and low hold-up volumes. Reuse Line modules provide a linear and predictive scale-up process from laboratory to pilot-scale to manufacturing scale by using matching materials, fluid-path length, and performance characteristics. Reuse Line modules are fully scalable from batch volumes from 10 mL up to 1500 L with corresponding membrane surface areas from 0.056 ft<sup>2</sup> (0.0052 m<sup>2</sup>) up to 166.0 ft<sup>2</sup> (15.42 m<sup>2</sup>). Due to the inhouse production of the membrane the Reuse Line can offer you a high batch-to-batch consistency. Reuse Line hollow fiber modules can be sanitized and cleaned in 0.5 - 1.0 N NaOH, and stored in 0.1 N NaOH between uses.

## Unique Selling Points

For those applications where cleaning and Reuse is warranted, the Reuse Line hollow fiber modules offer:

- Self-contained
- Membrane offering from 3 kD to 750 kD MWCO and from 0.1 µm to 0.65 µm pore size
- Easy and reliable scale up and scale down
- Low hold up volume.
- High product flux and total capacity
- Robust, strong, multi-use hollow fiber membranes
- Low binding m-PES membrane - high yield and easy to clean



## Relevant Applications

- Concentration and purification of vaccines
- Concentration and diafiltration of gene therapy products
- Cell-harvest (e.g. excellent results have been achieved with both *E. Coli* whole cells and *E. Coli* lysates, as well as other microbial process streams.)
- Clarification of mammalian | CHO cell cultures and maximizing protein recovery
- Concentration | diafiltration of monoclonal antibodies, recombinant proteins, biological macromolecules and peptides.

## Technical Specifications

Attribute	Reuse Line Modules
Bio compatibility	The module as assembled is USP Class VI compliant
Membrane	Modified Polyethersulfone (m-PES)
Housing	White Polysulfone
Encapsulant	Medical Grade Epoxy

Name	Scale	Fiber Length	Dimensions in inches (cm)	Membrane surface area ft <sup>2</sup> (m <sup>2</sup> )*	Recommended batch volume per module	Connections	
						Feed   Retentate	Feed   Retentate
Discover	Lab	12 inch	0.38 × 12.88 (0.95 × 32.72)	0.056 (0.0052)	10 - 250 mL	Luer Lock	Luer Lock
		24 inch	0.38 × 24.38 (0.95 × 61.93)	0.115 (0.0107)	50 - 400 mL		
		41 inch	0.38 × 41.8 (0.95 × 106.2)	0.208 (0.0193)	80 - 850 mL		
Explorer	Lab	12 inch	0.5 × 12.3 (1.3 × 31.2)	0.17 (0.0155)	150 - 175 mL	½-inch TC	⅜-inch Hose Barb
		24 inch	0.5 × 23.8 (1.3 × 60.5)	0.35 (0.0321)	250 - 1,500 mL		
		41 inch	0.5 × 41.8 (1.3 × 106.2)	0.62 (0.0579)	300 - 3,000 mL		
Researcher	Lab   Pilot	12 inch	0.75 × 12.3 (1.91 × 31.2)	0.48 (0.0444)	400 - 2,000 mL	¾-inch TC	⅜-inch Hose Barb
		24 inch	0.75 × 23.66 (1.91 × 60.1)	1.01 (0.0940)	700 - 4,000 mL		
		41 inch	0.75 × 41.8 (1.91 × 106.1)	1.85 (0.1716)	1,000 - 8,000 mL		

<b>Researcher XL</b>	Lab   Pilot	12 inch	0.75 × 12.3 (1.91 × 31.2)	0.66 (0.0617)	0.56 - 2.8 L	¾-inch TC	⅜-inch Hose Barb
		24 inch	0.75 × 23.66 (1.91 × 60.1)	1.41 (0.1305)	0.98 - 5.6 L		
		41 inch	0.75 × 41.8 (1.91 × 106.1)	2.57 (0.2383)	1.4 - 11.2 L		
<b>Investigator</b>	Lab   Pilot	12 inch	1.32 × 12.01 (3.34 × 30.5)	1.4 (0.13)	1 - 6 L	1.5-inch TC	½-inch TC
		24 inch	1.32 × 23.5 (3.34 × 59.7)	3.0 (0.28)	2 - 12 L		
		41 inch	1.32 × 41.5 (3.34 × 105.4)	5.4 (0.51)	3 - 25 L		
<b>mini-Bioproducer</b>	Pilot   Production	12 inch	2.70 × 15.0 (6.86 × 38.1)	6.6 (0.61)	5 - 50 L	1.5-inch TC	0.75-inch TC
		24 inch	2.70 × 26.5 (6.86 × 67.3)	14.5 (1.35)	10 - 100 L		
		41 inch	2.70 × 44.5 (6.86 × 113)	26.9 (2.50)	25 - 250 L		
<b>Bioproducer</b>	Production	12 inch	3.50 × 15.0 (8.89 × 38.1)	13.1 (1.21)	25 - 100 L	1.5-inch TC	1-inch TC
		24 inch	3.50 × 28.5 (8.89 × 67.3)	29.0 (2.69)	50 - 250 L		
		41 inch	3.50 × 44.49 (8.89 × 113)	53.8 (5.00)	100 - 500 L		
<b>Maximizer</b>	Production	24 inch	4.62 × 28.8 (11.73 × 73.15)	58.4 (5.42)	100 - 500 L	2-inch TC	1-inch TC
		41 inch	4.62 × 46.5 (11.73 × 118.1)	104.9 (9.74)	250 - 1,000 L		
<b>Grand XL (only with 1.0 and 2.0 mm Fiber IDs)</b>	Production	24 inch	6.64 × 32.65 (16.88 × 82.93)	87.2 (8.10)	250 - 800 L	2.5-inch TC or 6-inch TC	1-inch TC
		43 inch	6.015 × 51.65 (16.88 × 131.19)	166.0 (15.42)	400 - 1,500 L		

\*Membrane surface valid for all fiber IDs. For number of fibers per module please consider the Validation Guide.

## Available MWCO and Pore Sizes for Different Fiber IDs


The tables stated below are valid for all sizes and lengths, if not stated otherwise.

MWCO									
Fiber ID	3 kD	5 kD	10 kD	30 kD	50 kD	100 kD	300 kD	500 kD	750 kD
0.5 mm**	■	■	■	■	■	■	■	■	■
1.0 mm	■	■	■	■	■	■	■	■	■
2.0 mm			■		■		■	■	■

Pore size				
Fiber ID	0.1 µm	0.2 µm	0.45 µm	0.65 µm
0.5 mm**	■	■	■	■
1.0 mm	■	■	■	■
2.0 mm	■	■	■	■

\*\* Not available for Grand XL

# Ordering Information




Please consult the specification table on page 3 to learn more about the possible combinations.

# Overview


## Discover

Exact length for different Fiber length: 12, 24 or 41 inch




## Explorer

Exact length for different Fiber length: 12, 24 or 41 inch




## Researcher & Researcher XL

Exact length for different Fiber length: 12, 24 or 41 inch




## Investigator

Exact length for different Fiber length: 12, 24 or 41 inch




## mini-Bioproducer

Exact length for different Fiber length: 12, 24 or 41 inch




## Bioproducer

Exact length for different Fiber length: 12, 24 or 41 inch




## Maximizer

Exact length for different Fiber length: 24 or 41 inch



## Grand XL

Exact length for different Fiber length: 24 or 43 inch





**Germany**

Sartorius Stedim Biotech GmbH  
August-Spindler-Strasse 11  
37079 Goettingen  
Phone +49 551 308 0

**USA**

Sartorius Stedim North America Inc.  
565 Johnson Avenue  
Bohemia, NY 11716  
Toll-Free +1 800 368 7178

 For further contacts, visit  
[www.sartorius.com](http://www.sartorius.com)