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In addition to sales of measurement devices, FUKUDA also offers advice on test conditions integrity and improvement in order for customers to test in a more accurate and safer environment.

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Rev. Feb.19 Printed Feb.19 1KJ Printed in Japan 10661-K-001-08

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Air Leak Test Systems For Packaging Containers

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The integrity testing of pharmaceuticals

Quantitative test methods that do not require wetting the product with water

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- Air leak testing device for sample inspections
- Barrier property checking (Oxygen, Steam, Microorganisms, Bacteria, etc.)
- Testing is also quantified for scientific quality management
- Test Results that are both quantitative and reproducible
- Equivalent to a minimum opening diameter of φ5μm
- GMP, DI, and Validation Compatible

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包装容器の密封試験を自動化、数値化、乾式化
The transition to automated, quantified, and dry packaging seal tests

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Sealed packaging is being used more and more in order to prevent the deterioration, degeneration etc. of medical products, cosmetics. As such, if there are any defects that impair the sealing properties, such as Pin-hole or seal abnormalities, these may cause serious problems that lead to product rejects. Air leak testing achieves the transition from previous water immersion testing to a quantitative method that also delivers reproducibility.

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Scientific validation is realized by air leak testing

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Differential Pressure Type Air Leak Testing

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For air leak testing of packaging containers, leakage tests of the object to be tested (hereinafter referred to test pieces) are performed using pressure variations from the test piece. As shown in the figure below, the test is performed by two steps of Large Leak Detection and Small Leak Detection.

- Large Leak Detection: Checking for sealing defects such as coming away, peeling, and tearing of the packaging.
- Small Leak Detection: Checking for small leak less than φ 20 μm in diameter.

Prior testing of the shape, volume, and test pressure, etc. of your products is required in order to perform optimal leak testing. Please do not hesitate to contact us to discuss your specific needs.

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Self Check Function

This function checks for malfunctions in various valves or
sensor faults.

Type Registration / Product Setting Function

Product type and setting conditions can be registered in 1-100
(100 types) groups.

Data Integrity, Data Falsification Preventing
(USB, Printer)

Log data such as test results, pressure values, and setting
values is encrypted, and can be recorded or converted to PDF
format on a PC.

ID Control

It is possible to identify the operator and the administrator
separately using ID's and passwords.

Barcode Function

You can select a type, and input an ID at login.

Sensitivity Check

Sensor sensitivity can be confirmed by built-in flow reference
units (FFM series).

Simple Setup

Automatic settings of testing conditions using good / defective
products.

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
Customized Ordering (Tabletop units / Total inspection units)

Please contact us through the website or directly through pharmaceutical
packaging specialist sales office for special order made or total inspection
units options.

E.g. For sample sizes that do not fit into the advertised model here.

➡ We will check the size of the product to be tested and design a model that
can handle it.

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Number		Item	Description
-	W	Measurement Method	Tank Pressurization Type/ Tank Partial Pressure Type
-	W	Pressure Range	Compound Pressure (Equipped with APU)
-	N	Additional Tank	None
①	V	Validation	With Validation
	N		Without Validation

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Item	Description
Measurement Method	Tank Pressurization Type/ Tank Partial Pressure Type
Test Pressure	1 ～ 200 kPa abs
Products Setting	1 ～ 100 Group（100 Types）
Operation Temperature/ Humidity Range	0 ～ 40℃ 45 ～ 85%RH（Without condensation）
Power Source/ Applied Current	100 ～ 240VAC 50/ 60Hz 100VA
External Input/ Output	Serial Communication 3Port (PC, Printer, Barcode reader) USB 1Port Ethernet 1Port
External Dimensions	MSP-0101 : W380 × D560 × H290 mm Dedicated Cradle : W380 × D560 × H700 mm
Mass	MSP-0101 25kg Dedicated Cradle 55kg
Maximum Test Parts Size	W130 × D180 × H65mm or below ※ Please consult FUKUDA if testing is required for sizes larger than that shown here.
Pneumatic Pressure Source	Pneumatic Pressure Source 0.4MPa or above. Please supply clean, stable, and dry air. Recommended conditions according to ISO 8573-1: 2010. Compressed air purity classes 1, 3, 1
Utility Dedicated Cradle (Sold separately)	Air leak testing is test method using compressed air. If the condition of dry air used is low, the test result will be affected. In order to ensure testing accuracy, we provide a dedicated cradle that meets the recommended pneumatic source conditions. The dedicated cradle is made up of a compressor, vacuum pump, and a dryer.
Optional (Sold separately)	Pin-hole (φ 5,10,15,20μm etc.), Pneumatic unit FR-57 with Bourdon tube, Vacuum pump, Barcode reader, USB memory, Printer (thermal/ impact dot matrix type), Printer roll paper

● Core(Capsule), Pin-hole, Utility (Vacuum pump, Compressor etc.) are not included in
the main unit price shown.

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