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
5991-5124JAJP

Agilent Technologies

mAb Titer Analysis with the Agilent Bio-Monolith Protein A Column

Authors

Emmie Dumont, Isabel Vandenheede, Pat Sandra, and Koen Sandra

Research Institute for Chromatography (RIC)

President Kennedypark 26
B-8500 Kortrijk, Belgium

James Martosella, Phu Duong
and Maureen Joseph
Agilent Technologies, Inc

Abstract

Monoclonal antibodies (mAbs) are becoming increasingly important in the treatment of various diseases. During development of recombinant mAbs, protein titer and yield from various cell culture supernatants must be monitored. This application note describes how the Agilent Bio-Monolith Protein A column was successfully applied in the determination of mAb concentrations.

Introduction

Protein A from *Staphylococcus aureus* has a very strong affinity for the Fc domain of immunoglobulins (IgG), allowing its capture from complex matrixes such as cell-culture supernatants. Affinity chromatography making use of Protein A is the gold standard in therapeutic monoclonal antibody (mAb) purification, and typically represents the first chromatographic step in downstream processing. Protein A chromatography finds applications beyond this large-scale purification. At the analytical scale, it is used early in the development of mAbs for the high-throughput determination of mAb titer and yield directly from cell culture supernatants, and to purify μg amounts of material for further measurements, for example, by mass spectrometry (MS), ion-exchange (IEX), size exclusion chromatography (SEC), or hydrophobic interaction chromatography (HIC).

This application note describes the use of the Agilent Bio-Monolith Protein A column in mAb titer analysis. This HPLC column (Figure 1) has a 5.2 mm id, a column length of 4.95 mm, and is composed of a highly cross-linked poly(glycidyl methacrylate-co-ethylene dimethacrylate) monolithic disk coated with native Protein A from *S. aureus*. Its monolithic nature, characterized by well-defined channels of 1,200 to 1,500 nm, and by the absence of pores and voids, delivers fast and efficient separations with negligible carryover and excellent robustness. These are features typically expected from a column for mAb titer analysis, to successfully guide clone selection and cell-culture optimization. We present the best practice for use of the column in the determination of absolute mAb concentrations in Chinese hamster ovary (CHO) cell-culture supernatants. Data from a trastuzumab biosimilar project are used for illustration purposes. Trastuzumab, marketed as Herceptin since 1998, is used in the treatment of HER2 positive breast cancer, and comes out of patent in 2014 in Europe, and 2018 in the United States.

Experimental

Materials

Water was obtained from Biosolve (Valkenswaard, The Netherlands). Citric acid, acetic acid, NaH_2PO_4 and Na_2HPO_4 were purchased from Sigma-Aldrich, Corp. (St. Louis, MO, USA). Humanized monoclonal antibody trastuzumab, marketed as Herceptin, was obtained from Roche (Basel, Switzerland). Trastuzumab biosimilar CHO cell culture supernatants were obtained from a local biotechnology company.

Sample preparation

Herceptin stock solution present at 21 mg/mL was diluted in mobile phase A prior to injection. Cell supernatants were diluted 1:1 in 50 mM Na_2HPO_4 . Supernatants were centrifuged at 5,000 g for 5 minutes prior to injection.

Instrumentation

Agilent Bio-Monolith Protein A measurements were performed on:

- Agilent 1100 Series Quaternary Pump (G1311A)
- Agilent 1100 Series Autosampler (G1313A)
- Agilent 1100 Series Diode Array Detector (G1315A)

Software

- Agilent Technologies OpenLAB CDS ChemStation revision C01.05 (35)

Conditions, Bio-Monolith column

Parameter	Value	
Column:	Agilent Bio-Monolith Protein A (p/n5069-3639)	
Mobile phase:	A) 50 mM phosphate, pH 7.4 B) 100 mM citric acid, pH 2.8 mM acetic acid, pH 2.6	
Gradient:	Time (min)	% B
	0 to 0.5	0 (binding)
	0.6 to 1.7	100 (elution)
	1.8 to 3.5	0 (regeneration)
Flow rate:	1 mL/min	
Injection volume:	Variable (50 μL , optimized for CHO cell culture supernatants)	
Detection:	UV at 280 nm	


Figure 1. UV 280 nm Agilent Bio-Monolith Protein A chromatogram of a trastuzumab-producing CHO clone, clone 9 (red), and of a Herceptin originator diluted in 50 mM Na-phosphate pH 7.4 to 0.2 mg/mL (blue). Note that the supernatant was diluted 1:1 in phosphate buffer.

Results and Discussion

Buffer selection

Figure 2 shows a typical chromatogram from the Protein A column. The example chromatogram is one injection of the supernatant of a specific trastuzumab-producing CHO clone. The unbound material eluted in the flow-through while the mAb was retained at neutral pH (binding) and was only released (elution) after lowering the pH upon applying a step gradient. In this case, 50 mM Na-phosphate at pH 7.4 was used for binding/loading, and 100 mM citric acid at pH 2.8 for elution. This represents a good starting condition for any application.

When developing a new method for a Protein A column, both binding and elution buffers should be optimized. For binding buffers, 50 mM Na phosphate, pH 7.4, is a good starting point, and can be optimized between pH 7 and 8. For elution buffers, the 100 mM citric acid used here is a good starting point. Other possible elution buffers are 500 mM acetic acid, pH 2.6, 100 mM glycine HCl, pH 2.8, and 12 mM HCl, pH 1.9.

Figure 3 compares the elution of a Herceptin originator with acetic acid and citric acid. Very similar peak shape and area were observed, although peaks were slightly sharper using citric acid. In the case of this Herceptin originator, no material was seen in the flow-through, which was not surprising since this represented the marketed product and was devoid of host-cell proteins. In the chromatograms shown, the flow rate was set at 1 mL/min. The monolithic nature of the support, characterized by convective instead of diffusive mass transfer, allowed for near flow-rate independence and, hence, high-throughput separations. This is highly desirable in mAb titer determination, which typically requires the processing of a wide range of samples. The maximum flow rate that can be applied on the column is 2 mL/min, which allows fast, sub-2-minute separations.


Figure 2. UV 280 nm Protein A chromatogram showing the supernatant of a trastuzumab-producing CHO clone. Injection volume was 50 μ L. Peak width at half height was 0.10 minutes for the unbound material and 0.06 minutes for the retained mAb.


Figure 3. UV 280 nm Protein A chromatogram of Herceptin originator diluted in 50 mM Na-phosphate, pH 7.4, to 0.5 mg/mL (50- μ L injection, 25- μ g column load). Elution was achieved using citric acid (A) and acetic acid (B). Peak width at half height is 0.057 and 0.067 minutes for citric acid and acetic acid, respectively.

Precision, linearity, carryover, and injection size

Precision is critically important in the determination of the mAb titer. Table 1 shows the peak area and retention time repeatability that can typically be expected upon injecting a Herceptin originator 10 times. Chromatograms are shown in Figure 4. More than acceptable relative standard deviation (RSD) values were obtained for both citric acid and acetic acid as elution buffers. Carryover was simultaneously assessed by injecting a buffer blank after the mAb injection sequence (Figure 5). At a 10-fold column load of 5 µg, carryover appeared to be nonexistent, which can again be attributed to the use of a monolithic support. Carryover at 1% levels became apparent upon a single load of 500 µg of mAb onto the column. This represents the maximum column load and is one typically not encountered in real-life experiments. It is worth noting that carryover was eliminated after the injection of a second buffer blank.


Figure 4. Replicate (n = 10) UV 280 nm Protein A chromatograms of a Herceptin originator diluted in 50 mM Na phosphate, pH 7.4, to 0.5 mg/mL (injection volume 5 µL). Elution was achieved using acetic acid.

Table 1. Retention time and peak area RSD values obtained for the 10-fold analysis of a Herceptin originator at 0.5 mg/mL (5 µL injection volume).

	Acetic acid		Citric acid	
	Peak area	RT (min)	Peak area	RT (min)
1	361	1.669	383	1.666
2	362	1.668	372	1.666
3	373	1.668	365	1.665
4	365	1.669	389	1.667
5	370	1.669	383	1.666
6	373	1.669	378	1.666
7	367	1.671	379	1.678
8	365	1.668	377	1.666
9	366	1.670	376	1.667
10	360	1.670	377	1.667
Mean	366	1.669	378	1.667
S	4.64	0.001	6.52	0.001
% RSD	1.27	0.06	1.73	0.06

The limit of detection (LOD) was around column loads of 0.5 µg. This put some demands on injection volume. If samples have low mAb levels, high volume injections are required. Figure 6 shows the linearity obtained when increasing the injection volume from 5 to 50 µL for a 1 mg/mL Herceptin originator. With the knowledge that 50 µL injections are perfectly feasible and that the lowest detectable amount on-column is 0.5 µg, samples with mAb concentrations at 10 µg/mL are within reach.


Figure 5. UV 280 nm Protein A chromatograms of a Herceptin originator diluted in 50 mM Na phosphate, pH 7.4, to 0.5 mg/mL, and a blank buffer analyzed after a sequence of 10 Herceptin injections. Elution was achieved using acetic acid, and injection volumes were 5 µL.


Figure 6. Linearity obtained by increasing the injection volume of a Herceptin originator (0.5 mg/mL) from 5 to 50 µL.

In mAb titer determination, it is important to be able to assess absolute mAb concentrations. These can be found by linking the peak areas measured in cell-culture supernatants to an external calibration curve constructed by diluting a mAb standard. For the Herceptin biosimilar project, this standard was found in the originator product, which was accurately formulated at 21 mg/mL. The calibration curves of a dilution series of Herceptin originators using citric acid and acetic acid as elution buffers are shown in Figure 7. The corresponding chromatograms are shown in Figure 8. In both cases, linearity was excellent, between 0.02 mg/mL and 2 mg/mL, which is the typical mAb titer range in CHO cells.



Figure 7. Herceptin Protein A calibration curve (0.02 to 2 mg/mL) using citric acid (A) and acetic acid (B) as elution buffer.

Application in mAb titer determination

The method possesses all the characteristics for the determination of mAb titer in cell-culture supernatants. It is fast, precise, and linear in the expected mAb concentration range and does not suffer from carryover. To illustrate this, nine trastuzumab-producing clones, generated in the framework of a Herceptin biosimilar development program, were analyzed using the Bio-Monolith Protein A column to determine absolute mAb concentrations. Chromatograms are displayed in Figure 9, and Table 2 reports the obtained mAb titers using both citric acid and acetic acid as elution buffers. Very consistent data were generated using both elution buffers. These results allow clear decisions to be made early in the development of mAbs. High-producing clones can be readily selected and subjected to further development.


Figure 9. Overlaid UV 280 nm Protein A chromatograms of nine trastuzumab-producing CHO clones using citric acid (A) and acetic acid (B) as elution buffer.


Figure 8. Overlaid UV 280 nm Protein A chromatograms of Herceptin calibration points using citric acid (A) and acetic acid (B) as elution buffer.

Maximizing column lifetime

Column regeneration

A major benefit of using a monolithic disk is that the presence of channels instead of pores decreases the likelihood of column clogging when injecting cell-culture samples. This increases robustness and reduces cleaning efforts. Column contamination can be reduced by running a blank gradient injection after every 30 to 50 samples. If column deterioration is observed (tailing or broad peaks), the following cleaning procedure is recommended. Column regeneration is the first step. If performance is still suboptimal, the clean-in-place procedure can be used, which will reduce the amount of Protein A available.

Column regeneration

1. Wash with 2 mL (20 column volumes (CV)) of 100 mM phosphate buffer + 1 M NaCl, pH 7 to 8, at 0.5 to 1.0 mL/min.
2. Wash with 2 mL (20 CV) of low-pH solution (such as elution buffer).
3. Re-equilibrate with binding buffer.

Clean-in-place

1. Wash with 1 to 2 mL (10 to 20 CV) of 0.1 M NaOH (reverse flow direction) at 0.2 to 0.5 mL/min.
2. Wash with 1 to 2 mL (10 to 20 CV) of DI water at 0.5 to 1.0 mL/min.
3. Wash with 1 to 2 mL (10 to 20 CV) of concentrated buffer (0.1 to 0.5 M) to restore normal pH (7.0 to 7.4).
4. Re-equilibrate with 5 mL (50 CV) of binding buffer.

Conclusions

The Agilent Bio-Monolith Protein A column was successfully applied in the selection of trastuzumab-biosimilar-producing clones based on both titer and structural similarity to the originator. This clone selection process is of utmost importance early in the development of innovator and biosimilar mAbs

For More Information

These data represent typical results.

For more information on our products and services, visit our Web site at www.agilent.com/chem.

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Fast and High Resolution Analysis of Intact and Reduced Therapeutic Monoclonal Antibodies (mAbs)

The Agilent Bio-inert LC and AdvanceBio RP-mAb Columns

Authors

M. Sundaram Palaniswamy
Agilent Technologies Pvt Ltd
Bangalore, India

Abstract

Therapeutic monoclonal antibodies (mAbs) have become the most rapidly growing class of therapeutics in development for many diseases. Novel mAbs are entering clinical trials at a rate of 40 per year. There is also an urgent need for an analytical method that can be used for high-throughput analysis of large number of samples to support the growing biopharma development. This Application Note describes a fast and high-resolution method for the analysis of intact and reduced therapeutic Innovator and Biosimilar mAbs by reverse phase HPLC. Separation was achieved using an Agilent 1260 Infinity Bio-inert LC system with Agilent AdvanceBio RP-mAb C4 and Diphenyl columns. RP-mAb columns give the advantage of superficially porous 3.5 μm particles with 450 \AA wide pores for improved accuracy and short analysis time compared to fully porous particles of the same size. The bio-inertness of the system, together with high resolution, and short and reproducible methods makes it highly suitable for biopharma QA/QC analysis.

Introduction

Evaluating the molecular similarity of a biosimilar to the reference or the innovator molecule is crucial during biosimilar development. A number of physicochemical methods are required by regulatory agencies involving a wide range of comparability programs. The authorities want to see comparability data on platforms that the previous company or the innovator submitted, primarily high-performance liquid chromatography (HPLC), TOF, Q-TOF mass spectrometry, and capillary electrophoresis. HPLC is a well-established technique for the determination of intact protein by size exclusion or ion exchange. However, technological developments in the field of reverse phase (RP) chromatographic stationary phases (a large pore size of 300 Å or fused core particles with short alkyl chains) have made them promising tools for analyzing intact proteins¹.

Historically, mAbs and their fragments are analyzed with limited success using widepore, totally porous particle RP HPLC. Due to their large size and limited diversity, analysis times are often unacceptably long, and mAb peaks can elute as broad bands, compromising resolution. In contrast, high efficiency superficially porous columns easily separate mAbs and their fragments in minutes with high efficiency.

In this work, we have demonstrated the suitability of the Agilent 1260 Bio-inert Quaternary LC system and Agilent AdvanceBio RP-mAb columns to achieve high resolution and rapid separation of intact and fragmented therapeutic innovator and biosimilar rituximab. The unique design of the AdvanceBio RP-mAb column offers unique selectivity due to its superficially porous particles (3.5 µm) with wide pores (450 Å). The column delivers a significant speed and resolution advantage while maintaining compatibility with all instruments.

Experimental

Equipment

A completely biocompatible Agilent 1260 Infinity Bio-inert Quaternary LC System with a maximum pressure of 600 bar consisting of the following modules was used:

- Agilent 1260 Infinity Bio-inert Quaternary LC Pump (G5611A)
- Agilent 1260 Infinity Bio-inert High Performance Autosampler (G5667A)
- Agilent 1200 Infinity Series Thermostat (G1330B)
- Agilent 1260 Infinity Thermostatted Column Compartment containing bio-inert click-in heating elements (G1316C, option 19)
- Agilent 1260 Infinity Diode Array Detector with 60 mm Max-Light high sensitivity flow cell (G4212B option 33)
- Agilent AdvanceBio RP-mAb Diphenyl, 2.1 × 50 mm, 3.5 µm (p/n 799775-944)
- Agilent AdvanceBio RP-mAb C4, 2.1 × 50 mm, 3.5 µm (p/n 799775-904)

The complete sample flow path is free of any metal components, therefore, the sample never contacts metal surfaces. Solvent delivery is free of any stainless steel or iron components.

Software

Agilent ChemStation B.04.03 (or higher)

Reversed-phase HPLC parameters

Chromatographic parameters for intact and reduced mAb analysis using AdvanceBio RP-mAb columns are shown in Table 1.

Table 1. Chromatographic parameters used for intact and reduced analysis.

Parameter	HPLC (intact and reduced mAbs)	
Mobile phase	A) Water + 0.1 % TFA B) IPA:ACN:Water (70:20:10) + 0.09 % TFA	
Columns	Agilent AdvanceBio RP-mAb Diphenyl, 2.1 × 50 mm, 3.5 µm Agilent AdvanceBio RP-mAb C4, 2.1 × 50 mm, 3.5 µm	
Gradient	Time (min)	% B
	0	15
	0.5	25
	1.5	35
	1.51	35
	3.0	60
4.0	60	
Post time	2 minutes	
Injection volume	1 µL	
Flow rate	1.0 mL/min	
TCC	80 °C	
UV detection	220 and 280 nm	

Reagents, samples, and materials

Alnnovator and biosimilar rituximab were purchased from a local pharmacy and stored according to the manufacturer's instruction. PBS and tris(2-carboxyethyl) phosphine (TCEP) were purchased from Sigma-Aldrich. All chemicals and solvents were HPLC grade, and highly purified water from a Milli-Q water purification system (Millipore Elix 10 model, USA) was used. Acetonitrile and 2-propanol were purchased from Lab-Scan (Bangkok, Thailand). For intact and reduced analysis, rituximab samples were diluted to 2 mg/mL using PBS.

Sample preparation

Reduction of mAbs

For the separation of the light and heavy chains, an aliquot of 0.5 MTCEP stock was added to the mAb samples to obtain a final concentration of 10 mM. The mixture was held at 60 °C for 30 minutes.

Results and Discussion

The AdvanceBio RP-mAb column with superficially porous particles and wide pores delivers higher resolution and faster run times to provide accurate, reproducible results when analyzing monoclonal antibodies for biopharma discovery, development, and QA/QC applications. Combined with the Agilent 1260 Infinity Bio-inert Quaternary LC System with a power range up to 600 bar, it can be used for mAb separation. The mobile phase was a combination of isopropanol (IPA), acetonitrile (ACN), water, and trifluoroacetic acid (TFA). Figures 1 and 2 depict the optimized RP HPLC elution profile of intact innovator and biosimilar rituximab on an AdvanceBio RP-mAb Diphenyl, 2.1 × 50 mm, 3.5 μm and AdvanceBio RP-mAb C4, 2.1 × 50 mm, 3.5 μm column, respectively, demonstrating excellent peak shape and fast separation in 4 minutes. Comparing Figures 1 and 2 demonstrates that different selectivity can be obtained through the use of different bonded phases using the same chromatographic conditions, with the diphenyl phase resolving in finer detail.


Figure 1. RP-HPLC analysis of innovator rituximab (A) and biosimilar rituximab (B) separated on an Agilent AdvanceBio RP-mAb Diphenyl 2.1 × 50 mm, 3.5 μm column. C and D show zooms.


Figure 2. RP-HPLC profile of intact innovator rituximab (A) and biosimilar rituximab (B) on an Agilent AdvanceBio RP-mAb C4, 2.1 × 50 mm, 3.5 μm column. C and D show zooms.

Reduced mAb analysis

TCEP was used to separate free antibody light and heavy chains. AdvanceBio RP-mAb columns are very effective in providing fast and high-resolution separations of antibody fragments. The profiles in Figures 3 and 4 show a rapid reversed-phase analysis optimized for the separation of antibody fragments in approximately 4 minutes using C4 and diphenyl phases, respectively.

In both cases, due to reduction of the disulfide bonds, mAbs eluted as distinct light chain (LC) and heavy chain (HC) separations with high efficiency. The same gradient used for the intact analysis was able to resolve the LC and HC for the reduced samples. As we have seen with intact mAb analysis, the LC and HC show different selectivities with diphenyl and C4 columns. RP HPLC analysis of intact and reduced innovator and biosimilar using AdvanceBio RP-mAb diphenyl and C4 columns indicates that the mAb pair are highly similar.


Figure 3. RP-HPLC profiles of innovator rituximab (A) and biosimilar rituximab (B) separated on an Agilent AdvanceBio RP-mAb Diphenyl, 2.1 × 50 mm, 3.5 μm column. Mirror plot image overlays (C).


Figure 4. RP-HPLC profiles of innovator rituximab (A) and biosimilar rituximab (B) separated on an Agilent AdvanceBio RP-mAb C4, 2.1 × 50 mm, 3.5 μm column. Mirror plot image overlays (C).

Precision of retention time and area

Tables 2 and 3 present the average retention times and area RSDs from six replicates of intact and TCEP reduced innovator, and biosimilar rituximab for the diphenyl and C4 phases. The results show that both columns provide precision of RT and area within the acceptable limit of $\pm 3\%$ and $\pm 5\%$, respectively.

Conclusion

In this application note, we have demonstrated a simple LC-UV-based approach to define the molecular similarity between a biosimilar and its innovator reference. We first used the Agilent 1260 Bio-inert Quaternary LC system with Agilent AdvanceBio RP-mAb Diphenyl and C4 columns to develop a high-resolution and rapid separation of intact mAbs. Using the same method, we were also able to show the separation of light chain and heavy chain after TCEP reduction. Area and RT precision of intact and reduced analysis using AdvanceBio RP-mAb columns were excellent, and show the reliability of the method. Such fast, simple, and reproducible methods for intact and reduced analysis of mAbs, coupled with bio-inertness of the system makes this solution suitable for the comparability analysis of mAbs for the biopharma industry.

Reference

1. Navas, N; *et al.*, Anal. Bioanal. Chem. 2013, 405, pp 9351-9363.

Table 2. Retention time and peak area RSD (%), n = 6 for intact analysis

Samples	Retention time		Peak area	
	Mean (min)	RSD	Mean (mAU/min)	RSD
Agilent AdvanceBio RP-mAb, C4, 2.1 × 50 mm, 3.5 μm				
Innovator rituximab	1.96	0	71.61	1.98
Biosimilar rituximab	1.95	0.26	77.3	0.47
Agilent AdvanceBio RP-mAb, Diphenyl, 2.1 × 50 mm, 3.5 μm				
Innovator rituximab	2.51	0.20	66.7	0.458
Biosimilar rituximab	2.51	0	73.3	1.86

Table 3. Retention time and peak area RSD (%), n = 6 for reduced analysis

Samples	Retention time		Peak area	
	Mean (min)	RSD	Mean (mAU/min)	RSD
Agilent AdvanceBio RP-mAb, Diphenyl, 2.1 × 50 mm, 3.5 μm				
Innovator rituximab LC	2.32	0.60	19.71	4.24
Innovator rituximab HC	2.58	1.52	57.33	1.57
Biosimilar rituximab LC	2.32	0.07	23.56	3.25
Biosimilar rituximab HC	2.60	0.05	58.40	5.61
Agilent AdvanceBio RP-mAb, C4, 2.1 × 50 mm, 3.5 μm				
Innovator rituximab LC	1.83	0	21.5	1.4
Innovator rituximab HC	2.03	0.04	51.2	2.25
Biosimilar rituximab LC	1.83	0.03	24.47	3.84
Biosimilar rituximab HC	2.03	0.06	52.66	0.84


Disulfide Linkage Analysis of IgG1 using an Agilent 1260 Infinity Bio-inert LC System with an Agilent ZORBAX RRHD Diphenyl sub-2 μm Column

Authors

M. Sundaram Palaniswamy
Agilent Technologies, Inc.
Bangalore, India

Abstract

This Application Note describes a simple method for the analysis of disulfide linkages in monoclonal antibodies (mAbs) by reversed-phase HPLC. Separation was achieved using an Agilent 1260 Infinity Bio-inert LC System and an Agilent ZORBAX RRHD 300 Diphenyl sub-2 μm particle column. Diphenyl 1.8 μm columns deliver UHPLC performance for reversed-phase separations of intact proteins and peptide digests. Together with UHPLC instruments, these versatile columns enable higher order characterization with shorter analysis times. The 1260 Infinity Bio-inert LC System has a power range up to 600 bar and is capable of handling the higher pressures demanded by emerging column technologies with smaller particles down to 1.7 μm .

Introduction

Although recombinant mAb therapeutics have advanced enormously in recent years, little is known about their disulfide bond patterns. Complete disulfide bond assignment of IgG1 antibodies can be challenging due to their large size and substantial number of disulfide linkages. Disulphide bonding is critical to maintaining immunoglobulin (IgG) tertiary and quaternary structure for therapeutic monoclonal antibodies (mAb). Both inter- and intra-chain disulphide bonds are formed intracellularly in the expression host prior to secretion and purification during mAb production processes. Disulphide bond shuffling has previously been reported for IgG2 and disulphide-mediated arm-exchange for IgG4, reflecting innate behaviour of these IgG classes^{1, 2}. However, a typical and significant reduction in the number of disulphide bonds has been observed in IgG13 that present significant issues for manufacturing of therapeutic mAbs. This Application Note demonstrates the suitability of the 1260 Bio-inert Quaternary LC System for separating and analyzing the disulfide linkages of IgG1 by reversed-phase HPLC on ZORBAX RRHD 300 Diphenyl, 2.1 × 100 mm, 1.8 μm column. Ultrahigh performance liquid chromatography (UHPLC) separation using sub-2 μm particles improves resolution per time and sensitivity, shortens run times, and thus enables the analysis of IgG1, reduced IgG1, and the peptides resulting from digestion of IgG1.

Equipment

Instrumentation

A completely biocompatible Agilent 1260 Infinity Bio-inert Quaternary LC System with a maximum pressure of 600 bar consisting of the following modules was used:

- Agilent 1260 Infinity Bio-inert Quaternary LC Pump (G5611A)
- Agilent 1260 Infinity Bio-inert High Performance Autosampler (G5667A)
- Agilent 1200 Infinity Series Thermostat (G1330B)
- Agilent 1260 Infinity Thermostatted Column Compartment containing bio-inert click-in heating elements (G1316C, option 19)
- Agilent 1260 Infinity Diode Array Detector with 60 mm Max-Light high sensitivity flow cell (G4212B option 33)
- Agilent ZORBAX RRHD 300 Diphenyl, 2.1 × 100 mm, 1.8 μm column (p/n858750-944).

The complete sample flow path is free of any metal components such that the sample never gets in contact with metal surfaces. Solvent delivery is free of any stainless steel or iron components.

Software

Agilent OpenLAB CDS ChemStation Edition for LC & LC MS Systems, Rev. C.01.04

HPLC analysis

Table 1 Chromatographic parameters used for RP HPLC.

Parameter	Conditions	
Mobile phase A	Water + 0.1% TFA	
Mobile phase B	Acetonitrile + 0.09% TFA	
Gradient	Time (min)	Mobile phase B (%)
	0 minutes	25
	3 minutes	35
	4 minutes	40
	5 minutes	40
	15 minutes	90
	16 minutes	25
	Post time	5 minutes
Injection volume	3 μL (Needle with wash, flush port active for 7 seconds)	
Flow rate	0.3 mL/min	
Data acquisition	214 and 280 nm	
Acquisition rate	20 Hz	
Flow cell	60 mm path	
Column oven	50 °C	
Sample thermostat	5 °C	

Reagents, Samples and Materials

The human monoclonal antibody IgG1 was a proprietary pharmaceutical molecule. DL-Dithiothreitol (DTT), iodoacetamide, trizma base, and Endoproteinase Lys C were purchased from Sigma Aldrich. All chemicals and solvents used were HPLC grade and highly purified water from Milli Q water purification system (Millipore Elix 10 model, USA) was used. Acetonitrile was of gradient grade and purchased from Lab-Scan (Bangkok, Thailand).

Reduction and alkylation of Intact IgG1

IgG1 was diluted to 2 mg/mL using 100 mM Tris HCl and 4 M Gu HCl, pH 8.0. An aliquot of 10 μ L of 0.5 M DTT stock was added to obtain a final concentration of 5 mM. The mixture was held at 37 $^{\circ}$ C for 30 minutes. The reaction mixture was cooled briefly to room temperature (RT). An aliquot of 26 μ L of 0.5 M Iodoacetamide stock was added for a final concentration of 13 mM. It was allowed to stand for 45 minutes. Once removed, the solution was quenched with 20 μ L of DTT for a final concentration of 10 mM.

Lys C digestion of IgG1 and reduction

IgG1 was diluted to 1 mg/mL using 100 mM Tris HCl, pH 8.0. Endoproteinase Lys C in 100 mM Tris HCl, pH 8.0 was added at an enzyme protein ratio of 1:100 (w/w). The mixture was incubated overnight at 37 $^{\circ}$ C. The reaction was stopped by lowering the pH to 6.0 by adding 10 % TFA. Later, the reduction of Lys C digested IgG1 was carried out as described earlier in this Application Note.

Results and Discussion

Separation and Detection

A ZORBAX RRHD 300 Diphenyl 1.8 μ m column has the advantage of low pH and temperature stability, and, combined with the 1260 Infinity Bio-inert Quaternary LC System with a power range up to 600 bar and capabilities of handling the higher pressures, can be used for protein separation. Figure 1 **A** depicts the optimized RP HPLC elution profile of intact IgG1 on a ZORBAX RRHD 300 Diphenyl, 2.1 \times 100 mm, 1.8 μ m column demonstrating excellent retention of IgG1 in 15 minutes. The reproducibility of analysis was tested with six replicates. Figure 1 **B** shows the overlay of six replicates.


Figure 1. RP HPLC profile of intact IgG1 on an Agilent ZORBAX RRHD 300 Diphenyl, 2.1 \times 100 mm, 1.8 μ m column (A), and an overlay of six replicates (B).

Table 2 Retention time and area RSD (%), n = 6 for intact IgG1

Retention time		Peak area	
Mean (min)	RSD (limit: \pm 3.0 %)	Mean (mAU/min)	RSD (limit: \pm 5.0 %)
8.838	0.086	1,170	0.461

The effect of reduction and alkylation of the disulfide bonds in intact IgG1 was tested. Figure 2 shows the reversed-phase chromatogram of A reduced and alkylated IgG1 B overlay with reduction/alkylation buffer blank and C overlay of six replicates showing separation reproducibility. Due to the reduction of the disulfide bonds, the IgG1 is separated into its light and heavy chains. The IgG1 eluted as distinct light chain (LC) and heavy chain (HC) as indicated in Figure 2; however, this was not confirmed by mass spectral analysis.


Figure 2. RP HPLC profile of (A) Reduced and alkylated IgG1, (B) overlaid with buffer blank, (C) overlay of six replicates.

Table 3. Retention time and area RSD (%), n = 6 for Light chain

Retention time		Peak area	
Mean (min)	RSD (limit: ± 3.0 %)	Mean (mAU/min)	RSD (limit: ± 5.0 %)
8.638	0.091	504.33	2.780

Table 4. Retention time and area RSD (%), n = 6 for Heavy chain

Retention time		Peak area	
Mean (min)	RSD (limit: ± 3.0 %)	Mean (mAU/min)	RSD (limit: ± 5.0 %)
9.084	0.152	1,520	0.390

Peptide maps resulting from Lys C digestion of intact IgG1 under nonreducing conditions resulted in a less intricate RP HPLC profile. A representative chromatogram of the IgG1 digest (Figure 3 A) displays the two (baseline separated) peaks that were selected for area and RT precision. The overlay results in sharp peaks with good resolution and excellent separation reproducibility (Figure 3 B).

Further, we wanted to compare the reversed-phase profile of IgG1 under nonreduced and reducing conditions to determine the peptides bound through disulfide linkages. The overlay of Lys C peptide maps of nonreduced IgG1 (red trace) and reduced IgG1 (blue trace) is depicted in Figure 4. The appearance of additional peaks (indicated by an asterix) after reduction of Lys C digested IgG1 confirms they are bound through disulfide linkages.

Precision of retention time and area

The precision of the retention time and area for intact IgG1, reduced IgG1 and endoproteinase Lys C digested IgG1 under nonreduced conditions are given in Tables 2, 3, and 4. The results show that the ZORBAX diphenyl sub-2 μm column shows precision of RT and area to be within 3% and 5 % respectively.


Figure 3. RP HPLC profile of (A) Lys C digested IgG1 and (B) overlay of six replicates. Peaks selected for RT and area RSDs are indicated by an asterix.

Table 5 Retention time and area RSD (%), n=6 for Lys C digested IgG1

	Retention time		Peak area	
	Mean (min)	RSD (limit: $\pm 3.0\%$)	Mean (mAU/min)	RSD (limit: $\pm 5.0\%$)
Peak 1	5.525	0.307	60	0.544
Peak 2	7.444	0.140	132.45	1.113


Figure 4. Comparison of peptide maps of Lys C digested IgG1 under nonreducing condition (red trace) followed by reduction (blue trace). Peptides bound through disulfide linkages are indicated by an asterix.

Conclusion

Disulfide linkage analysis is important to study some of the post-translational modifications of proteins for biopharma process development and monitoring. We have shown the combination of an Agilent 1260 Infinity Bio-inert Quarternary LC System and an Agilent ZORBAX RRHD 300 Diphenyl, 2.1 \times 100 mm, 1.8 μ m column to perform reproducible and high resolution analysis of disulfide linkage analysis of monoclonal antibodies for biopharma process development and monitoring. Area and RT precision of the method were excellent and proved reliability. Further, the 1260 Infinity Bio-inert LC has a power range up to 600 bar and is capable of handling the higher pressures demanded by emerging column technologies with smaller particles down to 1.7 μ m. The bio-inertness and corrosion resistance of the instrument coupled with a simple and reproducible method make this solution particularly suitable for the QA/QC analysis of monoclonal antibody for the biopharmaceutical industry.

References

1. R. Mhatre, J. Woodard, C. Zeng, Strategies for locating disulfide bonds in a monoclonal antibody via mass spectrometry, *Rapid Commun. Mass Spectrom.* 13 (1999) 2503–2510.
2. T.-Y. Yen, H. Yan, B.A. Macher, Characterizing closely spaced, complex disulfide bond patterns in peptides and proteins by liquid chromatography/electrospray ionization tandem mass spectrometry, *J. Mass. Spectrom.* 37 (2002) 15–30.
3. Mullan et al. *BMC Proceedings* 2011, 5 (Suppl 8):P110

Determination of Drug-to-Antibody Distribution in Cysteine-Linked ADCs.


An Analysis of ADCs of IgG1 and IgG2 Subclasses

Authors

Lily Liu
Agensys, Inc.
Santa Monica, CA, US
Andy Coffey
Agilent Technologies, Inc

Abstract

This Application Note shows the advantages of using sub-2 μm reversed-phase columns to gain extra resolution and accuracy in the determination of drug-to-antibody ratios. These ratios were determined in antibody-drug conjugates derived from antibody intermediates of the IgG1 and IgG2 subclasses in this study.

Introduction

Antibody-drug conjugates (ADCs) represent a growing class of anticancer therapeutics that combine the specificity of an antibody with the potency of chemotherapeutic agents using covalent and chemically stable linkages. The ADC field is expanding with an increasing number of conjugation technologies being developed. One dominant class of ADCs includes conjugation to cysteine residues that are involved in the formation of interchain disulfide bonds through maleimide linkages. One of the principal critical quality attributes for ADCs that directly correlates with potency is the drug-to-antibody ratio (DAR). Up to eight or 12 drug-linkers may be conjugated per antibody, depending on the IgG antibody subclass¹.

Hydrophobic interaction chromatography (HIC) is a common approach for determining conjugate distribution, and calculating DAR for ADCs manufactured from IgG1 mAbs. However, monitoring conjugate distribution and DAR for ADCs manufactured from IgG2 mAbs by HIC is challenging, due to incomplete resolution between positional isomers and variably conjugated species. Reversed-phase (RP) chromatography can be used as an alternative or orthogonal technique for determining the DAR of ADCs following reduction of interchain disulfide bonds. Using this technique, the DAR may be calculated experimentally from the distribution of unconjugated and conjugated light and heavy chains. For ADCs manufactured from IgG2 antibodies, RP is a more suitable method. This is because elution between unconjugated and variably conjugated light and heavy chains is dictated by the number of conjugated drug-linkers, regardless of the site of conjugation.

A limited variety of suitable HIC and RP columns are available for these applications. This Application Note describes the use of the Agilent ZORBAX RRHD SB300-C8 column for characterizing the distribution of unconjugated and variably conjugated light and heavy chains, and for determining the average DAR. Here, we describe RP UHPLC methods suitable for ADCs manufactured from both IgG1 and IgG2 antibodies. Compared to methods using common HPLC columns, the ZORBAX RRHD SB300-C8 column offers improved peak resolution, and yields similar distributions of unconjugated and conjugated light and heavy chains and DARs.

Materials and Methods

Reagents, samples, and materials

ADCs manufactured from fully human IgG1 and IgG2 antibody intermediates are proprietary. DL-dithiothreitol (DTT) was purchased from Thermo Scientific (Pierce NoWeigh DTT). All solvents used were HPLC grade, and were purchased from either VWR or Fisher Scientific.

Samples in their respective formulation buffers (pH 5–6) were diluted to 5 mg/mL, and the pH was adjusted to approximately pH 8 with 1 M Tris pH 9. Partial reduction of the interchain disulfide bonds was achieved by incubation in 40 mM DTT at 37 °C for 15 minutes. After cooling to room temperature, reduced samples were diluted 1:1 with 2 % formic acid in 50 % acetonitrile to quench the reduction reaction.

UHPLC method

Parameter	Value
Column	Agilent ZORBAX RRHD SB300-C8, 50 mm × 2.1, 1.8 μm
Other columns	Vydac 214MS, C4, 2.1 × 50 mm, 5 μm, 300 Å Agilent PLRP-S, 2.1 × 50 mm, 5 μm, 1,000 Å
Mobile phases	A) 0.1 % TFA in H ₂ O B) 0.08 % TFA in 90 % ACN
Column temperature	80 °C (IgG1) 70 °C (IgG2)
Post-column cooler	35 °C
Injection volume	2 μL (IgG1) 3 μL (IgG2)
Flow rate	1 mL/min (IgG1) 0.8 mL/min (IgG2)
Detection	UV at 214 and 280 nm
Autosampler temperature	10 °C
IgG1 Gradient	Time (min) %B 0 34.5 3 38.0 5.5 38.5 25 55.0 25.1 75.0 26 75.0 26.1 34.5 Post time: 4 minutes
IgG2 Gradient	Time (min) %B 0 30.0 3 30.0 21 45.0 21.1 75.0 22 75.0 22.1 30.0 Post time: 2 minutes

Instruments

UHPLC with DAD detection system from an external vendor.

The equivalent Agilent instrument for UHPLC analysis is the Agilent 1290 Infinity II LC system, which is expected to deliver comparable, or better performance

Peak assignments

Peak identities were confirmed by coupling the UHPLC with in-line mass spectrometry (data not shown). The major peaks corresponded to unconjugated and variably conjugated light and heavy chains. Peaks eluting as trailing shoulders from the major peaks were identified as having one or more intrachain disulfide bonds reduced. Multiple peaks corresponding to conjugated heavy chains were observed having the same mass, and were identified as being positional isomers, where the drug-linker was conjugated at different cysteine residues.

DAR calculation

The DAR value was calculated from the analysis of the UV chromatogram, using Equation 1.

$$\text{DAR} = 2 \left(\sum_{n=0}^1 \frac{\text{LC peak area} \times n_{\text{drug}}}{\text{Total LC peak area}} + \sum_{n=0}^j \frac{\text{HC peak area} \times n_{\text{drug}}}{\text{Total HC peak area}} \right)$$

Figure 1.


Figure 1. Comparison of reduced RP profiles obtained using a C4 HPLC column (A) and an Agilent ZORBAX RRHD SB300-C8 (B) for a cysteine-conjugated IgG2 ADC.

Results and Discussion

Figure 1 compares the reversed-phase chromatographic profiles of an ADC manufactured from an IgG2 mAb intermediate using a conventional C4 HPLC column and the ZORBAX RRHD SB300-C8 column. The peaks observed in the chromatogram corresponding to the UHPLC method using the ZORBAX RRHD SB300-C8 column are sharper, and show better resolution compared to the C4 HPLC column. Peak separation and resolution achieved using the ZORBAX RRHD SB300-C8 column enabled improved peak integration accuracy, and the DAR value was calculated to be 0.1 higher as a result. resolution achieved using the ZORBAX RRHD SB300-C8 column enabled improved peak integration accuracy, and the DAR value was calculated to be 0.1 higher as a result.

Figure 2 shows RP chromatograms of a reduced ADC manufactured from an IgG1 mAb resulting from analyses using both the Agilent PLRP-S column and the ZORBAX RRHD SB300-C8 column. Separation of the unconjugated and conjugated light and heavy chains on the ZORBAX RRHD 300SB-C8 column resulted in sharper peaks and improved resolution of the minor species in comparison to the PLRP-S HPLC column. The DAR value calculated from the results using the ZORBAX RRHD column was 0.1 higher than the result from the PLRP-S column. This difference can be attributed to more accurate peak integration.


Figure 2. Comparison of PLRP-S HPLC (A) and Agilent ZORBAX RRHD SB300-C8. 1.8 μm (bottom) profiles for a cysteine-conjugated IgG2 ADC.

Distributions and DAR value

Table 2 shows the distribution of unconjugated and conjugated light and heavy chains, and the calculated average DAR for the ZORBAX RRHD column and the two different HPLC columns. Relative peak areas and concomitant DAR values calculated using results obtained with the ZORBAX RRHD column are similar to those determined using the HPLC columns. The major advantage of using a sub-2 μm column is the improved peak resolution. This improved resolution results in better separation of minor species that were coeluting with the main peaks using the HPLC columns. The gradients presented in this Application Note can also be used to detect minor changes in sample stability and characterization of minor peaks by LC/MS. Also, the increased peak resolution and the higher pressure compatibility of the ZORBAX RRHD columns allow for development of shorter gradients for high-throughput testing.

Conclusion

Reversed-phase methods can determine the distribution of unconjugated and conjugated light and heavy chains and for calculation of the DAR for ADCs.

The Agilent ZORBAX RRHD SB300-C8 column has been shown to be suitable for these purposes to support characterization of ADCs derived from both IgG1 and IgG2 antibody intermediates. Peaks eluting from the ZORBAX RRHD SB300-C8 column were observed to be sharper and better resolved in comparison to the HPLC columns included in the comparison. This improvement enabled more accurate peak integrations and concomitant DAR values. The improvement in peak resolution from the ZORBAX RRHD SB300-C8 column also allowed detection of conjugation site positional isomers. RP methods using the ZORBAX RRHD SB300-C8 column have been shown to be suitable for characterizing the conjugate distribution of partially reduced ADCs in place of, or orthogonal to, hydrophobic interaction chromatography.

Reference

1. Wiggins, B.; *et al.* Characterization of Cysteine-Linked Conjugation Profiles of Immunoglobulin G1 and Immunoglobulin G2 Antibody–Drug Conjugates. *J. Pharm. Sci.* **2015**, 104(4), 1362-1372.

Table 2. Relative peak areas of each species and average DAR value using an Agilent ZORBAX RRHD SB300-C8 compared with HPLC columns.

Peak ID	IgG2 ADC		IgG1 ADC	
	% by HPLC (C4)	% by UHPLC (Agilent ZORBAX RRHD)	% by HPLC (PLRP-S)	% by UHPLC (Agilent ZORBAX RRHD)
L0	22.6	23.7	14.3	13.3
L1	2.3	3.4	15.6	17.7
H0	17.3	15.6	16.8	15.4
H1	12.5	13.1	26.2	27.8
H2	32.4	30.9	16.8	15.8
H3	7.4	8.1	10.5	10.0
H4	2.1	3.2	–	–
H5	3.3	2.2	–	–
DAR	3.5	3.6	3.6	3.7

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PLRP-S

資料番号	タイトル
5991-6263EN	Drug-to-Antibody Ratio (DAR) Calculation of Antibody-Drug Conjugates (ADCs) Using Automated Sample Preparation and Novel DAR Calculator Software
5991-6559EN	Measuring Drug-to-Antibody Ratio (DAR) for Antibody-Drug Conjugates (ADCs) with UHPLC/Q-TOF

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
High-Resolution Mapping of Drug Conjugated Peptides in an ADC Digest

Peptide map comparison of mAb and drug conjugated mAb

Authors

Suresh Babu C.V.
Agilent Technologies Inc.

Introduction

Currently, antibody drug conjugates (ADCs) are prime protein drugs for biotherapeutic use. When a cytotoxic drug is conjugated to a biotherapeutic monoclonal antibody (mAb), there are several options for the conjugation site. As part of the characterization of ADCs, it is important to be able to identify these conjugation sites. This can be done using peptide mapping. The specificity of the enzyme to cleave the mAb into peptide fragments results in different cleavage patterns, and, hence, peptide fragments, around the conjugation site. High-resolution peptide mapping can be used to identify peptides that are produced as a result of conjugation of the cytotoxic drug. This Application Note demonstrates the use of the Agilent AdvanceBio Peptide Mapping Column and an Agilent 1290 Infinity LC system for ADC peptide mapping analysis. For analysis using UV, it is important to have high resolution to identify the individual peptides, therefore, the method was developed with an optimized flow rate and gradient time for increased peak capacity. Comparison of a peptide map of Trastuzumab biotherapeutic mAb and its cytotoxic drug conjugate, ADC, revealed the peptide map differences corresponded to drug-conjugated peptides. These hydrophobic peptides were resolved on the AdvanceBio Peptide Mapping Column.

Materials and Methods

Therapeutic proteins, ADC, and Trastuzumab were purchased from a local pharmacy. All chemicals and solvents were HPLC grade. Tryptic digestion of mAbs was carried out as described elsewhere¹. Before the digestion of the mAbs with trypsin, the disulfides were reduced and alkylated under denaturing conditions.

An Agilent 1290 Infinity LC system with the following configuration was used for the study:

- Agilent 1290 Infinity Binary Pump with integrated vacuum degasser (G4220A)
- Agilent 1290 Infinity Autosampler (G4226A)
- Agilent 1290 Infinity Thermostatted Column Compartment (G1316C)
- Agilent 1290 Infinity Diode Array Detector (DAD) (G4212A) with 10 mm Max-Light flow cell (G4212-60008)
- AdvanceBio Peptide Mapping Column (p/n 651750-902)

Conditions

Parameter	Value
Column	Agilent AdvanceBio Peptide Mapping, 2.1 × 250 mm, 2.7 μm (p/n 651750-902)
Mobile phase	A) 0.1 % TFA in water B) 0.08 % TFA in 90 % ACN
Gradient	Time (min) %B 0 3 60 50 65 90 66 90 70 3
Injection volume	5 μL (10 μg/μL)
Flow rate	0.5 mL/min
Data acquisition	210 nm/4 nm, 252 nm/4 nm
Thermostatted column compartment	60 °C
Sample thermostat	5 °C
Postrun time	10 minutes

Results and Discussion

Peak capacity is often used as an evaluation criterion to measure the performance of a column under given chromatographic conditions. mAb digestion produces many peptides for analysis. Therefore, it is necessary to develop a method that can increase peak capacity. Also, peak capacity is essential in a peptide mapping study so small impurity peaks or sample heterogeneity can be addressed. Due to the heterogeneous nature of ADC with glycosylation and cytotoxic drug conjugates, tryptic-digested ADC will generate more complex peptides than unconjugated mAb. To monitor the tryptic-digested ADC peptide mixture, optimization of gradient time and flow rates is critical to achieve high peak capacity.

The peak capacity values were calculated by dividing the gradient time by the average peak width of five peptide peaks at baseline (5σ). Figure 1 depicts the effect of gradient time and flow rate on peak capacity. The results suggested that a 0.5 mL/min flow rate and 60 minutes gradient time gave the highest peak capacity values for the 2.1 × 250 mm, 2.7 μm column. These would, therefore, be the optimum conditions for identifying the peptides that have cytotoxic drug conjugation with high resolution.


Figure 1. Effect of flow rate and gradient time on peak capacity.

Figure 2 shows the peptide map of the tryptic digested ADC. The peptide map shows excellent performance with baseline separation and resolution across the entire gradient profile. There was a significant improvement in separation time with the 250 mm column (60 minutes) compared to a 150 mm column (220 minutes) as previously reported in the literature for peptide mapping of ADC2. A peak capacity value of 354 was obtained and the RSD values demonstrate the excellent reproducibility of retention time and peak area and, thus, the precision of the system (Table1).

To identify peptides that have the cytotoxic drug attached, the peptide digests of the mAb and its conjugate, ADC, were analyzed by monitoring the UV trace at 252 nm (Figure 3). Peptide maps of ADC are different from those of Trastuzumab. It is clearly evident that the more hydrophobic drug-bonded peptides in ADC are eluted later (~ 40 to 60 minutes). Comparing the two peptide maps shows a group of later-eluting peptides identified in the ADC digest that are not present in the digest of the mAb. These hydrophobic peptides are the ones with the cytotoxic drug conjugation.

Table 1. RSD of retention time and area (n = 5) of peaks shown in Figure 2.

	Mean RT (min)	RSD RT (%)	Mean area (mAU/min)	RSD Area (%)
Peak 1	5.37	0.13	369.2	0.76
Peak 2	14.27	0.06	106.1	1.66
Peak 3	28.84	0.02	202.61	0.09
Peak 4	35.86	0.02	193.83	0.58


Figure 2. Peptide map of tryptic-digested ADC separated on an Agilent AdvancedBio Peptide Mapping column (*peaks selected for RSD calculations).


Figure 3. Overlay of peptide map of tryptic-digested ADC and Trastuzumab separated on an Agilent AdvancedBio Peptide Mapping Column.

Conclusions

High-resolution peptide maps are obtained when the 250 mm Agilent AdvanceBio Peptide Mapping Column is used with an Agilent 1290 Infinity LC System. By comparing the peptide maps of the mAb and its conjugate, ADC, it is possible to identify the peptides conjugated with a cytotoxic drug in the ADC digest. Additionally, we demonstrated that the AdvanceBio Peptide Mapping Column provided resolution across the range of peptide types. Good peak shapes and reproducibility were obtained for the analysis of the more hydrophobic conjugated peptides, enabling identification and quantitation.

References

1. Gudihal, R.; Waddell, K. Identification of Oxidation Sites on a Monoclonal Antibody Using an Agilent 1260 Infinity HPLC-Chip/MS System Coupled to an Accurate-Mass 6520 Q-TOF LC/MS; *Agilent Technologies Application Note*, publication number 5990-8768EN, **2011**.
2. Kim, M. T.; *et al.* Statistical Modeling of the Drug Load Distribution on Trastuzumab Emtansine (Kadcyla), a Lysine-Linked Antibody Drug Conjugate. *Bioconjugate Chem.* **2014**, *25*, pp 1223-1232.

Analysis of Tryptic Digests of a Monoclonal Antibody and an Antibody-Drug Conjugate with the Agilent 1290 Infinity II LC

Authors

Gerd Vanhoenacker,
Mieke Steenbeke,
Isabel Vandenheede,
Pat Sandra, and Koen Sandra

Research Institute for
Chromatography

President Kennedypark 26
B-8500 Kortrijk Belgium

Udo Huber and
Sonja Schneider

Agilent Technologies, Inc.
Waldbronn, Germany.

Abstract

An Agilent 1290 Infinity II LC in combination with an Agilent AdvanceBio Peptide Mapping column was used to analyze tryptic digests of the monoclonal antibody (mAb) trastuzumab, or trade name Herceptin, and the antibody-drug conjugate (ADC) ado-trastuzumab emtansine, or trade name Kadcyla. The use of the highly efficient Agilent V380 Jet Weaver mobile phase mixer reduced the baseline noise caused by the trifluoroacetic acid (TFA) modifier significantly, enabling detection of low abundant peptides by UV at 214 nm. By modifying the gradient steepness, peak capacities of approximately 300, 450, and 900 could be obtained for total analysis times of 25, 45, and 205 minutes, respectively.

Introduction

Monoclonal antibodies (mAbs) have emerged as important therapeutics for the treatment of cancer and autoimmune diseases, among others^{1,2}. The successes of mAbs have triggered the development of various next-generation formats including antibody-drug conjugates (ADCs), which combine a specific mAb and a cytotoxic drug by a stable linker^{1,2}. The promise of ADCs is that highly toxic drugs can selectively be delivered to tumor cells, thereby substantially lowering side effects typically experienced with classical chemotherapy. Peptide mapping is an important methodology in the analysis and characterization of these molecules. Hundreds of peptides with varying physicochemical properties present in a wide dynamic concentration range exist in mAb and ADC tryptic digests, demanding the best in terms of separating power. This Application Note describes how an Agilent 1290 Infinity II LC in combination with an Agilent AdvanceBio Peptide Map column is successfully applied to tackle these challenging separations.

Experimental

Instrumentation

An Agilent 1290 Infinity II LC was used, comprising:

- Agilent 1290 Infinity II High Speed Pump (G7120A)
- Agilent 1290 Infinity II Multisampler (G7167B)
- Agilent 1290 Infinity II Multicolumn Thermostat (G7116B)
- Agilent 1290 Infinity II Diode Array Detector (G7117B) equipped with a 10 mm flow cell

An Agilent 1290 Infinity II LC was used, comprising:

- Agilent Jet Weaver mixer, both 35 and 100 μ L configurations were tested (G4220-60006)
- Agilent High-Performance Jet Weaver mixer, 380 μ L (G4220-60012)

Samples and sample preparation

Trastuzumab (Herceptin) and ado-trastuzumab emtansine (Kadcyla) were obtained from Roche (Basel, Switzerland).

A 100 μ g amount of protein, diluted in 0.05 % Rapigest/100 mM Tris-HCl, pH 8, was reduced at 60 °C for 30 minutes by the addition of 5 mM dithiothreitol, and alkylated at 37 °C for 1 hour by adding 10 mM iodoacetamide. Trypsin was subsequently added at an enzyme-to-substrate ratio of 1:25 (w:w). Digestion proceeded for 16 hours at 37 °C. The final concentration was 0.5 μ g/ μ L.

Method parameters

Parameter	Value
Column	Agilent AdvanceBio Peptide Mapping column, 2.1 \times 250 mm, 2.7 μ m (p/n 651750-902)
Mobile phase	A) 0.05 % TFA in water/acetonitrile 99:1 (v:v) B) 0.045 % TFA in acetonitrile
Flow rate	0.35 mL/min
Gradient	0 to 60 %B in various gradient slopes (see peak capacity) 60 to 90 %B in 0.5 minutes and hold for 4.5 minutes 7 minutes post time at 0 %B Example for 40 minutes gradient: 0 to 40 minutes – 0 to 60 %B 40 to 40.5 minutes – 60 to 90 %B 40.5 to 45 minutes – 90 %B
Temperature	60 °C
Injection	5 μ L Needle wash flush port, 5 seconds, 0.05 % TFA in water/acetonitrile 20:80 (v:v)
Detection	Signal 214/4 nm, Reference 360/60 nm Signal 252/4 nm, Reference 360/60 nm (for ADC) >0.025 minutes (0.5 seconds response time) (10 Hz)

Results and Discussion

Agilent AdvanceBio Peptide Mapping column

The Agilent AdvanceBio Peptide Mapping column is packed with 2.7 μm superficially porous C18 particles with 120 \AA pore size. It is a state-of-the-art column for peptide mapping, enabling high-resolution separations in short analysis times. Figure 1 shows the results for a fast analysis (25 minutes total analysis time) of a tryptic digest of Herceptin.

Agilent Jet Weaver mixers

For comprehensive peptide mapping, the system should be able to detect high as well as low abundant peptides. Peptide mapping is generally carried out with UV detection at 214 nm and a water/acetonitrile mobile phase containing trifluoroacetic acid (TFA) because of its beneficial effects on peptide retention and peak shape. It is known that these conditions contribute to increased baseline noise due to the UV absorption of the TFA modifier³. The noise will depend on column dimensions and flow, gradient slope, and system gradient formation.

When mobile phase mixing is inadequate, small variations present in the mobile phase composition may persist after passage through the column, and reach the detector. In the case of TFA, which has significant UV absorbance at low wavelengths (for example, 214 nm), small fluctuations in TFA concentration and water/acetonitrile ratio will be visible with UV or DAD. Excessive noise can be tackled by increasing the mobile phase mixing performance to stabilize solvent composition as much as possible.

The impact of mixing has been tested with the selected column and samples by comparison of three state-of-the-art Jet Weaver mobile phase mixers.


Figure 1. Analysis of a Herceptin tryptic digest, 20-minute gradient, UV detection at 214 nm.

- Agilent V35 Jet Weaver: internal volume of 35 μL , most commonly used because of the small impact on delay volume
- Agilent V100 Jet Weaver: internal volume of 100 μL , used for applications that require higher mixing performance and low delay volume
- Agilent V380 High-Performance Jet Weaver: internal volume of 380 μL , used for applications that need best mixing

Figure 2 shows the results of a blank (mobile phase A) injection with a 40-minute gradient carried out with the various Jet Weaver mixers. It is clear that despite the excellent behavior of the V35 Jet Weaver for most analyses, its performance for the detection of low-level compounds under these particular conditions is poor. Increasing the mixing volume to 100 μL significantly improves baseline stability, and by installing the 380 μL mixer, the noise caused by TFA is nearly eliminated. It is striking how the small system peaks present in the blank analysis (retention time 23 to 27 minutes) are not detectable with the V35 Jet Weaver whereas they are easily visible with the V380 mixer. Therefore, the V380 High-Performance Jet Weaver was selected for further analyses.


Figure 2. Analysis of a blank solution with an Agilent V35 (green), Agilent V100 (red), and Agilent V380 (blue) Jet Weaver mixer, 40-minute gradient, UV detection at 214 nm.

Switching between a low-volume and high-volume mixer will change the delay volume of the setup, and will affect retention times and potentially selectivity. To maintain the separation, some delay time should be added to the gradient to correct for this. This time shift for the gradient enabled nearly identical retention times and selectivity to be obtained with all three mixers (Figure 3).

Considerations on peak capacity

As shown in Figure 3, the tryptic digests are complex samples, and high peak capacity is mandatory for adequately characterizing mAbs and ADCs. With the current system and column configuration, the peak capacity can be adjusted according to the desired performance by changing only the gradient steepness and gradient time⁴.

As illustrated, the Kadcylya tryptic digest was analyzed with 12 different gradient times between 5 and 200 minutes. Four peptides were selected to calculate the peak capacity at 4 sigma (= 13.4 %) peak height. Figure 4 and Table 1 show the results. Short gradients with peak capacity below 250 can be used for fast (high productivity) analyses, whereas long gradients will result in peak capacities close to 900. The Herceptin tryptic digest was injected with some selected gradient conditions as a control, and peak capacities were in accordance to expectations. Figure 5 shows the results for these analyses. Note that the slope of the curve flattens strongly from 120 minutes onward, and that, in fact, working at slower slopes is pointless.


Figure 3. Analysis of Kadcylya tryptic digest with an Agilent V35 (green), Agilent V100 (red), and Agilent V380 (blue) Jet Weaver mixer, 40-minute gradient, UV detection at 214 nm.


Figure 4. Peak capacity as a function of gradient time. UV detection at 214 nm (detector set at 40 Hz to ensure sufficient data points in fast analyses).

Table 1. Peak capacity according to gradient for Kadcylya and Herceptin.

Gradient time (min)	Peak capacity	
	Kadcylya	Herceptin
5	95	134
10	177	
15	238	
20	280	342
30	359	
40	438	475
50	506	
60	596	
80	635	646
120	822	
160	836	
200	893	891

Further comments on the analysis of the Herceptin and Kadcylla digests

A 40-minute gradient resulted in a peak capacity of approximately 450, and these conditions, a good compromise between analysis time and separation performance, were used to highlight the difference between the mAb and the ADC. The protein sequence of Kadcylla is identical to Herceptin; the difference is in the conjugation of the cytotoxic agent emtansine to lysine residues. Figure 6 shows an overlay of both samples, with detection at 214 nm and 252 nm. Overall, the chromatograms are similar except for the cluster of peaks eluting late in the gradient (between 25 and 35 minutes). These are peptide-drug conjugates all containing emtansine. The complexity stems from the fact that a high number of lysine residues are available for conjugation. Figure 7 shows a detail of the chromatograms recorded at 252 nm. The addition of the drug to the peptide increases hydrophobicity and, therefore, retention. Since this cytotoxic agent has a UV absorbance at 252 nm, the conjugates are better observed with the detector set at this wavelength. The repeatability of the developed method was evaluated by five replicate injections of the Kadcylla tryptic digest. The overlay in Figure 8 shows that the injection and retention time precision are excellent.


Figure 5. Analysis of a Herceptin tryptic digest with a 20, 40, 80, and 120-minute gradient, UV detection at 214 nm. Note, for comparison, that the y-scale was increased with the Δ -gradient factor.


Figure 6. Comparison of analysis of a Herceptin (blue) and Kadcylla (red) tryptic digest with a 40-minute gradient, UV detection 214 nm and 252 nm.

Conclusions

High peak capacities can be achieved for the analysis of tryptic digests of mAbs and ADCs with an Agilent 1290 Infinity II LC in combination with an Agilent AdvanceBio Peptide Map column. The Agilent V380 Jet Weaver mobile phase mixer is effective in reducing the noise caused by the TFA modifier. This opens opportunities to quantify low abundant species in complex mixtures. The excellent precision obtained with the current setup suggests that this is a valuable tool for detailed analysis of protein biopharmaceuticals such as mAbs and ADCs. A next level of detail can be obtained when analyzing these digests on an LC×LC setup as described recently^{5,6}.

References

1. Sandra, K.; Vandenhede, I.; Sandra, P. Modern chromatographic and mass spectrometric techniques for protein biopharmaceutical characterization. *J. Chromatogr. A* **2014**, *1335*, 81–103.
2. Fekete, S.; *et al.* Chromatographic, electrophoretic, and mass spectrometric methods for the analytical characterization of protein biopharmaceuticals. *Anal. Chem.* **2016**, *88*, 480–507.
3. Huesgen, A. G. Performance of the Agilent 1290 Infinity Quaternary Pump using Trifluoroacetic Acid (TFA) as the mobile phase modifier, *Agilent Technologies Technical Overview*, publication number 5990-4031EN, **2009**.
4. Vanhoenacker, G.; *et al.* Tryptic digest analysis using the Agilent 1290 Infinity LC System, *Agilent Technologies Application Note*, publication number 5990-4031EN, **2009**.
5. Vanhoenacker, G.; *et al.* Comprehensive two-dimensional liquid chromatography of therapeutic monoclonal antibody digests, *Anal. Bioanal. Chem.* **2015**, *407*, 355-366.
6. Sandra, K.; *et al.* Multiple heart-cutting and comprehensive two-dimensional liquid chromatography hyphenated to mass spectrometry for the characterization of the antibody-drug conjugate ado-trastuzumab emtansine, *J. Chromatogr. B*, doi: 10.1016/j.jchromb.2016.04.040.


Figure 7. Detail of the comparison of an analysis of a Herceptin and Kadcyra tryptic digest with a 40-minute gradient, UV detection at 252 nm.


Figure 8. Overlay of five replicate injections of a Kadcyra tryptic digest with a 40-minute gradient, UV detection at 214 nm.

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 **Agilent**
 Trusted Answers

Determination of Amino Acid Composition of Cell Culture Media and Protein Hydrolysate Standard

The Agilent AdvanceBio Amino Acid Solution

Authors

M. Sundaram Palaniswam
Agilent Technologies, Ltd

Abstract

This study presents a method for analyzing primary amino acids in cell culture media using the Agilent AdvanceBio Amino Acid Analysis (AAA) solution with absorbance detection. Derivatization using an online injector program with OPA and FMOC decreases sample preparation time, and increases reproducibility over traditional offline methods. The effectiveness of this solution for routine analysis was confirmed using a system suitability test and retention time and area precision studies. The AdvanceBio AAA solution provides sensitive and high-resolution separation of all amino acids in cell culture media. The limit of detection (LOD), limit of quantification (LOQ), and linearity for selected amino acids for qualitative assays are also reported.

Introduction

Amino acids are the basic building blocks of proteins. They constitute all proteinaceous material of the cell including the cytoskeleton and the protein component of enzymes, receptors, and signaling molecules. In addition, amino acids are used for the growth and maintenance of cells. Cell culture media plays a key role in the biopharma industry. A large proportion of the amino acids supplied from cell culture media are diverted to pathways that could influence the fate of the cells in a culture. The identification of the optimal concentration of amino acids is important in fed batch and perfusion culture. Therefore, the design of an amino acid supplementation strategy might be streamlined by identifying the amino acid demands of a cell culture due to host cell growth and product production.

HPLC with precolumn derivatization is commonly used for the analysis of amino acids. Precolumn derivatization of free amino acids in solution for HPLC separations with UV or fluorescence detection is sometimes done manually, offline. Some immediate drawbacks to offline derivatization are sources of error due to operator skill, competence, and laboratory technique. Other drawbacks include extra sample manipulation, extra time required, and increased risk of contamination. Automated online derivatization minimizes these error sources, immediately improves precision, and saves time. Thus, a rugged high-resolution HPLC method including online derivatization, can increase productivity compared to offline methods. Consistent automated OPA derivatization, using the injector programming of the HPLC's autosampler and highly efficient Agilent AdvanceBio AAA columns, generate a rapid-reproducible amino acid method ideal for cell culture media. This method is convenient because the cell media samples are simply transferred to autosampler vials and analyzed. The selectivity of the AdvanceBio AAA column and the mobile phase gradient provides high resolution of 23 amino acids.

Materials and Methods

Instrumentation

Analyses were performed using an Agilent 1290 Infinity LC, which was equipped with an Agilent 1290 Infinity binary pump delivery system (G4220A), Agilent 1290 Infinity autosampler (G4226A), Agilent 1290 Infinity thermostatted column compartment (G1316C), and Agilent 1290 Infinity DAD (G4212A).

Reagents, samples, and materials

Cell culture media for compositional analysis, Minimum Essential Medium Eagle (M4655), Non-Essential Amino Acid (M7145), RPMI 1640 (R0083), Na₂HPO₄, and Na₂B₄O₇·10H₂O, were bought from Sigma-Aldrich. Protein hydrolysate was obtained from Fisher Scientific. Acetonitrile and methanol used were bought from Lab-Scan (Bangkok, Thailand). HPLC grade and highly purified water from a Milli-Q water purification system (Millipore Elix 10 model, USA) was used.

Column

Agilent AdvanceBio AAA, C18, 4.6 × 100 mm, 2.7 μm (p/n 655950-802)

Preparation of HPLC mobile phase

Mobile phase A contained 10 mM Na₂HPO₄ and 10 mM Na₂B₄O₇, pH 8.2. Mobile phase B contained acetonitrile, methanol, and water (45:45:10, v:v:v). Since mobile phase A is consumed at a faster rate than B, it is convenient to make 2 L of A for every 1 L of B produced.

Injection diluent

The injection diluent was 100 mL mobile phase A, plus 0.4 mL concentrated H₃PO₄ in a 100 mL bottle, stored at 4 °C. To prepare 0.1 N HCl, 4.2 mL concentrated HCl (36%) was added to a 500 mL volumetric flask that was partially filled with water, mixed, then filled to the mark with water. This solution was then used for making extended amino acid and internal standard stock solutions. It was stored at 4 °C.

Agilent AdvanceBio AAA standards and reagents kit, p/n 5190-9426, includes:

Part number	Component
5061-3339	Borate buffer: 0.4 M in water, pH 10.2, 100 mL
5061-3337	FMOC reagent, 2.5 mg/mL in ACN, 10 × 1 mL ampules
5061-3335	OPA reagent, 10 mg/mL in 0.4 M borate buffer and 3-mercaptopropionic acid, 6 × 1 mL ampules
5062-2479	Dithiodipropionic acid (DTDPA) reagent, 5 g
5061-3330	AA standard, 1 nmol/μL, 10 × 1 mL
5061-3331	AA standard 250 pmol, 10/pk
5061-3332	AA standard, 100 pmol/μL, 10 × 1 mL
5061-3333	AA standard, 25 pmol/μL, 10 × 1 mL
5061-3334	AA standard, 10 pmol/μL, 10 × 1 mL
5062-2478	Amino acids supplement kit, 1 g each

Derivatization reagents

Derivatization reagents (borate buffers, OPA, and FMOC) are ready-made solutions supplied by Agilent, and were transferred from their containers to autosampler vials. Precautions included:

- OPA is shipped in ampules under inert gas to prevent oxidation. Once opened, the OPA is potent for approximately 7 to 10 days. Therefore, 100 μL aliquots of OPA were transferred in microvial inserts and refrigerated. The OPA autosampler microvial was then replaced daily. Each ampule lasted 10 days (one vial/day).
- FMOC is stable in dry air, but deteriorates in moisture. Therefore, FMOC was transferred in 100 μL aliquots to microvial inserts and refrigerated. An open FMOC ampule transferred to 10 microvial inserts should last 10 days.
- Borate buffer was transferred to a 1.5 mL autosampler vial without a vial insert, and replaced every 3 days.

Preparation of amino acid standards

- Solutions of 17 amino acids in five concentrations are available from Agilent (10 $\text{pmol}/\mu\text{L}$ to 1 $\text{nmol}/\mu\text{L}$) for calibration curves. Each 1 mL ampule of standards was divided into 100 μL portions in conical vial inserts, and stored at 4 $^{\circ}\text{C}$.
- The extended amino acid (EAA) stock solution was produced by weighing 59.45 mg asparagine, 59.00 mg hydroxyproline, 65.77 mg glutamine, and 91.95 mg tryptophan into a 25 mL volumetric flask. This flask was filled halfway with 0.1 N HCL, and shaken or sonicated until the amino acids were dissolved. It was then filled to mark with water to produce a total concentration of 18 $\text{nmol}/\mu\text{L}$ of each amino acid.
- For the high-sensitivity EAA stock solution, 5 mL of this standard-sensitivity solution was diluted with 45 mL water (1.8 $\text{nmol}/\mu\text{L}$). Solutions containing extended standards were unstable at room temperature, and were kept frozen, and discarded at the first signs of reduced intensity.

Internal Standard (ISTD) stock solution

For primary amino acid ISTD stock solutions, 58.58 mg norvaline was weighed into a 50 mL volumetric flask. For secondary amino acids, 44.54 mg sarcosine was weighed into the same 50 mL flask. This flask was filled halfway with 0.1 N HCL, and shaken or sonicated until dissolved, then filled to mark with water for a final concentration of 10 nmol each amino acid/ μL (standard sensitivity). For high-sensitivity ISTD stock solution, 5 mL of standard-sensitivity solution was diluted with 45 mL of water, and stored at 4 $^{\circ}\text{C}$.

Calibration curves may be made using two to five standards, depending on experimental need. Typically, 100 $\text{pmol}/\mu\text{L}$, 250 $\text{pmol}/\mu\text{L}$, and 1 $\text{nmol}/\mu\text{L}$ are used in a three-point calibration curve for standard analytical sensitivity analysis. The following tables should be followed if an internal standard or other amino acids (for example, the extended amino acids) are added. Table 1 describes standard analytical sensitivity concentrations typically used in UV analysis.

Table 1. Chromatographic parameters used for intact and reduced analysis.

	Concentration of Final AA solution ($\text{pmol}/\mu\text{L}$)		
	900	225	90
Take 5 mL of 18 nmol EAA	5 mL	5 mL	5 mL
Dilute with water	–	15 mL	45 mL
Diluted EAA mix	5 mL	20 mL	50 mL
Take 5 mL of diluted EAA mix	5 mL	5 mL	5 mL
Add 10 nmol ISTD solution	5 mL	5 mL	5 mL
EAA-ISTD mix	10 mL	10 mL	10 mL
Take 100 μL of EAA-ISTD mix	100 μL	100 μL	100 μL
For 1 nmol AA, add:	900 μL	–	–
For 250 pmol AA, add:	900 μL	–	–
For 100 pmol AA, add:	–	–	900 μL
Final AA solution with EAA and 500 $\text{pmol}/\mu\text{L}$ ISTD	1 mL	1 mL	1 mL

Online derivatization

Depending on the autosampler model, the automated online derivatization program differs slightly. For the Agilent G4226A well plate automatic liquid sampler (WPALS), the injection program was:

1. Draw 2.5 μL from borate vial (p/n 5061-3339).
2. Draw 1.0 μL from sample vial.
3. Mix 3.5 μL in wash port five times.
4. Wait 0.2 minutes.
5. Draw 0.5 μL from OPA vial (p/n 5061-3335).
6. Mix 4.0 μL in wash port 10 times default speed.
7. Draw 0.4 μL from FMOC vial (p/n 5061-3337).
8. Mix 4.4 μL in wash port 10 times default speed.
9. Draw 32 μL from injection diluent vial.
10. Mix 20 μL in wash port eight times.
11. Inject.
12. Wait 0.1 minutes.
13. Valve bypass.

The location of the derivatization reagents and samples is up to the analyst and the ALS tray configuration. Using the G1367C with a 2 \times 56-well plate tray (p/n G2258-44502), the locations were:

- Vial 1: Borate buffer
- Vial 2: OPA
- Vial 3: FMOC
- Vial 4: Injection diluent
- P1-A-1: Sample

Thermostatted column compartment (TCC)

Left and right temperatures were set to 40 $^{\circ}\text{C}$. Analysis was enabled when the temperature was within ± 0.8 $^{\circ}\text{C}$.

Diode array detector (DAD)

Signal A: 338 nm, 10 nm bandwidth, reference wavelength 390 nm, 20 nm bandwidth.

Signal B: 262 nm, 16 nm bandwidth, reference wavelength 324 nm, 8 nm bandwidth.

Signal C: 338 nm, 10 nm bandwidth, reference wavelength 390 nm, 20 nm bandwidth.

To detect both OPA and FMOC derivatized amino acids in a single chromatogram, it was necessary to switch detector wavelengths. This switch took place between the last eluting OPA derivatized amino acid, lysine (peak 20 in the standard), and the first eluting FMOC derivatized amino acid, hydroxyproline (peak 21 in the standard).

With the DAD, determining the appropriate transition point was possible by initially collecting two channels. Signal A, 338 nm, detected OPA derivatized amino acids, and signal B, 262 nm, detected FMOC derivatized amino acids. From this analysis, the ideal point at which to switch wavelength during the run was determined. Subsequent runs were then made using a single channel, with the detector timetable function used to program a wavelength switch from 338 to 262 nm at the appropriate time between the elution of OPA-lysine and FMOC-hydroxyproline. This switch allowed both OPA and FMOC-derivatized amino acids to be detected in a single chromatogram. Peak width settings of >0.01 minutes were used for all columns.:

Linearity, limit of detection (LOD), and limit of quantification (LOQ) determination

As an example extended amino acid (EAA) stock solution, asparagine (59.45 mg), glutamine (65.77 mg), and tryptophan (91.95 mg) were used for linearity, LOD, and LOQ determination. These standards were weighed into a 25 mL volumetric flask, which was filled halfway with 0.1 N HCl, and mixed or sonicated until they dissolved. The flask was then filled to mark with water for a total concentration of 18 nmol/ μL of each amino acid.

Linearity was studied in the range of 0.9–1,000 pmol/ μL of these standard amino acids. Appropriate AA standard solutions were prepared in triplicate and injected into the chromatograph. The LOD and LOQ were estimated from the calibration function. LOD and LOQ were calculated as $3 (SD(a)/b)$ and $10 (SD(a)/b)$, respectively, where $SD(a)$ is the standard deviation of the intercept, and b is the slope of the calibration function.

Gradient program	
Time (min)	%B
0	2
0.35	2
13.4	57
13.5	100
15.7	100
15.8	2
18	end

Flow rate: 1.5 mL/min for 4.6 mm id

Table 3. Retention time and area RSD precision for amino acids (1,000 pmol) separated on an Agilent AdvanceBio AAA, 4.6 × 100 mm, column (n = 6).

Amino acid	RT RSD (%)	Area RSD (%)	Amino acid	RT RSD (%)	Area RSD (%)
1. Aspartic acid	0.151	2.60	13. Valine	0.025	2.4
2. Glutamic acid	0.512	2.19	14. Methionine	0.025	1.78
3. Asparagine	0.124	2.13	15. Norvaline	0.019	1.77
4. Serine	0.114	1.74	16. Tryptophan	0.024	2.03
5. Glutamine	0.092	1.8	17. Phenylalanine	0.027	1.98
6. Histidine	0.077	1.39	18. Isoleucine	0.025	2.17
7. Glycine	0.068	1.48	19. Leucine	0.020	1.81
8. Threonine	0.059	2.26	20. Lysine	0.022	2
9. Arginine	0.027	1.66	21. Hydroxyproline	0.014	3.14
10. Alanine	0.031	1.87	22. Sarcosine	0.015	5.01
11. Tyrosine	0.034	2.04	23. Proline	0.011	2.58
12. Cysteine	0.030	2.22			

System suitability as per the European Pharmacopoeia (Ph. Eur.)

The European Pharmacopoeia (Ph. Eur.) defines requirements for the qualitative and quantitative composition of amino acids and mixtures of amino acids. The requirements for allowed impurities are also defined. Manufacturers of amino acids are legally bound to prove that their amino acids meet these specifications before they can distribute their products in Europe. Leucine (Leu) is a branched-chain α -amino acid that is produced by a fermentation process (Figure 2). During this process, isoleucine can be produced as a by-product. The Ph. Eur. states that leucine and isoleucine should have a resolution of not less than 1.5 [1].

Ten concentration points for three amino acids were assayed in triplicate. The three standard amino acids showed good linearity in the tested range. The area response obeyed the equation $y = mx + C$, where the intercept C was zero within 95 % confidence limits, and the square correlation coefficient (R^2) was always greater than 0.99. Figure 3 shows the linearity curve for asparagine, glutamine, and tryptophan in the concentration range evaluated.

The LOD and LOQ were approximately 0.9 pmol and 3.8 pmol, respectively, using UV detection, indicating that the method was sensitive. Table 5 shows the observed LOD and LOQ values of asparagine, glutamine, and tryptophan.


Figure 2. Isoleucine and leucine chemical relationship.

Table 4. System suitability testing using AdvanceBio AAA columns and AA standards.

System suitability	Agilent AdvanceBio AAA, C18, 4.6 × 100 mm, 2.7 μ m	Agilent AdvanceBio AAA, C18, 3 × 100 mm, 2.7 μ m
Resolution between leucine and isoleucine (≥ 1.5)	4.5	4.6

Table 5. LODs and LOQs for three amino acids.

Asparagine		Glutamine		Tryptophan	
Concentration (pmol)	S/N ratio	Concentration (pmol)	S/N ratio	Concentration (pmol)	S/N ratio
0.9 (LOD)	5.3	0.9 (LOD)	3.0	0.9 (LOD)	4.5
1.9 (LOQ)	10.8	3.8 (LOQ)	13.8	3.8 (LOQ)	20.5

Results and Discussion

High-throughput routine analysis

The chromatogram in Figure 1 illustrates the standard analytical sensitivity achieved in high-throughput separations of amino acids. This chromatogram was obtained using an Agilent 1290 Infinity LC with an Agilent AdvanceBio AAA, 4.6 × 100 mm, 2.7 μm column using the amino acid method with DAD detection. A single run was completed in 18 minutes (including re-equilibration) with adequate resolution. The primary amino acids (1-20, OPA-derivatized), shown in Figure 1, were monitored at 338 nm, while the secondary amino acids (21-23, Fmoc-derivatized) were monitored at 262 nm.

Precision of retention time and area (n = 6)

Tables 2 and 3 summarize the average retention times and area RSDs for all the amino acids for the 100 and 1,000 pmol from six replicates of an amino acid method. The retention time RSDs for all amino acid peaks, including the early eluting peak 1 were less than 1.2%, demonstrating excellent gradient reproducibility. Peak area RSDs were less than 5 %, indicating precise sample injection. The RSD values demonstrate the robustness and precision of the amino acid method.


Figure 1. Separation of amino acid standard (1 nmol) on an Agilent AdvanceBio AAA 4.6 × 100 mm, 2.7 μm column, using the amino acid method.

Table 2. Retention time and area RSD precision for amino acids (100 pmol) separated on an Agilent AdvanceBio AAA, 4.6 × 100 mm, column (n = 6).

Amino acid	RT RSD (%)	Area RSD (%)	Amino acid	RT RSD (%)	Area RSD (%)
1. Aspartic acid	1.270	1.066	13. Valine	0.084	2.47
2. Glutamic acid	0.973	1.85	14. Methionine	0.073	1.82
3. Asparagine	0.605	1.79	15. Norvaline	0.073	1.72
4. Serine	0.629	1.82	16. Tryptophan	0.054	1.57
5. Glutamine	0.470	1.56	17. Phenylalanine	0.051	1.66
6. Histidine	0.430	1.22	18. Isoleucine	0.047	1.72
7. Glycine	0.477	1.92	19. Leucine	0.03	1.7
8. Threonine	0.440	1.95	20. Lysine	0.028	1.66
9. Arginine	0.251	2.15	21. Hydroxyproline	0.021	4.13
10. Alanine	0.280	3.06	22. Sarcosine	0.026	1.15
11. Tyrosine	0.128	1.65	23. Proline	0.021	4.36
12. Cysteine	0.067	1.9			

Amino acid analysis of cell culture media and protein hydrolysate standard

We analyzed the amino acid composition of commonly used cell culture supplements. These standards included: Minimum Essential Medium Eagle (MEM), Non-Essential Amino Acid (NEAA), RPMI 1640 R0083, and protein hydrolysate standard. The results were then compared with the amino acid standards. Figures 4 to 7 show the overlays of amino acid composition of the media and the amino acid standards.

It is evident that the amino acid composition of cell culture supplements, as determined by this method, matches accurately with their theoretical composition. In addition, baseline resolution of isoleucine and leucine was observed with a resolution factor of 4.35 for the protein hydrolysate standard, meeting the regulatory requirements for these components significantly better than competitive columns. Such applications are useful in monitoring and adjusting amino acid composition, which is an essential part of optimizing the manufacturing process to ensure high quality and optimum yield of the final biopharmaceutical product.


Figure 3. Linearity curve with 10 standard concentrations of asparagine, glutamine, and tryptophan ranging from 0.9 to 1,000 pmol, showing excellent coefficient values.


Figure 4. Amino acid analysis of MEM media (blue trace) mirrored with AA standards (red trace).


Figure 5. Amino acid analysis of NEAA media (blue trace) and comparison with AA standards (red trace).


Figure 6. Amino acid analysis of RPMI 1650 media (blue trace) and comparison with AA standards (red trace).


Figure 7. Amino acid analysis of protein hydrolysate standard. The resolution between leucine and isoleucine with the Agilent AdvanceBio AAA, 4.6 × 100 mm, 2.7 μm column is much higher than the reported value for system suitability requirement for this pair.

Conclusion


Amino acid analysis is an important approach for the characterization of protein and peptide-based products. Studying the roles of amino acids during bioprocesses leads to better understanding the feeding strategy, and to improving the yield and quality of the product. In addition, the determined amino acid composition can confirm sample identity, and give a measure of sample purity. This Application Note demonstrates several Agilent tools for the analysis of amino acids. We first used the Agilent 1290 Infinity LC and the Agilent AdvanceBio AAA kit for the automated online derivatization of amino acids using OPA/FMOC chemistries. The derivitized amino acids were then separated on an AdvanceBio AAA LC column to achieve a fast, sensitive, and reproducible separation of amino acids. Area and RT precision of the method were excellent, and met the system suitability requirement. Linearity curves with 10 standard concentrations of three amino acids, ranging from 0.9 pmol to 1 nmol, had excellent coefficient of linearity values, indicating that the method was quantitative and accurate. The LOD and LOQ for the amino acids were 0.9 pmol and 3.8 pmol, respectively, indicating that the method was sensitive. In addition, this method was able to separate and detect amino acids from cell culture media and protein hydrolysate standard. The amino acid composition determined using this method correlated well with their theoretical compositions..

Reference

1. European Pharmacopoeia 9.0 (2.2.56) *Amino Acid Analysis*

For More Information

These data represent typical results. For more information on our products and services, visit www.agilent.com/chem.

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1. Wang, X.; Barber, W. E.; Long, W. J. Applications of superficially porous particles: High Speed, high efficiency or both? *J. Chromatogr. A.* **2012**, 1228, 72-88.
2. Schuster, R.; Apfel, A. A new technique for the analysis of primary and secondary amino acids; Application note, Hewlett-Packard Publication number 5954-6257, **1986**.
3. Schuster, R. Determination of amino acids in biological, pharmaceutical, plant and food samples by automated precolumn derivatization and high-performance liquid chromatography. *J. Chromatogr. B.* **1988**, 431, 271-284.
4. Henderson, Jr., J. W.; Ricker, R. D.; Bidlingmeyer, B. A.; Woodward, C. Rapid, Accurate, Sensitive and Reproducible HPLC Analysis of Amino Acids; Application note, Agilent Technologies, Inc. Publication number 5980-1193E, **2000**.
5. Woodward, C.; Henderson, Jr., J. W.; Todd Wielgos, T. High-Speed Amino Acid Analysis (AAA) on Sub-Two Micron Reversed-phase (RP) Columns; Application note, Agilent Technologies, Inc. Publication number 5989-6297EN, **2007**.
6. Gratzfeld-Huesgen, A. Sensitive and Reliable Amino Acid Analysis in Protein Hydrolysates using the Agilent 1100 Series HPLC; Application note, Agilent Technologies, Inc. Publication number 5968-5658EN, **1999**.
7. Greene, J.; Henderson, Jr., J. W.; Wikswo, J. P. Rapid and Precise Determination of Cellular Amino Acid Flux Rates Using HPLC with Automated Derivatization with Absorbance Detection; Application note, Agilent Technologies, Inc. Publication number 5990-3283EN, **2009**.
8. Henderson, Jr., J. W.; Brooks, A. Improved Amino Acid Methods using Agilent ZORBAX Eclipse Plus C18 Columns for a Variety of Agilent LC Instrumentation and Separation Goals; Application note, Agilent Technologies, Inc. Publication number 5990-4547EN, **2010**.
9. Kirkland, J. J.; van Straten, M. A.; Claessens, H. A. Reversed-phase high-performance liquid chromatography of basic compounds at pH 11 with silica-based column packings. *J. Chromatogr. A.* **1998**, 797, 111-120.
10. Tindall, G. W.; Perry, R. L. Explanation for the enhanced dissolution of silica column packing in high pH phosphate and carbonate buffers. *J. Chromatogr. A.* **2003**, 988, 309-312.
11. Anon. Extending Column Lifetime in Pharmaceutical Methods with High pH-Stable Poroshell HPH Chemistries; Technical overview, Agilent Technologies, Inc. Publication number 5991-5022EN, **2014**.

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Fast and Efficient HILIC Methods for Improved Analysis of Complex Glycan Structures

Authors

James Martosella,
Oscar Potter,
Danny Mancheno,
and Jia Liu
Agilent Technologies, Inc

Introduction

N-linked glycosylation is a critically important and very complex post-translational modification. It therefore needs to be controlled and monitored throughout development, processing, and manufacture of drug glycoproteins. Therapeutic protein characteristics, including safety, efficacy, and serum half-life, can be affected by differences in their glycosylation pattern, and so the analysis of these patterns is an important part of the characterization of therapeutic glycoproteins, particularly mAbs. Separation using HILIC with fluorescence detection is a robust method for glycan analysis, whereas HILIC/LC can also be coupled to mass spectrometry to obtain important mass and structure information.

One of the growing challenges in HILIC/LC, however, is achieving high-resolution separations with fast analysis times. With ever-increasing demands placed on biocharacterization for higher throughput, researchers are looking for improved separation (HILIC) methods, but not at the cost of lost separation performance. Since glycans include many closely related structures, it is critical to achieve the highest resolution possible, and preferably during a fast analysis time.

In this work, we used a sub-2 μm UHPLC HILIC column with amide chemistry for high-throughput glycosylation profiling. Specifically, we profiled 2-AB labeled human IgG and bovine fetuin N-linked glycans using a 1.8 μm , 2.1 \times 150 mm column with fluorescence detection. Rapid, sensitive and selective separations were achieved to provide ultra-high resolution of these complex glycans in run times as short as 9 minutes. In a run time comparison to a currently available UHPLC glycan column, we observed a 40 % reduction in analysis speed for human IgG N-linked glycans under identical conditions.

Materials and Methods

Conditions, recombinant human IgG1

Parameter	Value
Sample:	Agilent 2-AB labelled IgG N-linked glycan library, 200 pmol (p/n5190-6996)
Column:	Agilent AdvanceBio Glycan Mapping, 2.1 × 150 mm, 1.8 μm (p/n859700-913)
Mobile phase:	A, 100 mM NH ₄ formate, pH 4.5; B, ACN
Inj vol:	2 μL in 70/30 ACN:water
Column temp:	55 °C
Sample thermostat:	10 °C
Detection:	Fluorescence, excitation 260 nm, emission 430 nm
Instrument:	Agilent 1290 Infinity LC with 1260 Fluorescence Detector

Conditions, bovine fetuin

Parameter	Value
Sample:	ProZyme Glyko 2-AB bovine fetuinN-linked library
Column:	Agilent AdvanceBio Glycan Mapping, 2.1 × 150 mm, 1.8 μm (p/n859700-913)
Mobile phase:	A, 100 mM NH ₄ formate, pH 4.5; B, ACN
Inj vol:	5 μL (20 pmol) in 70/30 ACN:water
Column temp:	55 °C
Other conditions as above	

Conditions, dextran ladder

Parameter	Value
Sample:	Agilent 2-AB labelled dextran ladder standard (p/n5190-6998)
Column:	Agilent AdvanceBio Glycan Mapping, 2.1 × 150 mm, 1.8 μm (p/n859700-913)
Mobile phase:	A, 100 mM NH ₄ formate, pH 4.5; B, ACN
Inj vol:	2 L 75:25 ACN:water (10 pmol total glycan)
Other conditions as above	

Conditions, mass spectrometry, recombinant human IgG1

Parameter	Value
Instrument:	Agilent 6550 iFunnel Q-TOF LC/MS
Source:	Agilent Dual JetStream
Drying gas temp:	200 °C
Drying gas flow:	12 L/min
Sheath gas temp:	250 °C
Sheath gas flow:	12 L/min
Nebulizer pressure:	25 psi
Capillary voltage:	3,500 V
Nozzle voltage:	500 V
Fragmentor voltage:	250 V
Skimmer voltage:	45 V
Octupole RF Vpp:	550 V
MS range:	100 to 1,700 m/z
MS scan rate:	2 spectra/s

The workflow is shown in Figure 1.


Figure 1. Amino acid analysis of protein hydrolysate. The resolution between leucine and Isoleucine with the AdvanceBio AAA, 4.6 x 100 mm, 2.7 μm column is much higher than the reported value for system suitability requirement.

Results and Discussion

Dextran ladder

Figure 2 shows the separation of a 2-AB labeled dextran ladder. In this separation, a homopolymeric series of 20 glucose oligomers were efficiently baseline resolved in less than 15 minutes.


Figure 2. The Agilent AdvanceBio Glycan Mapping column separates a 2-AB labeled dextran ladder to baseline in less than 15 minutes

Human IgG glycans

Figure 3 demonstrates an ultra-high resolution separation of 2-AB labeled N-linked human IgG glycans.


Figure 3. Very high resolution of 2-AB labeled N-linked human IgG glycans on the Agilent AdvanceBio Glycan Mapping column.

Fast separation in less than 10 minutes was also achieved, as shown in Figure 4.


Figure 4. The Agilent AdvanceBio Glycan Mapping column separates 2-AB labeled N-linked human IgG glycans in less than 10 minutes.

MS N-linked glycans

Figure 5 shows eight representative mass spectra from the ultra-high resolution separation of human IgG glycans. The spectra were generated by Q-TOF analysis (experimental). All spectra matched theoretical masses to within 6 ppm.


Figure 5. Eight representative mass spectra from the ultra-high resolution separation of human IgG glycans (fluorescence detector separations are displayed on the panel to the left).

Sub-2 µm HILIC comparison

Figure 6 shows the results of a comparison of glycan amide columns. Using the same chromatographic conditions, the AdvanceBio Glycan Mapping column delivered better resolution and narrower bands, with higher peak capacity, at a 40 % faster separation time than another brand of sub-2 µm HILIC column in a 2.1 × 150 mm configuration.


Figure 6. The AdvanceBio Glycan Mapping column delivers better resolution and narrower bands, with higher peak capacity at a 40 % faster separation time.

Fetuin glycans

Finally, we revealed the fast and highly efficient performance of the AdvanceBio Glycan Mapping column in a separation of 2-AB labeled bovine fetuin N-linked glycans (Figure 7).



Figure 7. Fast and highly efficient separation of 2-AB labeled bovine fetuin N-linked glycans, 2.1 × 150 mm 1.8 µm AdvanceBio Glycan Mapping column.

Table 1. Glycan nomenclature and structure assignments for the optimized rapid separation of bovine fetuin 2-AB labeled N-linked glycans displayed in top chromatogram of Figure 7.

Peak	Retention	GU value	Glycan structure	Structures
1	6.70	9.4	A2G2S2	
2	7.27	9.8	A2G2S2	
3	8.73	10.8	A3G3S3, A3G3S2 (trace)	
4	9.22	11.2	A3G3S3, A3G3S2 (trace)	
5	9.70	11.6	A3G3S3, A3G3S4 (trace)	
6	10.20	12	A3G3S4, A3G3S3	
7	10.63	12.4	A3G3S4	

Conclusions

The Agilent 1.8 µm HILIC amide AdvanceBio Glycan Mapping column provided separation of N-linked glycans with high speed, excellent resolution and increased efficiency. In a 2.1 × 150 mm configuration and under identical chromatographic conditions, the column enabled a well-resolved separation of 2-AB labeled IgG N-linked glycans, with a 40 % reduction in elution time compared to another brand of sub-2 µm HILIC column. A separation of 2-AB labeled bovine fetuin N-linked glycans demonstrated the column's excellent analytical selectivity and resolving power for separating these complex biantennary and triantennary glycans.

Acknowledgement

This work was presented in a poster session at the 62nd ASMS Conference on Mass Spectrometry and Allied Topics, 15–19 June, 2014, Baltimore.

N-Glycan Analysis of mAbs and Other Glycoproteins with UHPLC and Fluorescence Detection


The Agilent 1290 Infinity Binary LC System with the Agilent AdvanceBio Glycan Mapping Column

Authors

Sonja Schneider and
Oscar Potter
Agilent Technologies, Inc.
Waldbronn, Germany

Abstract

This Application Note describes the analysis of N-linked glycans with hydrophilic interaction chromatography (HILIC) using the Agilent 1290 Infinity Binary LC together with the Agilent 1260 Infinity Fluorescence Detector and the Agilent 6530 Accurate-Mass Q-TOF LC/MS. Enzymatic glycan release with PNGase F followed by derivatization with 2-aminobenzamide (2-AB) was conducted on monoclonal antibodies (mAbs) and two other glycoproteins, fetuin and ovalbumin. The excellent resolution provided by the Agilent AdvanceBio Glycan Mapping column allowed detection and identification of all major N-glycans in the mAb sample. Furthermore, the highly complex N-glycans released from fetuin and ovalbumin were well resolved.

Introduction

Glycosylation is one of the most frequently observed post translational modifications. Mammalian glycoproteins contain three major types of glycans: N-linked, O-linked, and glycosylphosphatidylinositol (GPI) lipid anchors, which consist of one or more monosaccharide units. A single glycosylation site can generate considerable heterogeneity of the mass and charge of glycoproteins. These oligosaccharides are involved in many biological regulation and recognition processes, for example, protein sorting, immune and receptor recognition, inflammation, pathogenicity, metastasis, and other cellular processes^{1,2}. In addition, properties such as safety, efficacy, and the serum half-life of therapeutic proteins can be affected by their glycosylation pattern.

Recombinant monoclonal antibody therapeutics (mAbs) represent the largest group of therapeutic proteins. The efficacy of these therapeutics is highly dependent on the correct glycosylation pattern of the mAbs and, so far, all licensed therapeutic mAbs are immunoglobulins G (IgGs)³. Human IgG has a single conserved N-linked glycosylation site located on the Fc region of each heavy chain at Asn297 (Figure 1), resulting in the presence of two N-glycans per IgG. This typically consists of a handful of major structures and numerous minor structures⁵. The combination of glycans at each of the two glycosylation sites on the Fc region leads to large numbers of different glycoforms in each batch of mAb production.

The glycan structure plays a critical role in complement activation and receptor affinity⁶, which affect the efficacy of therapeutic mAbs. Moreover, non-human glycans are a safety issue due to induced immune responses. Therefore, analysis of the glycan pattern is an important part of the characterization of therapeutic glycoproteins, especially mAbs.

This Application Note uses symbolic glycan structures according to the Consortium for Functional Glycomics (CFG), as shown in Figure 2. Assigned glycans are also described by the Oxford glycan nomenclature and by another style of nomenclature, which is popular for mAb glycans, shown here in italics.

Figure 2A shows the general nomenclature used to describe sugar residues of different glycan structures on proteins. Figure 2B shows some predominant glycan structures present on the Asn-297 site in human IgG. In general, N-glycans have a core structure, containing two b-D-N-acetylglucosamine (GlcNAc) and three mannose (Man) units. IgG Fc N-glycans are predominantly biantennary complex-type structures, partially core-fucosylated (for example, FA2 or G0F).


Figure 1. IgG antibody structure.


Figure 2. Glycan structure and isoforms. A) Monosaccharide description after the Consortium for Functional Glycomics, B) predominant glycan structures of human IgGs.

Different strategies for the analysis of N-glycans have been described. Many methods are based on enzymatic release of N-glycans from the protein by PNGase F. Due to the lack of intrinsic chromophores, it is also common to derivatize the glycans with a fluorescent label prior to analysis⁷. Each N-glycan contains a single reducing end site that can be reacted with an excess of fluorescent label, such that each N-glycan will be attached to one fluorophore. The processed sample is, therefore, appropriate for relative quantification by separation with fluorescence detection without the need for any quantitation standards or calibration. 2-AB is a stable, neutral label that is popular for N-glycan analysis^{7,8,9}. Figure 3 illustrates 2-AB labeling by reductive amination (Schiff's base intermediate not shown).

Subsequent purification using hydrophilic interaction chromatography/solid phase extraction (HILIC/SPE) is performed to remove the large excess of 2-AB so that it does not interfere with the HILIC/FLD analysis.

Here, we show enzymatic release of N-glycans using PNGase F with subsequent derivatization with 2-AB prior to separation by HILIC UHPLC, with fluorescence detection and identification by on-line quadrupole time-of-flight mass spectrometry (Q-TOF/MS).


Figure 3. Labeling of a glycan with 2-aminobenzamide (2-AB).

Experimental

Reagents, samples, and materials

The Agilent 1290 Infinity Binary LC System consisted of the following modules:

- Agilent 1290 Infinity Binary Pump (G4220A) with 35 μ L Jet Weaver
- Agilent 1290 Infinity Autosampler (G4226A)
- Agilent 1290 Infinity Thermostat (G1330B)
- Agilent 1290 Infinity Thermostatted Column Compartment (G1316C)

- Agilent 1260 Infinity Fluorescence Detector (G1321B) with standard flow cell

Reagents, samples, and materials

Agilent 6530 Accurate-Mass Q-TOF LC/MS

Column

Agilent AdvanceBio Glycan Mapping, 2.1 \times 150 mm, 1.8 μ m (p/n 859700-913)

Software

Agilent OpenLAB CDS ChemStation Edition for LC & LC/MS Systems Rev. C.01.05 [38] and Agilent MassHunter Workstation Software, Version B.05.01, Build 4.0.479.0. Glycan structures were created with GlycoWorkbench, Version 2.1, stable (146).

Sample preparations

Deglycosylation procedure:

N-glycans were released from a monoclonal antibody, fetuin, and ovalbumin using PNGase F. This enzyme cleaves asparagine-linked high mannose as well as hybrid and complex oligosaccharides from the glycoproteins and leaves the glycans intact. Fetuin has three N-glycosylation sites (Asn-81, Asn-138, and Asn-158) and four O-linked sites (Ser-253, Thr-262, Ser-264, and Ser-323)¹¹. Ovalbumin has only one glycosylation site, whereas the mAb contains two glycosylation sites. The amount of PNGase F was adjusted to the amount of N-glycosylation sites. The proteins were deglycosylated according to instructions for 3 hours at 37 °C. The reaction was then stopped, and the sample was vacuum-dried for further processing.

2-AB-labeling for fluorescence detection and sample cleanup

The dried glycan samples were labeled with 2-aminobenzamide according to the protocol for 3 hours at 65 °C. After the labeling procedure, the samples were purified using the HILIC cleanup cartridges according to the instruction manual. After the cleanup procedure, the samples were vacuum-dried and reconstituted in ultrapure water:acetonitrile 30:70 (v/v) for analysis.

Solvents and samples

Buffer A was 100 mM ammonium formate in water, pH 4.5 and buffer B was acetonitrile. All solvents were LC grade. Fresh ultrapure water was obtained from a Milli-Q Integral system equipped with a 0.22 μ m membrane point-of-use cartridge (Millipak). The monoclonal antibody was RAT Anti-DYKDDDDK Tag Antibody. Ammonium formate, fetuin and ovalbumin, PNGase F from Elizabethkingia miricola, GlycoProfil 2-AB Labeling Kit, and GlycoProfil Glycan Cleanup Cartridges were purchased from Sigma-Aldrich Corp., St. Louis, USA.

Instrumental conditions

	Antibody standard gradient	Fetuin gradient	Ovalbumin gradient
Starting flow rate	0.5 mL/min	0.5 mL/min	0.5 mL/min
Gradient	0 minutes 85 % B 5 minutes 75 % B 35 minutes 64 % B 40 minutes 50 % B 42 minutes, flow 0.5 mL/min 42.01 minutes, flow 0.25 mL/min 43 minutes 0 % B 48 minutes 0 % B 50 minutes 85 % B 50.01 minutes, flow 0.25 mL/min 51 minutes, flow 0.5 mL/min	0 minutes 75 % B 45 minutes 50 % B 47 minutes 40 % B, flow 0.5 mL/min 47.01 minutes, flow 0.25 mL/min 49 minutes 0 % B 51 minutes 0 % B 51.01 minutes 75 % B, flow 0.25 mL/min 52.00 minutes, flow 0.5 mL/min	0 to 6 minutes 85 % B 10 minutes 80 % B 60 minutes 70 % B 65 minutes 50 % B, flow 0.5 mL/min 65.01 minutes, flow 0.25 mL/min 68 minutes 0 % B 73 minutes 0 % B 74 minutes 85 % B, flow 0.25 mL/min 75.00 minutes, flow 0.5 mL/min
Stop time	51 minutes	52 minutes	75 minutes
Post time	20 minutes	20 minutes	20 minutes
Injection volume	5 µL	1 µL	1 µL
Thermostat autosampler	5 °C		
Column temperature	60 °C		
FLD	Ex. 260 nm, em. 430 nm		
Peak width	> 0.013 minutes (0.25 seconds resp. time) (37.04 Hz)		

MS parameters	
Gas temperature	250 °C
Sheath gas temperature	250 °C
Gas flow	8 L/min
Sheath gas flow	8 L/min
Nebulizer	25 psi
Vcap	3,500 V
Nozzle	1,000 V
Fragmentor	200 V
Skimmer	45 V
Oct 1 RF Vpp	550
Collision energies	15 and 30 V
Mode	MS and targeted MS/MS

Results and Discussion

Analysis of N-glycans from monoclonal antibodies

Figure 4 shows the separation of the mAb N-glycans. The mAb glycan pattern was optimally resolved, allowing separation and integration of all major N-glycans. Relative quantification was made based on the calculation of the peak area percentage. High intensity of the labeled glycans was achieved by setting the optimal wavelengths for glycan detection on the Agilent 1260 Infinity Fluorescence Detector, using 260 nm as excitation wavelength and 430 nm as emission wavelength¹⁰.

The resulting HILIC glycan profile was assigned to the corresponding glycan structures based on the parent ion mass observed and the related MS/MS spectra. The parent masses were entered into the GlycoMod tool from ExPASy to find related glycan structures.


Figure 4. Separation of mAb N-glycans with fluorescence detection with 260 nm as the excitation wavelength.

GlycoMod predicts possible glycan structures (labeled or unlabeled) from the experimentally determined masses (<http://web.expasy.org/glycomod/>). Another helpful tool for glycan assignment and glycan structure design is GlycoWorkbench15, which was used in this work to prepare glycan structure cartoons.

As an example of the workflow, the N-glycan FA2G1Sg1 with a parent mass of 1026.88 [z = 2] (peak 7 and 8) was chosen. The glycan databases revealed the two most likely glycan structures for this mass (Figure 5).

MS/MS data was then used to distinguish between these two potential structures. Figure 6 shows the collision-induced dissociation (CID) MS/MS spectrum of the N-glycan FA2G1Sg1. The MS/MS data confirm the presence of a type of sialic acid, N-glycolylneuraminic acid (NeuGc), which results in strong signals for fragment ions at m/z 308 (NeuGc) and m/z 673 (NeuGc attached to galactose and N-acetylglucosamine). Meanwhile, there are no signals at m/z 292 or m/z 657, which would have indicated the presence of N-acetylneuraminic acid (NeuAc). Therefore, these results provided strong evidence that the structure was FA2G1Sg1 (containing NeuGc) rather than A2G2S1 (which contains NeuAc). The decision to assign a structure with a core fucose was also supported by the MS/MS data, based on the lack of strong fragment ion signals at m/z 512 or m/z 350, which would have been present if the fucose had instead been attached in the outer arm region.


Figure 5. Two most likely glycan structures for the parent mass of 1,026.88 [z = 2], FA2G1Sg1 and A2G2S1.

All other peaks were similarly assigned using their MS and MS/MS spectra. The assigned structures are given in Table 1. The results show that the mAb mainly contains core fucosylated complex glycans, including several structures with NeuGc. These findings are typical for an IgG antibody produced by rat cells. NeuGc does not normally occur in human glycoproteins and is undesirable in therapeutic proteins¹². The sialic acids occurring in human glycoproteins are typically N-acetylneuraminic acids.


Figure 6. MS/MS spectra of FA2G1Sg1 - 1026.88 [z = 2] - 1931.6876 Da.

Table 1. Overview of masses and assigned 2AB- glycan structures of rat monoclonal antibody.

Peak	Oxford	Biopharma mAb style	Mass detector (Da)	Structure
1	FA1	G0-GlcNAc	A2G2S2	
2	FA2	G0F	792.3130 [z = 2]	
3	M5	Man5	1,355.5 [z = 1]	
4,5	FA2G1	G1F	873.34 [z = 2]	
6	FA1G1Sg1	G1FSg1-GlcNAc	925.34 [z = 2]	
7,8	FA2G1Sg1	G1FSg1	1,026.88 [z = 2]	
9	FA2G2Sg1	Ag1F	1,107.9135 [z = 2]	
10	FA2G2Sg2	Ag2F	1,261.446 [z = 2]	

Analysis of antibodies from fetuin and ovalbumin

N-glycans from two more proteins, fetuin and ovalbumin, were released by PNGase F, derivatized with 2-AB and analyzed using HILIC/UHPLC with online MS. Figure 7 shows the separation of bovine fetuin N-glycans. This glycosylation profile was dominated by complex non-fucosylated biantennary and triantennary glycans containing NeuAc. Nine major peaks could be assigned using Q-TOF/MS detection. Table 2 shows the assigned glycan structures.


Figure 7. Separation of 2-AB-labeled fetuin.

Table 2. Detailed information of N-glycan ovalbumins.

Peak	Oxford	Structure
1	A2G2S1	
2,3	A2G2S2	
4	A3G3S2	
5	A3G3S3, A3G3S2 (trace)	
6	A3G3S3, A3G3S2 (trace)	
7	A3G3S3, A3G3S4 (trace)	
8	FA3G3S4, A3G3S3	
9	A3G3S4	

GlycoMod Figure 8 shows the separation of ovalbumin glycans. Ovalbumin is N-glycosylated only at one site (Asn292), but a complex glycosylation pattern can be associated to this site¹³. Due to the complexity of the glycan profile, the gradient had to be adjusted to a longer separation time to achieve higher resolution. The high performance of the AdvanceBio Glycan Mapping column allowed resolution of over 50 peaks with a good signal-to-noise (S/N) ratio. Twenty major peaks were assigned based on the parent ion-mass data (Table 3). Detailed structural conclusions were not achievable due to the high chance of isobaric structures occurring, several of which cannot necessarily be distinguished from the MS/MS data. Instead, the N-glycans are described in terms of their monosaccharide composition.

Compared to the relatively simple glycan pattern of the mAb, the two other glycoproteins had a greater variety of glycan structures. No fucosylated glycans were detected in ovalbumin in contrast to the mAb glycans, which matches previously reported findings that avian egg glycoproteins are non-fucosylated¹⁴.


Figure 8. Separation of N-glycans released from ovalbumin.

Table 3. Assigned masses and monosaccharide composition of ovalbumin N-glycans; H = hexoses, i.e. galactose or mannose; N = N-acetylglucosamine.

Peak	Mass + 2AB (Da)	Calculated mass (Da)	PComposition (short form)	PComposition (long form)
1,2	1,234.48 [M]	1,114.48	H3N3	(HexNAc) ₁ +(Man) ₃ (GlcNAc) ₂
3	1,193.45 [M]	1,073.45	H4N2	(Hex) ₄ (HexNAc) ₂
4	1,437.56 [M]	1,317.56	H3N4	(HexNAc) ₂ +(Man) ₃ (GlcNAc) ₂
5	1,396.53 [M]	1,276.53	H4N3	(Hex) ₁ (HexNAc) ₁ +(Man) ₃ (GlcNAc) ₂
6, 7	1,640.64 [M]	1,520.64	H3N5	(HexNAc) ₃ +(Man) ₃ (GlcNAc) ₂
8	1,355.506 [M]	1,235.51	H3N5	(Hex) ₂ +(Man) ₃ (GlcNAc) ₂
9	1,599.61 [M]	1,479.61	H4N4	(Hex) ₁ (HexNAc) ₂ +(Man) ₃ (GlcNAc) ₂
10	1,843.73 [M]	1,723.73	H3N6	(HexNAc) ₄ +(Man) ₃ (GlcNAc) ₂
11	1,802.74 [M]	1,682.74	H4N5	(Hex) ₁ (HexNAc) ₃ +(Man) ₃ (GlcNAc) ₂
12, 14	2,046.884 [M]	1,926.88	H3N7	(HexNAc) ₅ +(Man) ₃ (GlcNAc) ₂
13	1,517.56 [M]	1,397.56	H6N2	(Hex) ₃ +(Man) ₃ (GlcNAc) ₂
15	1,761.6574 [M]	1,641.66	H5N4	(Hex) ₂ (HexNAc) ₂ +(Man) ₃ (GlcNAc) ₂
16	2,005.8098 [M]	1,885.81	H4N6	(Hex) ₁ (HexNAc) ₄ +(Man) ₃ (GlcNAc) ₂
17	2,249.9728 [M]	2,129.97	H3N8	(HexNAc) ₆ +(Man) ₃ (GlcNAc) ₂
18	1,964.82 [M]	1,844.82	H5N5	(Hex) ₂ (HexNAc) ₃ +(Man) ₃ (GlcNAc) ₂
19, 20	2,208.87 [M]	2,088.87	H4N7	(Hex) ₁ (HexNAc) ₅ +(Man) ₃ (GlcNAc) ₂

Conclusions

The Agilent 1290 Infinity Binary LC System, together with the Agilent 1260 Infinity Fluorescence Detector and Agilent 6530 Accurate-Mass Q-TOF/MS, was an ideal combination for the analysis of released N-glycans that were derivatized with 2-aminobenzamide. Sample preparation using PNGase F for the release of N-linked glycans followed by 2-AB derivatization with subsequent HILIC sample cleanup was demonstrated for one monoclonal antibody and two other glycoproteins.

The Agilent AdvanceBio Glycan Mapping column demonstrated excellent resolving power, allowing separation and identification of all major N-glycans in a rat mAb sample. Complex biantennary and triantennary N-glycans from fetuin and ovalbumin were also analyzed with very high resolution. Optimized fluorescence excitation and emission wavelengths of 260 and 430 nm provided better S/N ratios. Electrospray ionization Q-TOF MS analysis allowed assignment of different glycan structures or monosaccharide compositions, depending on the complexity of the sample.

References

1. Rademacher, T. W; Williams, P; DwekMark, R. A. Agalactosyl glycoforms of IgG autoantibodies are pathogenic. *P. Natl. Acad. Sci.* **1994**, 91, pp 6123-6127.
2. Peracaula, R; *et al.* Glycosylation of human pancreatic ribonuclease: differences between normal and tumor states. *Glycobiology* **2003**, 13, pp 227-244.
3. Jefferis, R. Glycosylation of recombinant antibody therapeutics. *Biotechnol. Progr.* 2005, 21, pp 11-16. 4. Arnold, J. N; *et al.* Human immunoglobulin glycosylation and the lectin pathway of complement activation. *Adv. Exp. Med. Biol.* **2005**, 564, pp 27-43.
4. Fernandes, D. Demonstrating comparability of antibody glycosylation during biomanufacturing. *Euro. Biopharm. Re.* **2005**, Summer, pp 106-110.
5. Abès, R; Teillaud, J. L. Impact of Glycosylation on Effector Functions of Therapeutic IgG. *Pharmaceuticals* **2010**, 3, pp 146-157.
6. Ruhaak, L. R; *et al.* Glycan labeling strategies and their use in identification and quantification. *Anal. Bioanal. Chem.* **2009**, 397, pp 3457-3481.
7. Royle, L; *et al.* HPLC-based analysis of serum N-glycans on a 96-well plate platform with dedicated database software. *Anal. Biochem.* **2008**, 376, pp 1-12.
8. Huhn, C; *et al.* IgG glycosylation analysis. *Proteomics* **2009**, 9, pp 882-913.
9. Melmer, M; *et al.* HILIC analysis of fluorescence-labeled N-glycans from recombinant biopharmaceuticals. *Anal. Bioanal. Chem.* **2010**, 398, pp 905-914.
10. Ding, W; *et al.* Identification and Quantification of Glycoproteins Using Ion-Pairing Normal-phase Liquid Chromatography and Mass Spectrometry. *Mol. Cell. Proteom.* **2009**, 8, pp 2170-2185.
11. Leibiger, H; *et al.* Variable domain-linked oligosaccharides of a human monoclonal IgG: structure and influence on antigen binding. *Biochem. J.* **1999**, 338, pp 529-538.
12. Anumula, K. R. Advances in fluorescence derivatization methods for high-performance liquid chromatographic analysis of glycoprotein carbohydrates. *Anal. Biochem.* **2006**, 350, pp 1-23.
13. Montreuil, J; Vliegthart, J. F. G; Schachter, H., Eds.; Glycoproteins II. Elsevier B. V., Amsterdam, **1997**.
14. Ceroni, A; *et al.* GlycoWorkbench: a tool for the computer-assisted annotation of mass spectra of glycans. *J. Proteome Res.* **2008**, 7, pp 1650-1659.

Agilent Solutions for High-throughput N-linked Glycan Profiling from Biotherapeutics

Authors

Arunkumar Padmanaban,
Sreelakshmy Menon, and
Suresh Babu CV

Agilent Technologies India Pvt
Ltd Bangalore India

Abstract

This Application Note presents an automated high-throughput sample preparation for 2-AB labeled N-linked glycans from innovator and biosimilar monoclonal antibody (mAb) followed by liquid chromatography (LC) analysis.

Agilent AssayMAP Bravo liquid handling platform was used for the automated sample preparation involving glycan cleavage and release from mAb, 2-AB labeling, and purification. The downstream ready samples were then analyzed using an Agilent 1290 Infinity LC system with Agilent AdvanceBio Glycan Mapping columns. The glycan distribution and heterogeneity between the samples were deduced by comparing the chromatogram from both innovator and biosimilar mAb.

The study highlights the high-throughput application of the AssayMAP Bravo platform for automated and reproducible sample preparation for glycan profiling, followed by fast chromatographic separation using a 1290 Infinity LC system.

Introduction

Monoclonal antibodies (mAbs) are leading the race in biotherapeutics, and have revolutionized the way diseases are treated and intervened. Patents for most of the 20 clinically approved¹ first-generation mAbs have either expired, or are about to expire. This has increased the opportunity for generating generic versions, referred to as biosimilars. Regulatory bodies such as the United States Food and Drug Administration (FDA) and European Medicines Agency (EMA) have published guidelines for the characterization of biosimilar protein therapeutics including post-translation modifications². Glycosylation modification results in structural heterogeneity that changes an mAb's target binding capacity, stability, charge, and mass³. During the course of developing the antibodies from the cell line, the glycans are monitored for cell line alteration and culture optimization. This requires a high-throughput sample preparation to screen several cell culture conditions in parallel. The Agilent Bravo Automated Liquid Handling Platform is a precise liquid-handling system designed for high-throughput applications like glycan profiling. Agilent Bravo with Agilent AssayMAP technology combines the automation with miniature 5 μ L pack bed cartridges for sample preparations, enabling high-throughput chromatography and sample preparations in a shorter time⁴.

This Application Note presents a workflow solution for profiling N-glycans from mAb using an Agilent AssayMAP Bravo platform. The system automates the N-glycan purification and derivatization using 2-aminobenzamide (2-AB) dye, which is then followed by Agilent HPLC analysis using Agilent Glycan Mapping columns. We showcase the ease-of-use of an AssayMAP Bravo for the automated sample preparation required for high-throughput profiling studies.

Experimental

Materials

ProZyme's GlykoPrep-plus Rapid N-Glycan Sample Preparation with 2-AB kit, and 2-AB labeled Oligomannose 5 and 6 were purchased from ProZyme. An Agilent AdvanceBio Glycan Mapping column, 2.1 \times 150 mm, 1.8 μ m (p/n 859700-913), and Agilent 2-AB labeled human IgG N-linked glycans standards (p/n 5190-6996) were purchased from Agilent Technologies. Innovator rituximab and biosimilar product were purchased commercially from a local pharmacy. All other chemicals were procured as HPLC grade from Sigma.

Sample preparation

Reagents from a ProZyme GlykoPrep plus kit with an AssayMAP protocol were used for the sample preparation. An innovator and biosimilar sample was diluted to 1 mg/mL, then loaded onto three columns of the 96-well plate (24 replicates each). After placing the sample plates and reagents plates as specified in the guideline⁵, the samples were processed by launching the N-Glycan Sample Prep: RX digestion & 2-AB labeling module from Agilent VWorks software. The protocol consists of five modules that are performed in sequential order to immobilize the samples, digest the glycans, elute, label with 2-AB, and complete a final cleanup to remove the excess dyes. The final purified labeled glycan from each well were then transferred to HPLC vials and analyzed immediately, or stored at -80°C .

Figure 1 presents a schematic diagram of the complete workflow.


Figure 1. Schematic workflow for the glycan isolation, labelling, purification followed by LC analysis using Agilent solutions. Components of the Agilent AssayMAP steps are shown in blue; the ProZyme kits are shown in orange.

Instrumentation

Agilent 1290 Infinity LC System including:

- Agilent 1290 Infinity Binary Pump G4220A
- Agilent 1290 Infinity Autosampler G4226A
- Agilent 1290 Infinity TCC G1316 C
- Agilent 1260 Fluorescence Detector G1321 B

The LC method described earlier was adopted for this study⁶. Every sample was analyzed in quadruplet injections, followed by a blank injection.

Software

- Agilent VWorks Automation Control 11.4.0.1233
- Agilent AssayMAP Launch Pad 3.0
- Agilent N-Glycan Sample Prep: RX digestion & 2-AB labeling protocols 1.0
- Agilent ChemStation C.01.06

Results and Discussion

Agilent AssayMAP sample preparation

The samples were processed following a protocol consisting of five modules, as presented in Figure 2. The software suite consists of a deck layout, labware table, and application settings tab for each module to be performed (Figure 3). The user was prompted to place the appropriate consumables and reagents listed in the labware table in the specified deck positions. After setting up all labware, the protocol was executed, and AssayMAP Bravo completed the protocol, and a confirmation message was displayed to proceed to the next module.

The final Cleanup Protocol module eluted the labeled and purified glycans in an aqueous buffer into a clean 96-well plate. The samples were then analyzed in quadruplets, along with blanks, using the Agilent AdvanceBio Glycan Mapping column for the downstream LC analysis.


Figure 2. Agilent VWorks modules of an Agilent AssayMAP N-Glycan Sample Preparation.


Figure 3. The deck layout, labware table, and application settings of the Immobilization and Digestion Protocol, one of the five sample preparation modules in Agilent VWorks software.

HPLC Analysis

The N-glycan profiles were compared between the innovator and biosimilar mAbs using a fast analysis method with the Agilent 1290 Infinity system. The Agilent 2-AB labeled IgG N-linked glycan library was used as the standard to assess the column performance, and to annotate the peaks in mAb samples. The 1290 Infinity system with high backpressure enabled the analysis of the samples in less than 6 minutes, saving analysis time per sample, and increasing sample throughput.

Each sample was HPLC analyzed in replicates of four. Figure 4 presents the peak area of four major glycans species from 96 replicates, demonstrating the excellent column-to-column reproducibility of the AssayMAP micro chromatography pipette tips. The Reproducibility Standard Deviation (RSD) calculated for peak area and peak height showed a coefficient of variation (CV) of less than 6 % for all glycan species. This demonstrates the very robust and reproducible sample preparation capability of the AssayMAP Bravo system.

Glycan profiling and comparison of innovator and biosimilar rituximab

The chromatograms of the innovator and biosimilar rituximab were compared with the standard N-linked IgG glycan library, and the peaks corresponding to glycans were annotated. Separate standards comprising 2-AB labeled oligomannose 5 and 6 were also used to annotate additional peaks. The glycosylation pattern of the major abundant glycans, such as the G0F, G1F isoforms, and G2F was comparable between the innovator and biosimilar product (Figure 5). Small differences in the low abundant glycans were observed; the biosimilar sample contained lower amounts of mannose (Man5), as shown in the zoomed view of Figure 6. Despite minor differences in some low abundant glycans, the glycan profile of the innovator and biosimilar rituximab can be concluded to be comparable.


Figure 4. Peak area of major glycan species from 96 replicates for each innovator and biosimilar sample.


Figure 5. Glycosylation pattern of innovator and biosimilar rituximab.

The area percentage for each N-glycan core was calculated and compared between both samples. Table 1 shows the ratio of N-glycans between the innovator and biosimilar rituximab. The table shows a similar N-glycan fingerprint profile for both innovator and biosimilar.

Both samples showed a similar trend in the distribution of N-glycans, with G0F as the major glycoform followed by G1F, G1F', and G2F. Apart from the annotated glycans, there were a few unknown glycan forms that may be assigned with an orthogonal detection system. A correlation graph (Figure 7) plotted for the area percentage of both samples shows high similarity, with an R2 of 0.973.

Table 1. Glycan area percentage between innovator and biosimilar mAbs.

Percentage of N-glycans		
N-glycans	Innovator	Biosimilar
G0	1.5	0.8
G0F	42.2	50.5
Man5	1.7	1.2
G1F	28.6	26.7
G1F'	13.0	11.6
G2	0.3	0.2
G2F	8.7	5.5
G1FS1	0.9	0.6
A1F	0.2	0.5

Conclusions

This study highlights the versatility of the Agilent AssayMAP Bravo system for a high-throughput sample preparation.

- The Agilent VWorks Automation Control software suite simplifies the sample preparation with ready-to-go protocols, resulting in downstream-compatible samples with minimal hands-on operation.
- Monoclonal antibody samples from a 96-well plate were processed, in parallel, for enzymatic glycan cleavage, separation, derivatization with 2-AB, and purification.
- The AssayMAP demonstrated excellent reproducibility in the glycan purification, and performed robustly.
- The purified samples were then analyzed using an Agilent 1290 Infinity LC system with the Glycan Mapping column.
- The glycan species were well-resolved in a shorter time, and were annotated using standards.
- The distribution of the glycan species between the innovator and biosimilar were assessed, and the data suggest comparable glycan profiles for the innovator and biosimilar rituximab used in this study.


Figure 6. Zoomed view of Figure 5: the glycosylation pattern of innovator and biosimilar rituximab, showing minor differences in the low abundant glycans.


Figure 7. Linear graph showing the comparison of the area % of glycans released from innovator and biosimilar rituximab.

References

1. Chames, P.; Baty, D. Bispecific antibodies for cancer therapy. *mAbs* **2009**, 1:6, 539–547.
2. Scientific Considerations in Demonstrating Biosimilarity to a Reference Product: Guidance for Industry. Food and Drug Administration. Center for Biologics Evaluation and Research (CBER). April 2015. <http://www.fda.gov/downloads/Drugs/GuidanceComplianceRegulatoryInformation/Guidances/UCM291128.pdf>
3. Zhang, Q.; et al. A Perspective on the Maillard Reaction and the Analysis of Protein Glycation by Mass Spectrometry: Probing the Pathogenesis of Chronic Disease. *Journal of Proteome Research* **2009**, 8(2), 754–769.
4. Agilent AssayMAP Cartridges Automated Sample Preparation, Agilent Technologies Data Sheet, publication number 5991-1774EN, **2013**.
5. AssayMAP with ProZyme GlykoPrep-plus Rapid N-Glycan Sample Preparation with 2-AB Protocol Guide, Agilent Technologies Manual, publication number G5409-90012EN, Rev May **2013**.
6. Martosella, J.; et al. A Novel HILIC Column for High Speed N-linked Glycan Analysis, Agilent Technologies Application Note, publication number 5991-4886EN, **2014**.

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1. Guidance for Industry Immunogenicity Assessment for Therapeutic Protein Products, U.S. Department of Health and Human Services Food and Drug Administration Center for Drug Evaluation and Research (CDER) Center for Biologics Evaluation and Research (CBER), **2014**.
2. Mahler, H-C.; et al. Protein Aggregation: Pathways, Induction Factors and Analysis, *J. Pharm. Sci.* 2008, *98*(9).
3. Schneider, S. 2D-LC/MS Characterization of Charge Variants Using Ion Exchange and Reversed-Phase Chromatography, *Agilent Technologies Application Note* publication number 5991-6673EN **2016**.

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5991-7476JAJP

Agilent Technologies

Size Exclusion Chromatography in the Presence of an Anionic Surfactant

Intact Protein Profiling

Authors

Andy Coffey
Agilent Technologies, Inc

Abstract

Sodium dodecyl sulfate (SDS, or SLS) is a well known anionic detergent, frequently used to denature proteins. It is commonly used in polyacrylamide gel electrophoresis (SDS-PAGE), where a remarkably consistent level of binding across a wide range of proteins imparts a reliable charge-to-mass ratio. This allows separation of denatured proteins based on relative size due to their relative ion mobility. Conversely, size exclusion chromatography (SEC) for size-based separation of proteins is normally performed under nondenaturing conditions using predominantly aqueous buffers as mobile phase.

This Application Note used SEC with light scattering detection to investigate the impact of varying SDS concentration on the protein molecule. This was achieved by studying the changes in RT, apparent molecular weight, and hydrodynamic radius.

Introduction

To denature proteins, sodium dodecyl sulfate (SDS) is used above its critical micelle concentration (CMC). The CMC of SDS in pure water is 8.2 mM (approximately 0.2 % w/v). However, the presence of buffer salts or changes in pH and temperature may reduce the CMC significantly. CMC values of around 1.0 mM (approximately 0.04 % w/v) are observed in phosphate buffered saline (PBS), for example. A micelle forms when approximately 62 SDS molecules coalesce into a spherical shape with a hydrophobic core surrounded by an anionic surface, as depicted in Figure 1. Typically, SDS is used at significantly higher concentrations of 2 % w/v for sample preparation for SDS-PAGE, but also with a reducing agent to cleave disulfide bonds within the protein molecule. The resulting saturated SDS-protein complex has, on average, 1.4 g SDS per gram of protein [1]. Such denatured proteins are considered to adopt a rigid cylindrical shape resulting in the observation that ion mobility in gel electrophoresis is proportional to molecular weight. By omitting the reduction of disulfide bonds, it has been observed that the ratio of SDS to protein decreases significantly and the time to reach equilibrium may be doubled [2].

In contrast, adding SDS (approximately 10:1 mol ratio) has been found to arrest heat denaturation of BSA solution, which otherwise leads to creation of high molecular weight oligomers [3], or reduces the level of noncovalent aggregation. Historically, it was recommended to use 0.1 % w/v SDS in the mobile phase to reduce retention time (RT), and improve peak shape in protein size exclusion chromatography (SEC) [4].

The mechanism by which SDS denatures a protein is not fully understood. Since both the detergent molecule and the protein possess complementary ionic and hydrophobic regions, it is clear that a combination of mechanisms may be involved. Thermodynamic studies have suggested that interactions are not ionic alone. It is conceivable that SDS first begins to bind through an ionic interaction with positively charged amino acids side chains (Lys and Arg) on the surface of the protein. Increasing SDS concentration leads to higher levels of incorporation, and it is postulated that ultimately, the flexible protein chain is decorated with SDS micelles rather than the rigid cylinder model originally proposed [5].


Figure 1. Cartoon depicting SDS micelle formed at concentrations above CMC (aggregation number ~ 62).

Bovine serum albumin (BSA) is a well studied protein due to its capability of acting as a carrier of smaller molecules between tissues and cells. The structure and size (hydrodynamic radius) of the BSA monomer has been determined using multiple analytical techniques. Its primary structure is 583 amino acids with a molecular weight of 66,463 Da, Figure 2. Commercial samples frequently contain varying amounts of dimer and higher-order aggregates, and different isolation techniques may result in varying amounts of these oligomers. Since BSA contains 35 cysteine residues, there are 17 disulfide bridges and one free thiol group. It is believed that BSA oligomers are predominantly covalently linked through disulfide bonds involving the extra thiol group. Therefore, in the absence of reducing agents, denaturation with SDS will not result in such covalently linked dimers and higher aggregates reverting to a monomeric form.


Figure 2. Native heart-shaped conformation of bovine serum albumin.

To explore the effect of SDS denaturation on the characteristics of BSA, a series of experiments were performed. Using SEC of the nonreduced BSA protein, the amount of dimer and higher-order aggregates could be determined. The addition of a light scattering detector with dynamic light scattering (DLS) capability allowed the measurement of both molecular weight and hydrodynamic radius (Rh). The eluent chosen for this series of experiments was phosphate buffered saline (PBS), containing 10 mM phosphate buffer, and 140 mM NaCl, pH 7.4 (Eluent A). A second mobile phase was prepared in an identical manner, but with the addition of 2 % w/v SDS (Eluent B). The HPLC instrument was then run with increasing levels of SDS incorporated in the mobile phase by varying the proportions of Eluent A and Eluent B. Throughout the experiment, the same BSA sample was used, prepared at a concentration of 10 mg/mL in PBS alone.

The experiment was designed to determine:

- If SDS can denature BSA oligomers
- What effect is observed on the molecular weight measured using inline light scattering detection
- What change in hydrodynamic radius is observed

Conditions

Parameter	Value
Column:	Agilent AdvanceBio SEC 300 Å, 2.7 µm, 7.8 × 300 mm (p/n PL1180-5301)
Samples:	BSA (Sigma-Aldrich) Sample prepared 10 mg/mL in Eluent A
Eluent A:	PBS, pH 7.4 (10 mM phosphate, 140 mM NaCl)
Eluent B:	PBS, pH 7.4 (10 mM phosphate, 140 mM NaCl) + 2.0 % w/v SDS
Gradient:	Isocratic elution at 0 % B; 10 % B; 20 % B; 30 % B; 40 % B; 50 % B; 60 % B; 70 % B
Flow rate:	0.8 mL/min
Detector:	UV, 280 nm; LS 15°, LS 90°
Instrument:	Agilent 1260 Infinity Bio-inert Quaternary LC with Agilent 1260 Infinity Bio-SEC multidetector system.

Results and Discussion

Following two blank injections to allow for the column to become conditioned in the mobile phase, duplicate injections of BSA 10 mg/mL were made.

Figure 3 shows overlaid chromatograms (UV 280 nm signal) of four individual experiments at different mobile phase compositions. Table 1 contains RT information relating to the entire experimental series.

From Table 1, it is apparent that there is little change in total peak area for different runs. This indicates that increasing SDS concentration does not result in the absorption of the protein, or loss of protein through induced aggregation. RTs become shorter but stabilize from 0.6 to 1.4 % SDS (w/v), as seen in Figure 4. This shortening of RT also results in lower resolution between the peaks, so it proved difficult to accurately determine how much dimer and higher-order aggregates became denatured. It is clear that these multimer peaks are still abundant and must, therefore, be predominantly covalently linked.

Table 1. RT and total peak area.

% SDS in mobile phase	RT (min)			Total peak area
	Higher-order aggregates	Dimer	Monomer	
0.0 %	6.00	6.46	7.29	3,632
0.2 %	5.22	5.59	6.45	3,633
0.4 %	4.96	5.23	5.96	3,668
0.6 %	4.89	5.15	5.83	3,588
0.8 %	4.89	5.13	5.80	3,571
1.0 %	4.89	5.11	5.77	3,572
1.2 %	4.88	5.11	5.76	3,536
1.4 %	4.90	5.10	5.75	3,566


Figure 3. SEC chromatograms of BSA with increasing SDS concentration in mobile phase.


Figure 4. Effect of % SDS in mobile phase on RT.

In SEC, shorter RT is indicative of an increase in the size of a molecule in solution, but does not necessarily infer an increase in molecular weight. To ascertain the molecular weight of the three major peaks of the chromatogram, the BioSEC software was used to calculate the molecular weight using the relationship described in the following formulae [6]:

Where:

(LS) = Light scattering detector signal

K_{LS} = LS detector constant

dn/dc = Specific refractive index increment*

(UV) = UV detector signal (280 nm)

ϵ = Extinction coefficient

* Using the approximation of $dn/dc \approx 0.186$ mL/g for nonglycosylated proteins

Figure 5 shows the regions of each peak chosen to perform the light scattering molecular weight determination to try to minimize interference due to anticipated peak overlap. Table 2 shows the results for the molecular weight measurements. This table includes columns where the relative number of monomer units has been calculated. It is important to recognize that the molecular weight of the higher-order aggregates peak is nearly always three times larger than the molecular weight of the monomer peak, indicating it is a trimer. The molecular weight of the dimer peak is two times larger than the molecular weight of the monomer peak, as expected. In addition, the molecular weight of each of the peaks increases with increasing SDS concentration until a plateau is reached at 0.6% w/v SDS, in agreement with the stabilization of RTs at this point.


Figure 5. Molecular weight calculation regions for (1) higher-order aggregates, (2) dimer, and (3) monomer peaks of BSA.

Table 2. Molecular weight analysis results from light scattering detection.

SDS w/v%	Mass gain			SDS:BSA (g/g)		
	Higher-order aggregates	Dimer	Monomer	Higher-order aggregates	Dimer	Monomer
0.0	202,197	136,383	67,838	3.0	2.0	1.0
0.2	242,541	194,876	107,905	2.2	1.8	1.0
0.4	560,154	359,239	203,737	2.7	1.8	1.0
0.6	667,318	423,418	212,596	3.1	2.0	1.0
0.8	615,726	404,898	207,028	3.0	2.0	1.0
1.0	624,648	404,634	205,868	3.0	2.0	1.0
1.2	617,385	405,618	206,376	3.0	2.0	1.0
1.4	612,031	406,699	206,614	3.0	2.0	1.0

Table 3. Mass gain for BSA with increasing surfactant concentration, showing steady state at concentrations >0.6%.

SDS w/v%	Mass gain			SDS:BSA (g/g)		
	Higher-order aggregates	Dimer	Monomer	Higher-order aggregates	Dimer	Monomer
0.0	0	0	0	0.0	0.0	0.0
0.2	40,344	58,493	40,067	0.2	0.4	0.6
0.4	357,957	222,856	135,899	1.8	1.6	2.0
0.6	465,121	287,035	144,758	2.3	2.1	2.1
0.8	413,529	268,515	139,190	2.0	2.0	2.1
1.0	422,451	268,251	138,030	2.1	2.0	2.0
1.2	415,188	269,235	138,538	2.1	2.0	2.0
1.4	409,834	270,316	138,776	2.0	2.0	2.0

Since the observed molecular weight increases in line with the concentration of SDS present in the mobile phase, it may be inferred that each molecular species is actually gaining mass through accumulation of associated SDS. However, the observed mass gain is considerable and higher than expected, particularly since this BSA sample has not been reduced (Table 3). The steady state mass gain under the analysis conditions is 2 g SDS per 1 g protein, suggesting that SDS micelles may be accumulating along the protein molecule regardless of the number of oligomers it contains.

Closer inspection of the concentration-dependant UV trace of the chromatograms obtained at 0.0 % SDS and 1.4 % SDS shows further differences. There is clearly a reduction in resolution for the higher molecular weight species, however some smaller peaks eluting after the monomer peak have become evident in the 1.4 % SDS mobile phase composition (Figure 6). It was not possible to identify these lower molecular weight species.

Figure 7A shows the signals from the LS 90° detector for BSA at 0.0 % SDS and 1.4 % SDS concentrations. Since the light scattering detector is mass sensitive rather than concentration-dependant, it is not surprising that the signal obtained under the higher surfactant concentration conditions is much larger than the signal under native conditions.

In addition to the increased response seen in Figure 7A at the highest SDS concentration, there is also a vertical offset. The reason for this becomes apparent when the DLS data are analyzed to determine the Rh. Not only do we see an increase in the Rh of the monomer peak (from around 3.8 to 5.7 nm, close to literature values) [7,8,9], there is a continual background of particles detected of approximately 2.4–2.8 nm in size, almost certainly due to SDS micelles in solution.


Figure 6. Close up of the UV 280 nm signal of BSA in PBS mobile phase containing 0.0 % SDS and 1.4 % SDS.


Figure 7A. Overview of the LS 90° signal of BSA in PBS mobile phase containing 0.0 % SDS and 1.4 % SDS.


Figure 7B. Overview of the LS 90° signal of BSA in PBS mobile phase containing 0.0 % SDS and 1.4 % SDS with hydrodynamic radius (Rh) data added.

Conclusions

At the outset, this experiment was designed to determine if SDS can denature BSA oligomers using SEC with light scattering detection. It was found that accurate quantification of oligomer content was not possible through loss in resolution as the individual peaks corresponding to monomer, dimer, and higher-order aggregates eluted closer together. However, it is apparent the three individual species remain intact and become saturated with SDS to the same extent.

Simply denaturing the protein, making the molecule unfold, would be expected to give a larger structure in solution, leading to a shorter RT in SEC. However, as the SDS concentration was increased from 0.0 to 0.4 % w/v, the molecular weight of the monomer, dimer, and higher-order aggregates increased proportionally. At higher SDS concentrations (0.6–1.4 % w/v) the RTs (and, therefore, size in solution) and the molecular weight determined by light scattering stabilized as the molecular species were saturated by SDS.

Using DLS capability to determine hydrodynamic radius shows an increase in size comparable to literature values for a reduced BSA sample. The size increase does not reflect a change in conformation of the BSA sample run under nonreducing conditions, but replicates the observed increase in molecular weight as the molecule complexes with and becomes saturated by SDS micelles.

It is clear that, although the Agilent AdvanceBio SEC 300 Å column is able to tolerate SDS in the mobile phase, as seen by the good peak shape and reproducible recovery, the SEC separation is badly affected by the dramatic increases in size and molecular weight of the protein species being analyzed. The use of surfactants in the mobile phase for SEC should be avoided wherever possible.

References

1. Reynolds, J. A.; Tanford, C. The gross conformation of protein-sodium dodecyl sulfate complexes. *J. Biol. Chem.* **1970**, *245*(19), 5161-5165.
2. Pitt-Rivers, R.; Impiombato, F. S. A. The binding of sodium dodecyl sulphate to various proteins. *The Biochemical Journal* **1968**, *109*, 825-830.
3. Aoki, K.; Hiramatsu, K.; Kimura, K.; Kaneshina, S.; Nakamura, Y.; Sato, K. Heat Denaturation of Bovine Serum Albumin. I : Analysis by Acrylamide-gel Electrophoresis. *Bulletin of the Institute for Chemical Research, Kyoto University* **1969**, *47*(4), 274-282
4. Hagarová, D.; Horváthová, M.; Žúbor, V.; Breier, A. Optimization of conditions for size-exclusion chromatography of proteins. *Chemical Papers* **1991**, *45*(3), 341-348.
5. Shirahama, K.; Tsujii, K.; Takagi, T. Free-boundary electrophoresis of sodium dodecyl sulfate-protein polypeptide complexes with special reference to SDS-polyacrylamide gel electrophoresis. *J. Biochem.* **1974**, *75*, 309-319.
6. Wen, J.; Arakawa, T.; Philo, J. S. Size-Exclusion Chromatography with On-Line Light-Scattering, Absorbance, and Refractive Index Detectors for Studying Proteins and Their Interactions. *Anal. Biochem.* **1996**, *240*, 155-166.
7. Lorber, B.; Fischer, F.; Bailly, M.; Roy, H.; Kern, D. Protein Analysis by Dynamic Light Scattering: Methods and Techniques for Students. *Biochem. and Mol. Biol. Ed.* **2012**, *40*(6), 372-382.
8. Valstar, A.; Almgren, M.; Brown, W. The Interaction of Bovine Serum Albumin with Surfactants Studied by Light Scattering. *Langmuir* **2000**, *16*(3), 922-927.
9. Tanner, R. E.; Herpigny, B.; Chen, S.-H.; Rha, C. K. Conformational change of protein sodium dodecylsulfate complexes in solution: A study of dynamic light scattering. *J. Chem. Phys.* **1982**, *76*(8), 3866-3872.

Calibrating your Agilent AdvanceBio SEC Columns

Technical Overview

Authors

Andy Coffey
Agilent Technologies, Inc

Introduction

Agilent AdvanceBio SEC columns are specifically designed for aqueous size exclusion chromatography (SEC) of biomolecules. The particles have been manufactured using proprietary technology to combine optimum pore size and pore volume for separating molecules such as proteins, polysaccharides, or aqueous soluble polymers.

Even if the intention is to use the AdvanceBio SEC column for quantification of monomer and dimer content, it is still good practice to regularly perform a calibration using appropriate molecular weight standards. Regular calibration ensures better reproducibility, leading to improved accuracy, and enables earlier detection of potential problems, reducing system downtime and troubleshooting.

For protein separations, the standards should be a range of well characterized proteins covering the entire operating range of the column. For polysaccharide analysis, pullulan polysaccharide molecular weight standards should be used, and for synthetic polyethylene glycol or polyethylene oxide, PEG/PEO standards are suitable.

SEC is a relatively straightforward technique. It relies on simple diffusion into the pore structure of the stationary phase; larger molecules cannot permeate the particles, and elute first, while smaller molecules diffuse readily into the pores, and elute later. This provides two key aspects for the successful use of SEC:

- There should be no secondary interactions between the analyte and the stationary phase
- The pore size should be chosen to match the size of molecules being analyzed

Proteins are complex molecules containing numerous side chain functionalities: acid, basic, neutral, and hydrophobic. Finding the optimum conditions to avoid secondary interactions can be challenging, however the AdvanceBio SEC product range has a polymeric surface coating applied to the silica particle that overcomes many of these issues.

The mechanism of separation relies on differences in size of molecules in solution. Protein structures are often compact and globular in nature, and proteins often aggregate under stress conditions such as extremes of temperature, pH, or salt composition and for dimers and larger units [1]. This is a particular issue for molecules such as monoclonal antibodies, where the presence of aggregated proteins can lead to adverse effects if administered. SEC provides the ideal tool for quantifying and monitoring protein aggregation. In comparison, polysaccharides and synthetic polymers such as polyethylene oxide do not possess a wide range of functional groups, but are often made up of a series of closely related oligomers. In this case, SEC can provide the analysis of molecular weight, molecular weight distribution, and branching information [2].


Figure 1. Resolving ranges for Agilent AdvanceBio SEC columns.

Calibrating your Agilent AdvanceBio SEC column

Choose the eluent for analysis. The samples need to be soluble in the mobile phase, and there should be no secondary interactions that lead to molecules eluting earlier or later than expected. For proteins, we recommend 150 mM sodium phosphate, pH 7.0 as this eliminates sodium chloride, which may lead to corrosion, from your HPLC system. However, other mobile phases including phosphate buffered saline (PBS), or high ionic strength mobile phases such as 0.1 M sodium phosphate/0.1 M sodium sulphate may still be used (see Pages 10–11).

Dissolve the standards in the mobile phase. Ensure that the solution is clear and there is no insoluble material present. If necessary, filter the sample. Consider changing the mobile phase to ensure that the samples fully dissolve.

Record the chromatogram for each of the standards, and plot the retention time (x-axis) versus the molecular weight (y-axis) to create the calibration curve. Notice that the y-axis is plotted on a logarithmic scale.

Conditions

Parameter	Value
Column:	Agilent AdvanceBio SEC 130 Å, 2.7 µm, 7.8 × 300 mm (p/nPL1180-5350) Agilent AdvanceBio SEC 300 Å, 2.7 µm, 7.8 × 300 mm (p/nPL1180-5301)
Samples:	Agilent polyethylene glycol calibration kit, PEG-10, 10 × 0.2 g (p/nPL2070-0100) Agilent polyethylene oxide calibration kit, PEO-10, 10 × 0.2 g (p/nPL2080-0101) AdvanceBio SEC 130 Å Protein Standard, lyophilized, 1.5 mL (p/n5190-9416) AdvanceBio SEC 300 Å Protein Standard, lyophilized, 1.5 mL (p/n5190-9417) Protein standards (Sigma-Aldrich) Samples prepare at 0.5 – 1.0 mg/mL in mobile phase
Eluent A:	150 mM Sodium phosphate buffer, pH 7.0, PBS, pH 7.4 (10 mM phosphate, 140 mM NaCl) or 100 mM Sodium phosphate + 100 mM sodium sulfate, pH 7.0
Flow rate:	1.0 mL/min
Detector:	RI for polysaccharides and PEG / PEO standards UV, 220 nm for peptides and proteins
System:	Agilent 1260 Infinity Bio-inert LC (with additional refractive index detector)

Peptide and Protein Calibration

Proteins and peptides are generally very compact molecules containing multiple hydrogen bonds, electrostatic interactions, or covalent bonds such as disulphide bridges. Choosing a range of standards covering the operating range of the column for calibration is preferred. It is also possible to use prepared standard mixtures, as peaks are likely to be clearly defined and well resolved (Figures 3A and 3B).

The presence of a diverse range of side chain functionality from different amino acids may mean undesirable secondary interactions can occur, and so different mobile phase compositions may need to be tested.

The peaks shown are generally very sharp and well defined since they represent a single molecular species. It is quite common for proteins to contain aggregates, and these too are typically well resolved – evident as an earlier eluting peak or series of peaks.

Protein/Peptide	MW	Retention time (min)	
		130 Å	300 Å
Thyroglobulin	670,000	4.60	5.14
γ-Globulin	150,000	4.90	6.53
BSA	66,000	5.53	7.57
Ovalbumin	44,300	6.04	8.13
Myoglobin	17,600	6.77	8.79
Cytochrome C	12,327	6.95	8.92
Aprotinin	6,511	7.56	9.38
Neurotensin	1,672	9.42	10.54
Angiotensin-II	1,040	9.94	10.82


Figure 2. Peptide and protein calibration on Agilent AdvanceBio SEC columns.

Representative chromatograms for protein and peptide standards


Figure 3A. Individual peptide and proteins contained in Agilent AdvanceBio SEC 300 Å Protein Standard (p/n5190–9417) on an Agilent AdvanceBio SEC 300 Å column.


Figure 3B. Individual peptides and proteins contained in Agilent AdvanceBio SEC 130 Å Protein Standard (p/n5190–9416) on an Agilent AdvanceBio SEC 130 Å column.

Polysaccharide Calibration

Calibration with polysaccharide standards should be used when your analyte is also a polysaccharide, for example starch or cellulose. There are many types of polysaccharide, including linear, branched, or sulfated molecules. Unlike peptides or proteins, they comprise of a range of chain lengths and can extend to very large sizes. Therefore, it is important to ensure that the column you are using is capable of resolving the molecular weight range of polysaccharide you are investigating.

The mass difference of 162 for each additional hexose means that separating individual oligomers is very difficult unless the sample is sufficiently small (Figure 5C), achieved using two columns in series for additional resolution.

The combination of two columns run in series can increase resolution to the extent that individual oligomers are becoming visible in the 5,800 Da polysaccharide sample.

This sample contains oligomers of approximately 30–40 units, with a mass difference of 162 Da each.

Polysaccharide	MW	Retention time (min)	
		130 Å	300 Å
788K	788,000	4.48	4.74
380K	380,000	4.58	4.84
100K	100,000	4.73	5.33
48K	48,000	4.89	6.19
23.7K	23,700	5.39	7.42
5.8K	5,800	7.16	9.22
Maltotriose	504	9.46	10.72
Maltose	342	9.73	10.82
Glucose	180	10.01	11.09


Figure 4. Polysaccharide calibration on Agilent AdvanceBio SEC columns.

Representative chromatograms for polysaccharide standards


Figure 5A. Polysaccharide standards on an Agilent AdvanceBio SEC 300 Å column.


Figure 5B. Polysaccharide standards on an Agilent AdvanceBio SEC 130 Å column.


Figure 5C. Polysaccharide 5.8K run on Agilent AdvanceBio SEC columns 300 Å + 130 Å in series.

PEG/PEO Calibration

Polyethylene glycol (PEG) and polyethylene oxide (PEO) are synthetic polyethers, available in a range of sizes covering a diverse molecular weight range. Similar to polysaccharides, they comprise of a range of chain lengths, however the mass difference is just 44 Da.

PEG/PEO	MW	Retention time (min)	
		130 Å	300 Å
905K	905,000	4.41	4.73
692K	692,000	4.41	4.73
498.6K	498,600	4.58	4.83
305.5K	305,500	4.55	4.85
135.8K	135,800	4.65	4.93
77.35K	77,350	4.71	5.20
46.47K	46,470	4.77	5.67
21.3K	21,300	5.11	6.81
12.14K	12,140	5.65	7.74
8.73K	8,730	6.10	8.28
3.87K	3,870	7.26	9.32
1.48K	1,480	8.59	10.26
420	420	9.79	11.03
106	106	10.64	11.44


Figure 6. PEG/PEO Calibration on Agilent AdvanceBio SEC columns.

Representative chromatograms for PEG/PEO standards


Figure 7A. PEG/PEO standards on an Agilent AdvanceBio SEC 300 Å column.


Figure 7B. PEG/PEO standards on an Agilent AdvanceBio SEC 130 Å column.

Effect of Mobile Phase Composition

The impact of mobile phase on the size of the molecule in solution, or in overcoming secondary interactions needs to be taken into consideration, as this will influence the retention time and alter the shape of the calibration curve. This is clearly visible in the comparison of three different mobile phases: 150 mM sodium phosphate, pH 7.0; PBS, pH 7.4 (approximately 10 mM sodium phosphate, 140 mM NaCl); 0.1 M sodium phosphate + 0.1 M sodium sulfate, pH 7.0 (Figure 8).

For molecules that do not present secondary interactions, such as polysaccharides, the effect of mobile phase composition is minimal (Figure 9).


Figure 8. Effect of mobile phase composition on column calibration with protein standards.


Figure 9. Effect of mobile phase composition on column calibration with polysaccharide standards.

When choosing a mobile phase for peptide or protein analysis by SEC, it is important to determine the effect that differences in pH or ionic strength may have on the sample of interest. Figures 10A, 10B, and 10C show a noticeable shift in the retention time of BSA under different mobile phase conditions. Resolution factors between dimer and monomer peaks may also be affected, so method optimization and method robustness should be fully explored.


Figure 10A. γ -Globulin and BSA on an Agilent AdvanceBio SEC 300 Å column using 150 mM sodium phosphate, pH 7.0.


Figure 10B. γ -Globulin and BSA on an Agilent AdvanceBio SEC 300 Å column using PBS, pH 7.4.


Figure 10C. γ -Globulin and BSA on an Agilent AdvanceBio SEC 300 Å column using 100 mM sodium phosphate + 100 mM sodium sulfate, pH 7.0.

Conclusion

Calibrating your Agilent AdvanceBio SEC size exclusion column with appropriate standards ensures you understand the correct working range, and can allow you to use calibration curves to estimate the molecular size of unknown molecules. However, regular calibration with a selection of standards is beneficial, and can be used to monitor column performance over a period of time, allowing early detection of potential problems. In turn, corrective action can be taken, ultimately reducing system downtime and improving productivity.

Agilent offers a range of column dimensions and different pore sizes suitable for differently sized molecules. Together with the extensive range of Agilent LC instrumentation, consumables, and technical support, we can provide the ultimate workflow solution for your needs.

Ordering Information

Part Number	Description
PL1180-5301	Agilent AdvanceBio SEC 300 Å, 2.7 μ m, 7.8 \times 300 mm
PL1180-3301	Agilent AdvanceBio SEC 300 Å, 2.7 μ m, 7.8 \times 150 mm
PL1180-1301	Agilent AdvanceBio SEC 300 Å, 2.7 μ m, 7.8 \times 50 mm guard
PL1580-5301	Agilent AdvanceBio SEC 300 Å, 2.7 μ m, 4.6 \times 300 mm
PL1580-3301	Agilent AdvanceBio SEC 300 Å, 2.7 μ m, 4.6 \times 150 mm
PL1580-1301	Agilent AdvanceBio SEC 300 Å, 2.7 μ m, 4.6 \times 50 mm guard
5190-9417	Agilent AdvanceBio SEC 300 Å Protein Standard, lyophilized, 1.5 mL
PL1180-5350	Agilent AdvanceBio SEC 130 Å, 2.7 μ m, 7.8 \times 300 mm
PL1180-3350	Agilent AdvanceBio SEC 130 Å, 2.7 μ m, 7.8 \times 150 mm
PL1180-1350	Agilent AdvanceBio SEC 130 Å, 2.7 μ m, 7.8 \times 50 mm guard
PL1580-5350	Agilent AdvanceBio SEC 130 Å, 2.7 μ m, 4.6 \times 300 mm
PL1580-3350	Agilent AdvanceBio SEC 130 Å, 2.7 μ m, 4.6 \times 150 mm
PL1580-1350	Agilent AdvanceBio SEC 130 Å, 2.7 μ m, 4.6 \times 50 mm guard
5190-9416	Agilent AdvanceBio SEC 130 Å Protein Standard, lyophilized, 1.5 mL

See Also

1. Size exclusion chromatography for Biomolecule analysis: A "How-To" Guide; 5991-3651EN.
2. GPC/SEC standards: Product guide; 5990-7996EN.
3. Calibrating GPC columns: A Guide to Best Practice; 5991-2720EN.

References

1. Critical Reviews in Therapeutic Drug Carrier Systems **1993**, *10(4)*, 307-377.

www.agilent.com/cs/library/slidepresentation/Public/1-Conventional_GPC_-_Polymers_ans_Molecular_Weight.pdf

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Bio SEC-3

資料番号	タイトル
5991-2463EN	Choosing the Right Calibration for the Agilent Bio SEC-3
5991-3954EN	Detailed Aggregation Characterization of Monoclonal Antibodies Using the Agilent 1260 Infinity Multi-Detector Bio-SEC Solution with Advanced Light Scattering Detection
5991-3955EN	Determination of Protein Molecular Weight and Size Using the Agilent 1260 Infinity Multi-Detector Bio-SEC Solution with Advanced Light Scattering Detection
5990-9894EN	Optimum Pore Size for Characterizing Biomolecules with Agilent Bio SEC Columns
5991-0835EN	Development and partial validation of a SEC method for high-resolution separation and quantification of monoclonal antibodies

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Charge Heterogeneity Analysis of Rituximab Innovator and Biosimilar mAbs

Authors

Suresh Babu C.V.
Agilent Technologies India Pvt.
Ltd,
Bangalore, India

Abstract

This Application Note describes the high-resolution separation of charge variants of innovator and biosimilar rituximab using an Agilent 1260 Infinity Bio-inert Quaternary LC, biocolumns, and an Agilent OpenLAB ChemStation Software tool. An Agilent Bio MAb, 4.6 × 250 mm, 5 μm PEEK ion exchange column features a unique resin designed for the charge-based separation of monoclonal antibodies (mAbs). The optimized salt-gradient showed the differences in acidic and basic charge variant profiles between innovator and biosimilar rituximab. Precision of retention time, height, and area of charge isoforms were well within the acceptable range. C-terminal digestion by Carboxypeptidase B (CPB) revealed the major lysine variant peaks in biosimilar rituximab.

Introduction

Recombinant monoclonal antibodies (mAbs) are important biotherapeutics with a wide range of diagnostic and clinical applications. Recently, biosimilar products are increasing in popularity in biopharmaceuticals. mAbs can undergo various post-translational modifications (PTMs) including lysine truncation, deamidation, oxidation, glycosylation, and so forth, becoming heterogeneous in their biochemical and biophysical properties. Due to these modifications, charge variants can affect the efficacy, activity, and stability of mAbs as biotherapeutics. Hence, it is very important to characterize the charge heterogeneity in drug development that will serve as a quality control (QC) step in the biopharmaceutical industry. In addition, precise bioanalytical methods are necessary to demonstrate the similarity between a biosimilar and the innovator product.

Cation exchange chromatography (CEX) is the gold standard for charge-sensitive antibody analysis. In CEX, method parameters often need to be optimized for each protein, as ion exchange depends upon the reversible adsorption of charged protein molecules to immobilized ion exchange groups. This Application Note describes the salt-gradient method for separating the charge variants of innovator and biosimilar rituximab using an Agilent 1260 Infinity Bio-inert Quaternary LC and an Agilent Bio MAb NP5, 4.6 × 250 mm, PEEK ion exchange column. The method compares the CEX profiles of innovator and a rituximab biosimilar. Precision of retention time, height, area, and quantification of acidic, basic, and main forms was determined. Carboxypeptidase B (CPB) digestion was performed to study the contribution of C-terminal lysine variants.

Experimental

Instrumentation

An Agilent 1260 Infinity Bio-inert Quaternary LC, operating to a maximum pressure of 600 bar, was used for the experiments. The entire sample flow path was free of any metal components so that the sample did not come in contact with metal surfaces. Solvent delivery was free of any stainless steel or iron components.

Systems

- Agilent 1260 Infinity Bio-inert Quaternary LC Pump (G5611A)
- Agilent 1260 Infinity Bio-inert High Performance Autosampler (G5667A)
- Agilent 1200 Infinity Series Thermostat (G1330B)
- Agilent 1260 Infinity Thermostatted Column Compartment (TCC) containing bio-inert click-in heating elements (G1316C option 19)
- Agilent 1260 Infinity Diode Array Detector with with 10 mm bio-inert standard flow cell (G1315D)
- Agilent Bio MAb NP5, 4.6 × 250 mm, PEEK (p/n5190-2407)

Software

- Agilent OpenLAB CDS ChemStation Edition, revision C.01.062
- Agilent Buffer Advisor, Rev. A.01.01

Table 1. Chromatographic parameters used for IEX chromatography.

Parameter	Conditions				
Mobile phase A	Water				
Mobile phase B	NaCl (850.0 mM)				
Mobile phase C	NaH ₂ PO ₄ (41.0 mM)				
Mobile phase D	Na ₂ HPO ₄ (55.0 mM)				
Gradient	Time (min)	Mobile phase A (%)	Mobile phase B (%)	Mobile phase C (%)	Mobile phase D (%)
	0	30.3	0	59.6	10.1
	2	26.0	5.0	56.9	12.1
	8	21.5	10.0	54.9	13.6
	20	13.3	19.0	51.9	15.8
	21	30.3	0	59.6	10.1
Injection volume	5 µL				
Flow rate	0.75 mL/min				
Data acquisition	280 nm/4 nm, Ref.: 360 nm/100 nm				
Acquisition rate	5 Hz				
TCC	Room temperature				
Sample thermostat	5 °C				
Post run time	10 minutes				

Reagents, samples, and procedure

Innovator and biosimilar rituximab were purchased from a local pharmacy and stored according to the manufacturer's instructions. Sodium phosphate dibasic, sodium phosphate monobasic, sodium chloride, hydrochloric acid (HCl), and sodium hydroxide were purchased from Sigma-Aldrich. All the chemicals and solvents were HPLC grade, and highly purified water was from a Milli Q water purification system (Millipore Elix 10 model, USA). Carboxypeptidase B (C9584) was purchased from Sigma-Aldrich.

Ion exchange chromatography parameters

Table 1 shows the chromatographic parameters for ion exchange chromatography using a 1260 Infinity Bio-inert Quaternary LC. Rituximab (innovator and biosimilar) were diluted to 1 mg/mL in water, and the elution was monitored at 280 nm. Retention time (RT), area, and percent area were used to calculate standard deviation (SD) and relative standard deviation (RSD %) values. Relative percent area was used to quantify the charge variants of mAbs.

Carboxypeptidase B digestion

Biosimilar and innovator rituximab were diluted to 1 mg/mL using 10 mM sodium phosphate buffer, pH 7.5. To these, 0.25 units of CPB was added and incubated at 37 °C. At various time points, the reaction mixture was aliquoted and quenched with acetic acid before analysis.

Results and Discussion

The Agilent Buffer Advisor Software is an ideal tool to generate pH or ionic strength gradients for protein charge variant separation. It reduces the time required for method development. In this study, a series of method development scouting runs were carried out using the Buffer Advisor Software for optimal mAb charge variant separation. Figure 1 shows the charge variant profiles of innovator and biosimilar rituximab on a Bio MAb PEEK column, demonstrating high-resolution separation of charge variants in 20 minutes with three distinct peaks in biosimilar (Buffer: 30 mM, pH: 6.3, and NaCl: 0–161.5 mM). The Agilent Bio MAb columns contained a highly uniform, densely packed, weak cation exchange resin. Early and late-eluting peaks were called acidic and basic variants, respectively. The peak at 11.4 minutes was designated as the main peak. The overlay of five replicates of innovator and biosimilar rituximab shows excellent separation reproducibility (Figure 2). The average RTs and area RSDs for main peak are shown in the figure. The RSDs are within the acceptable range, which demonstrates the precision of the system.


Figure 1. Charge variant profiles of innovator (A) and biosimilar (B) rituximab using an Agilent Bio MAb 5 μ m column. C) Overlay of innovator and biosimilar rituximab. Peaks 1 and 2: acidic variants; 3: main form; 4, 5, 6 and 7: basic variants.

The high-resolution separation of mAbs facilitated the quantification of charge variants using peak areas. Table 2 summarizes the area percent of charge variants of five consecutive analyses. There was a significant difference in the area percent of the charge variants between two mAbs. The main form in the innovator rituximab was found to be 93.21 % and 29.78 % in biosimilar rituximab. The major charge variant in biosimilar rituximab was 69.46 % basic variants as compared to the innovator product (3.22 %).

Table 2. Charge variants quantification by area %, n = 5.

Innovator – Ristova	RT (min)	Area %
Acidic variant	10.84, 11.21	3.56
Main peak	11.44	93.21
Basic variant	11.9, 12.7	3.22
Biosimilar – Ristova	RT (min)	Area %
Acidic variant	10.73, 11.22	0.76
Main peak	11.45	29.78
Basic variant	11.87, 12.15, 12.59, 13.1, and 13.77	69.46


Figure 2. Overlay of five replicates of innovator (A) and biosimilar (B) rituximab on an Agilent 1260 Infinity Bio-inert Quaternary LC using an Agilent Bio Mab, 4.6 × 250 mm, 5 μm, PEEK column. Insert table shows the precision of retention time and area for main peak, n = 5.

To further characterize the basic variant mAU peaks, both mAbs were subjected to carboxypeptidase B digestion. Figures 3A and 3B show the overlay 50 of the IEX profiles before and after C-terminal cleavage of innovator and biosimilar rituximab, respectively. The disappearance of basic variant peaks after carboxypeptidase B treatment confirmed that the peaks correspond to lysine 20 variants. Figure 4 shows the overlay of the IEX profiles of biosimilar rituximab after CPB treatment and innovator rituximab without CPB treatment, revealing the charge variant similarity between the mAbs.


Figure 3. Characterization of basic charge variants. Separation of CPB treated (overnight) and untreated of innovator (A) and biosimilar (B) rituximab on an Agilent 1260 Infinity Bio-inert Quaternary LC using an Agilent Bio Mab, 4.6 × 250 mm, 5 μm, PEEK column.


Figure 4. Overlay of innovator rituximab without CPB treatment (red) and biosimilar rituximab after CPB treatment (blue).

Conclusion

The salt-gradient method described in this Application Note demonstrates the high-resolution separation of charge variant profiles of mAbs on an Agilent Bio MAb, 4.6 × 250 mm, 5 μm PEEK column. The innovator and biosimilar rituximab had different separation profiles with different degrees of acidic and basic variants. Carboxypeptidase B digestion confirmed that the major basic variant peaks in biosimilar correspond to lysine variants. The Agilent 1260 Infinity Bio-inert Quaternary LC with Bio MAb PEEK columns and reproducible method make this solution particularly suitable for the QA/QC analysis of mAbs for the biopharmaceutical industry.

References

Yan, He; et al. *J. Sep. Sci.* **2011**, 34, 548–555

Agilent publication number 5991-0895EN

Agilent publication number 5990-6844EN

Agilent publication number 5991-0565EN

Optimizing Protein Separations with Cation Exchange Chromatography Using Agilent Buffer Advisor


Protein separation with the Agilent 1260 Infinity Bio-inert Quaternary LC System

Authors

Sonja Schneider
Jochen Strassner
Agilent Technologies, Inc.
Waldbronn, Germany

Abstract

This Application Note shows that the Agilent Buffer Advisor software in combination with the Agilent 1260 Infinity Bio-inert Quaternary LC System is an ideal solution for automated protein separation by ionic strength gradients. Usually, pH scouting using premixed two-component gradients is time-consuming and work-intensive. Dynamically mixed four-component gradients calculated by the Buffer Advisor software shorten and simplify the workflow for pH scouting. In addition, excellent retention time precision and pH consistency were gained using the gradients calculated by the Buffer Advisor software.

Introduction

Proteins consist of many different amino acids comprising weak acidic (carboxylic) and basic (amine) groups. Therefore, proteins are amphoteric molecules that exist mostly as zwitterions in a certain pH range. The pH where the protein has no net charge and does not interact with a charged medium is the isoelectric point (pI). In ion exchange chromatography (IEX), the unique relationship between net surface and pH can be used for optimal protein separation. The pH defines the number of charges on the protein and also helps to stabilize the native structure of the protein in the buffer used during analysis.

To ensure optimal binding and elution characteristics of proteins of interest to the IEX column, pH and ionic strength of the deployed buffer are important factors. Even small changes in these two parameters can affect the separation. As a consequence, pH scouting is an important method to find the optimal separating conditions when working with ionic strength gradients. In contrast to pH-gradients, the pH is kept constant in ionic strength gradients. By increasing the ionic strength (salt concentration) of the mobile phase, the less strongly bound proteins are eluted earlier than the stronger bound proteins.

In general, a premixed two-component gradient is prepared for analysis with a starting buffer of low ionic strength and an elution buffer containing high ionic strength. This includes the following preparation steps:

- Dissolving the appropriate buffering compounds at defined concentration
- Titrating the pH with acid/base to the desired pH of the mobile phase
- Splitting the buffer and adding salt to one portion (elution buffer)
- Titrating the pH of the elution buffer with acid/base to the desired pH, if necessary

To perform pH scouting using premixed two-component gradients, prepare different bottles of buffer. To test, for example, six different pH values, it is necessary to prepare 12 bottles of premixed buffer. In contrast, with dynamically mixed four-component gradients, it is necessary to prepare only four bottles to generate various pH values. Further, dynamical mixing of a buffer eliminates the necessity to titrate the buffer solutions manually, which is typically time-consuming and errors prone. The application of dynamically mixed four-component gradients simplifies method development and reduces the time needed for buffer preparation to a large extent by just providing four bottles with stock solutions:

- Line A: Water
- Line B: Salt solution
- Line C: Acidic buffer component
- Line D: Basic buffer component

Using the four stock solutions, different buffers at different pH and salt concentration can be prepared. The Buffer Advisor software is a helpful tool to calculate the percentages of the stock solutions in order to achieve the desired pH, buffer concentration, and ionic strength (Figure 1).


Figure 1. Quaternary mixing to create a salt gradient with constant pH.

The Buffer Advisor software generates a timetable, which can be imported into the method of the 1260 Infinity Bio-inert Quaternary LC Pump using the Import Solvent Blending File function of the Agilent OpenLAB CDS ChemStation Edition software (Figure 2).


Figure 2. The generated Timetable can be imported into the method of the Agilent 1260 Infinity Bio-inert LC System through the Agilent OpenLAB CDS ChemStation Edition software.

The Buffer Advisor software can be applied for ionic strength or pH gradients in anion or cation exchange chromatography. The software provides a wide choice of different selectable buffers for single buffer (ionic strength gradients) or for composite buffer (pH gradients) applications. Depending on the proteins of interest and the used columns, the user can select buffers either for anion or cation exchange (Figure 3). To ensure optimal buffering capacity, recommended pH ranges and concentrations in which the buffers should be used are displayed.

A common issue in ion exchange chromatography with ionic strength gradients is the decrease in pH as an effect of added neutral salt like NaCl^{1,2}. The Buffer Advisor software counteracts this issue by recalculation of the overall mobile phase composition considering the concentration of acidic and basic buffer (Line C and D) to maintain the desired constant pH. In addition, if the pH deviation gets too large, the Buffer Advisor software automatically inserts additional time points into the pump timetable.

1. Select Buffer & Gradient Mode		2. Define Gradient Table				3. Compose Stock Solution	
<input type="radio"/> Single Buffer (pH / Salt Gradient) <input type="radio"/> Composite Buffer (Wide Range pH Gradient) <input checked="" type="radio"/> Cation Exchange <input type="radio"/> Anion Exchange		Time	Salt	pH	Buffer	A: Water B: NaCl C: MES D: MES - Na	
MES/Na (MES+MES/Na)		0	20	5.8	20		
		5	100	5.8	20		
		30	500	5.8	20		
		35	1000	5.8	20		
		36	20	5.8	20		
Sodium Citrate (Citric + NaOH)		pH 2.9-3.7, 3.7-6.2		7.5-15 mM			
Formic/Na (acid + Na salt)		pH 3.2-4.4		7.5-125 mM			
Formic/Na (acid + NaOH)		pH 3.3-4.6		10-50 mM			
Lactic/Na (acid + Na salt)		pH 3.2-4.5		7.5-125 mM			
Lactic/Na (acid + NaOH)		pH 3.4-4.7		7.5-50 mM			
Acetic/Na (Acetic+Acetate/Na)		pH 3.9-5.4		7.5-125 mM			
Acetic/Na (Acetic+NaOH)		pH 4.1-5.6		7.5-50 mM			
Succinic/Na (acid + Na salt)		pH 3.6-5.6		7.5-125 mM			
Succinic/Na (acid + NaOH)		pH 3.9-6.3		10-20 mM			
Malonic/Na (acid + Na salt)		pH 2.8-5.5		7.5-125 mM			
Malonic/Na (acid + NaOH)		pH 2.9-5.5		7.5-25 mM			
MES/Na (MES+NaOH)		pH 5.5-7.3		7.5-125 mM			
Maleic/Na (acid + Na salt)		pH 2.6-3.5, 5.0-6.2		7.5-125 mM			
Maleic/Na (acid + NaOH)		pH 2.6-3.6, 4.9-6.7		10-20 mM			
ACES/Na (acid + NaOH)		pH 6.1-7.7		7.5-40 mM			
MOPS/Na (acid + Na salt)		pH 6.2-8.1		7.5-125 mM			
MOPS/Na (acid + NaOH)		pH 6.5-8.3		7.5-40 mM			
HEPES/Na (HEPES + salt)		pH 6.6-8.5		7.5-125 mM			
HEPES/Na (HEPES + NaOH)		pH 6.9-8.7		7.5-40 mM			
BICINE/Na (BICINE + Na salt)		pH 7.3-9.1		7.5-125 mM			
BICINE/Na (BICINE + NaOH)		pH 7.6-9.3		7.5-50 mM			
TAPS (acid + Na salt)		pH 7.6-9.4		7.5-125 mM			
TAPS (acid + NaOH)		pH 7.9-9.8		7.5-40 mM			
Sodium borate (H3BO3 + Tetraborate)		pH 8.1-8.9		7.5-125 mM			
Sodium borate (Tetraborate+NaOH)		pH 9.4-10.6		7.5-70 mM			
Sodium borate (H3BO3 + NaOH)		pH 8.4-9.5		7.5-50 mM			
Bicarbonate (NaHCO3+Na2CO3)		pH 9.3-10.2		7.5-125 mM			

Figure 3. Buffer list for cation exchange chromatography, sorted by recommended pH range.

Experimental

Instrumentation

The Agilent 1260 Infinity Bio-inert Quaternary LC System consisted of the following modules:

- Agilent 1260 Infinity Bio-inert Quaternary Pump (G5611A)
- Agilent 1260 Infinity High performance Bio-inert Autosampler (G5667A)
- Agilent 1290 Infinity Thermostat (G1330B)
- Agilent 1290 Infinity Thermostatted Column Compartment (G1316C) with bio-inert solvent heat exchangers
- Agilent 1260 Infinity Diode Array Detector VL (G1315D) with bio-inert standard flow cell, 10 mm)
- Agilent 1260 Infinity Bio-inert Analytical-scale Fraction Collector (G5664A)

Column

- Agilent Bio MAb Column, PEEK, 4.6 × 250 mm, 5 μm

Software

- Agilent OpenLAB CDS ChemStation Edition for LC and LC/MS Systems, Rev. C.01.03 [32]
- Agilent Buffer Advisor, Rev. A.01.01

Solvents

Buffer A: H₂O

Buffer B: NaCl 3 M

Buffer C: MES (2-(N-morpholino) ethanesulfonic acid monohydrate) 60 mM

Buffer D: MES-Na (2-(N-morpholino) ethanesulfonic acid sodium salt) 35 mM

Sample

Mix of three proteins, solved in PBS (phosphate buffered saline), pH 7.4

Ribonuclease A: 13,700 Da pI 9.6

Cytochrom C: 1 2,384 Da pI 10–10.5

Lysozyme: 14,307 Da pI 11.35

All solvents used were LC grade. Fresh ultrapure water was obtained from a Milli-Q Integral system equipped with a 0.22 μm membrane point-of-use cartridge (Millipak). MES (2-(N-morpholino)ethanesulfonic acid monohydrate) and MES-Na (2-(N-morpholino)ethanesulfonic acid sodium salt) were purchased from Merck, Darmstadt, Germany. NaCl was purchased from VWR, Radnor, PA, USA.


Figure 4. Agilent Buffer Advisor software, showing the steps described in "Results and discussion".

Conditions

Parameter	Value
Flow rate:	1 mL/min
Gradient:	0 min – 20 mM NaCl 5 min – 20 mM NaCl 30 min – 500 mM NaCl 35 min – 1000 mM NaCl 36 min – 20 mM NaCl
Injection volume:	10 µL
Thermostat:	4 °C
Temperature TCC:	25 °C
DAD:	280 nm/4 nm
Ref.:	Off
Peak width:	> 0.05 min (1.0 s response time)(5 Hz)

Stock Solution Recipes	
Bottle B	NaCl: Sodium chloride Weight 175.32 g and fill up to 1 L.
Bottle C	MES: MES Weight 11.714 g and fill up to 1 L.
Bottle D	MES - Na: MES sodium salt Weight 7.6027 g and fill up to 1 L.

Figure 5. Stock solution recipes

Results and discussion

pH scouting was performed using calculations from the Agilent Buffer Advisor software for pH values from 5.8 to 6.8. A mix of three proteins (ribonuclease A, cytochrome C and lysozyme) was separated using a four-component salt gradient at six different pH values. Dynamically mixed four-component gradients were generated using the calculations from the software. The Buffer Advisor software simplifies the generation of different four-component gradients by calculating the percentage of the individual stock solutions in the mobile phase at defined time points (Figure 4).

1. After definition of the gradient parameters, such as time, maximum salt concentration, pH, and buffer concentration
2. The Buffer Advisor software calculates the needed stock concentrations. The Recipe button displays the absolute amount of needed chemicals for the preparation of the stock solutions (Figure 5). These proposed stock concentrations can be adjusted by the user.
3. Select the Process tab.
4. The Buffer Advisor software calculates the needed amount of each channel to maintain the correct pH during the complete chromatographic run. Furthermore, it calculates whether the pH, salt concentration and buffer concentration entered is suitable for the buffer system that was selected. The timetable displays also additional data, such as buffering capacity of the mobile phase.

The pH scouting for the three-protein mix of ribonuclease A (A), cytochrome C (B) and lysozyme (C) demonstrates the benefits of the Agilent Buffer Advisor software (Figure 6). Even small pH changes of 0.2 have a strong influence on the retention of the proteins on the weak cation exchange (WCX) column. Changes in the elution order become obvious when the pH is changed from 5.8 to 6.8.

Manual preparation of corresponding buffers for premixed two-component gradients includes several steps. For each pH and for each prepared bottle (one with low and one with high ionic strength), a manual titration of the buffers is necessary. pH scouting for six different pH values in order to achieve the optimal resolution results in preparation of 12 solvent bottles (including weighing chemicals, pH adjustment). This is a very time-consuming procedure and highly prone to error and variation.

In contrast, the Buffer Advisor software is capable of automatically and reproducibly mixing all six separation conditions out of four stock solutions without any manual interference. The optimal resolution was achieved at pH 5.8 (Figure 6).

Based on the results, the user has various options on how to proceed:

1. Fine-tuning of the resulting pH values and gradients
2. Transfer of the dynamically mixed four-component gradient to other instruments through the OpenLAB CDS ChemStation timetable
3. Implementation of pH scouting results into two-component gradients using premixed buffers

In the last case, however, deviations from correct pH are expected, due to the pH optimization procedure of the Buffer Advisor software.


Figure 6. pH scouting for the separation of a three-protein mix using dynamically mixed quaternary gradients.


Figure 7. Off line pH measurement

Conclusions

Using dynamically mixed four-component gradients, calculated by the Buffer Advisor software, shortens and simplifies the workflow for pH scouting. The employment of dynamically mixed gradients calculated with the Buffer Advisor software results in a significant decrease in buffer preparation time, particularly when compared to manual preparation of buffers for premixed two-component gradients. The Buffer Advisor software provides a wide range of prevalidated, user-selectable buffer systems for anion and cation exchange chromatography and delivers recipes for preparation of the most suitable stock solutions. Due to pH optimization of the software, resulting pH values are more accurate and precise than those resulting from premixed gradients formed out of manually prepared buffer solutions. The Buffer Advisor software counteracts this issue by the recalculation of the four-component gradient regarding the concentration of acidic and basic buffer to maintain the desired constant pH.

The Buffer Advisor software in combination with the Agilent 1260 Infinity Bio-inert Quaternary LC System is excellent for generating four-component gradients. The calculations of Buffer Advisor software lead to exact and reproducible protein analysis while providing an excellent tool for automated pH scouting and accurate ion exchange chromatography. The Buffer Advisor software is, therefore, an ideal tool for automatic development of analytical methods in ion-exchange chromatography, which can be seamlessly transferred to the corresponding QA/QC departments.

References

1. R. J. C. Brown & M. J. T. Milton. Observation of a combined dilution and salting effect in buffers under conditions of high dilution and high ionic strength, *Accred Qual Assur* 8(11): 505-510, **2003**.
2. A. E. Voinescu et al. Similarity of Salt Influences on the pH of Buffers, Polyelectrolytes, and Proteins, *J. Phys. Chem. B* 110: 8870-8876, **2006**.

Faster Separations Using Agilent Weak Cation Exchange Columns

Authors

Andrew Coffey
Agilent Technologies, Inc.

Abstract

Ion exchange is a commonly used technique for the separation of complex protein mixtures. Traditionally, such separations are performed using shallow gradients of increasing salt concentration with long column lengths providing the necessary resolution. The columns have often been packed using large diameter particles to minimize backpressure. This Application Note demonstrates how analysis times can be significantly reduced, increasing throughput without compromising analytical performance, by exploiting the benefits of small particle size, non-porous ion exchange sorbents.

Introduction

Proteins, polypeptides and oligonucleotides are often analyzed by ion exchange chromatography because they are complex molecules with multiple charges on their surfaces. The technique is ideally suited to the separation of charged biomolecules as it is non-denaturing and can provide good performance and resolution.

Traditionally, this has meant using highly porous particles to enable such large molecules to permeate the particles. In turn, columns of 15 cm or 25 cm in length, packed with 5 μm or 10 μm particles are commonly used.

The advent of non-porous sorbents such as Agilent's Bio IEX range, comprising a rigid polymeric core particle with a grafted hydrophilic layer containing the ion-exchange functionality, can improve resolution. This is because the diffusion-limited band broadening associated with a molecule penetrating the core of a large particle is eliminated. In turn, this means smaller particles and shorter column lengths can be used to significantly improve throughput, greatly reducing analysis times. The benefits for improved productivity for tasks such as fraction analysis are immediately evident.


Figure 1. Protein separation on Agilent Bio WCX 5 μm 4.6 \times 250 mm versus Agilent Bio WCX 3 μm , 4.6 \times 50 mm (flow rate 1.0 mL/min).

Materials and methods

Agilent Bio IEX columns are packed with rigid polymeric, nonporous particles grafted with a functionalized hydrophilic polymer layer. The resultant 1.7, 3, and 5 μm rigid particles provide high resolution and high separation efficiency by reducing the band broadening effects resulting from diffusion limitations with totally porous particles. The chemically bonded hydrophilic coating significantly reduces the effects of nonspecific binding and results in greater levels of recovery.

Conditions, Bio-Monolith column

Parameter	Value
Columns:	Agilent Bio WCX 5 μm , 4.6 \times 250 mm SS (p/n5190-2445) Agilent Bio WCX 3 μm , 4.6 \times 50 mm SS (p/n5190-2443) Agilent Bio WCX 1.7 μm , 4.6 \times 50 mm SS (p/n5190-2441)
Sample:	Ovalbumin (1), Ribonuclease A (2), Cytochrome c (3), Lysozyme (4)
Instrument:	Agilent 1260 Infinity Bio-inert Quaternary LC
Mobile Phase:	A: 20 mM sodium phosphate, pH 6.5 B: A + 1.6 M NaCl
Gradient:	0 to 50 % B
Temp.:	Ambient
Injection volume:	10 μL
Conc.:	0.5 mg/mL
Detection:	UV, 220 nm

Results and Discussion

The performance of a column, as measured by plate count, is dependent on particle size and column length. From this it may be inferred that a shorter column packed with smaller particles can be used to achieve the same level of performance when compared to a longer column packed with larger particles (Figure 1). This is commonly found in practice. However, for gradient elution, further modifications to the method need to be employed to provide the additional benefits of shorter run times and greater productivity.

Converting gradient times into column volumes is a useful way of calculating the shorter gradient program and can provide the desired outcome in terms of higher speed separations (Table 1). However, smaller particle sizes may require higher flow rates to attain maximum performance. This is illustrated by the van Deemter curves shown in Figure 2.

To maximize the separation efficiency using the Agilent Bio WCX 3 μm , 4.6 \times 50 mm column, the 4 minute gradient separation was carried out at 1.0, 1.5, 2.0, and 2.5 mL/min (Figure 3). As expected, the higher linear velocity created from higher flow rates improved the peak shape.


Figure 2. Typical van Deemter curves.

Table 1. Gradient Time to Column Volume Conversion

Time (minutes)	mM NaC	#CV
0	0	0.0
20	800	4.8
25	800	6.0
25.01	0	6.0
35	0	8.4

#CV = number of column volumes at 1.0 mL/min (4.6 \AA ~ 250 mm column)

Time (minutes)	mM NaCl	#CV
0	0	0.0
4	800	4.8
5	800	6.0
5.01	0	6.0
7	0	8.4

#CV = number of column volumes at 1.0 mL/min (4.6 \AA ~ 50 mm column)

In comparison, the Agilent Bio WCX 1.7 μm , 4.6 \times 50 mm column provided sharper peaks under identical conditions (Figure 4).

Increasing the flow rate should mean that it is possible to further reduce the gradient time. This was investigated using the Bio WCX 1.7 μm , 4.6 \times 50 mm column. The 0 to 800 mM NaCl gradient was reduced from 4 to 2 minutes.

It was found that at a flow rate of 1.7 mL/min the backpressure remained below 400 bar and still provided exceptional peak shape and resolution (Figure 5).


Figure 3. Effect of flow rate on chromatographic performance (Agilent Bio WCX 3 μm , 4.6 \times 50 mm).


Figure 4. Comparison of Agilent Bio WCX 3 μm , 4.6 \times 50 mm versus Agilent Bio WCX 1.7 μm


Figure 5. Agilent Bio WCX 1.7 μm , 4.6 \times 50 mm for protein separations under 3 minutes (flow rate 1.7 mL/min)

Conclusions

We have shown that by using shorter 5 cm columns packed with smaller particle size (3 μm and 1.7 μm), Agilent Bio WCX products can lead to significant reductions in run times from 20 or 30 minutes down to less than 3 minutes, and still retain excellent peak resolution. This enables much higher throughput in time-critical applications.

The backpressure of 400 bar shows that, by reducing the analysis time dramatically from over 30 minutes to less than four minutes for the entire gradient, a 600 bar system such as the Agilent 1260 Infinity Bio-inert LC is still sufficient.

For More Information

These data represent typical results. For more information on our products and services, visit our Web site at www.agilent.com/chem.

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資料番号	タイトル
5991-7442JAJP	Agilent1260 Infinity バイオイナート LC からAgilent1260 Infinity II バイオイナート LC へのシームレスなメソッド移管
5991-5273EN	Characterize mAb Charged Variants by Cation-exchange Chromatography
5991-5274EN	Characterize Fab and Fc Fragments by Cation-exchange Chromatography
5991-0895EN	Analysis of Intact and C-terminal Digested IgG1 on an Agilent Bio MAb 5 µm Column
5990-9629JAJP	pH グラジエント溶出によるモノクローナル抗体電荷変異体の分離の向上
5991-1407EN	High-resolution Analysis of Charge Heterogeneity in Monoclonal Antibodies Using pH-gradient Cation Exchange Chromatography
5991-1408EN	Protein Separation with pH Gradients Using Composite Buffer Systems Calculated by the Agilent Buffer Advisor Software
5991-4722EN	Reducing Cycle Time for Charge Variant Analysis of Monoclonal Antibodies
5991-3365EN	Simple Method Optimization in mAb Charge Variant Analysis using pH Gradients Generated from Buffer Advisor with Online pH and Conductivity Monitoring
5990-9270EN	Separation of Protein Standards on Agilent 3 µm Ion-Exchange Columns by Cation Exchange Chromatography
5990-9614EN	Analysis of proteins by anion exchange chromatography
5991-5221EN	Charge Profiling of 2AB-labelled N-linked Glycans

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© Agilent Technologies, Inc. 2014