




Claim Enhancements for Stripette® Serological Pipets

CORNING

We are pleased to announce three claim enhancements to our Stripette serological pipet portfolio. In the spirit of continuously evolving to support the changing needs of our customers, we are implementing upgrades in the form of increased sterility assurance level, a human DNA-free claim, and availability of extractables data.



Human DNA-free

The serological pipets will now feature a human DNA-free claim.

Pipets will be tested by DNA detection assay per ISO 18385 and found to be free of detectable human DNA. The limit of detection is 10^{-11} g of human genomic DNA.


Extractables data

Extractables analysis will be conducted in accordance with USP <665>. Data will be made available upon customer request.

Sterility Assurance Level (SAL) of 10^{-6}

The serological pipets are sterilized and dosimetrically released per the requirements of ANSI/AAMI/ISO 11137.

The sterility of the upgraded products will be increased from an SAL of 10^{-3} to 10^{-6} .



Ordering Information




Individually Wrapped, Clear Plastic Wrap

Cat. No.	Capacity (mL)	Graduations (mL)	Negative Grads. (mL)	Color-coded Stripe	Qty/Pk	Qty/Cs
4011	1	1/100	0.2	Yellow	50	1,000
4012	1	1/100	0.2	Yellow	50	200
4021	2	1/100	0.2	Green	50	1,000
4051	5	1/10	2.5	Blue	50	200
4101	10	1/10	3.0	Orange	50	200
4492*	10	1/10	3.0	Orange	50	200
4251	25	2/10	10.0	Red	50	200
4501	50	1/2	10.0	Purple	25	100
4484	100	1	n/a	Aqua	10	100



Individually Wrapped, Paper/Plastic Wrap

4485	1	1/100	0.2	Yellow	50	1,000
4486	2	1/100	0.2	Green	50	1,000
4487	5	1/10	2.5	Blue	50	200
4488	10	1/10	3.0	Orange	50	200
4489	25	2/10	10.0	Red	25	200
4490	50	1/2	10.0	Purple	25	100
4491	100	1	N/A	Aqua	10	100



Bulk Packed in Bags

4010	1	1/100	0.2	Yellow	50	1,000
4020	2	1/100	0.2	Green	50	1,000
4050	5	1/10	2.5	Blue	50	500
4100	10	1/10	3.0	Orange	50	500
4250	25	2/10	10.0	Red	25	200
4500	50	1/2	10.0	Purple	25	100

*Cat. No. 4492 features a wide tip for handling viscous fluids.

Related Products

Corning® Stripettor® Ultra Pipet Controller

- ▶ LCD display for continuous visualization of speed, modes, and battery status
- ▶ Colored nose pieces for user or application identification
- ▶ Well-positioned function buttons
- ▶ 3-position charging stand



Cat. No.	Description
4099	Corning Stripettor Ultra pipet controller, with universal power supply, three 0.2 µm filters, 3 colored nose pieces, wall hanger, and 3-position charging stand

Literature Available

- ▶ Corning Stripettor Ultra Pipet Controller (CLS-EQ-112)
- ▶ Liquid Handling Ergonomics in Every Step (CLS-LH-370)
- ▶ Top 10 Ways to Reduce Contamination Risk (MISC-Contamination Graphic-16)
- ▶ Controlling Cell Culture Contamination (PST-Cell-Contamination-16)

Warranty/Disclaimer: Unless otherwise specified, all products are for research use or general laboratory use only.* Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. These products are not intended to mitigate the presence of microorganisms on surfaces or in the environment, where such organisms can be deleterious to humans or the environment. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications. *For a listing of US medical devices, regulatory classifications or specific information on claims, visit www.corning.com/resources.

Corning's products are not specifically designed and tested for diagnostic testing. Many Corning products, though not specific for diagnostic testing, can be used in the workflow and preparation of the test at the customers discretion. Customers may use these products to support their claims. We cannot make any claims or statements that our products are approved for diagnostic testing either directly or indirectly. The customer is responsible for any testing, validation, and/or regulatory submissions that may be required to support the safety and efficacy of their intended application.

CORNING

Corning Incorporated
Life Sciences

www.corning.com/lifesciences

NORTH AMERICA

t 800.492.1110
t 978.442.2200

ASIA/PACIFIC

Australia/New Zealand

t 61 427286832

Chinese Mainland

t 86 21 3338 4338

India

t 91 124 4604000

Japan

t 81 3-3586 1996

Korea

t 82 2-796-9500

Singapore

t 65 6572-9740

Taiwan

t 886 2-2716-0338

EUROPE

CSEurope@corning.com

France

t 0800 916 882

Germany

t 0800 101 1153

The Netherlands

t 020 655 79 28

United Kingdom

t 0800 376 8660

All Other European Countries

t +31 (0) 206 59 60 51

LATIN AMERICA

grupoLA@corning.com

Brazil

t 55 (11) 3089-7400

Mexico

t (52-81) 8158-8400